



Feb 28, 2023

Whole-cell radioligand binding for receptor internalization

 [PLOS One](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.x54v9dz91g3e/v1>

Angus Li¹, Samuel Liu¹, Rennica Huang², Seungkirl Ahn¹, Robert J Lefkowitz^{1,2,3}

¹Department of Medicine, Duke University Medical Center, Durham, North Carolina, United States of America;

²Department of Biochemistry, Duke University Medical Center, Durham, North Carolina, United States of America;

³Howard Hughes Medical Institute, Duke University Medical Center, Durham, North Carolina, United States of America



Angus Li

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.x54v9dz91g3e/v1>

External link: <https://doi.org/10.1371/journal.pone.0283477>

Protocol Citation: Angus Li, Samuel Liu, Rennica Huang, Seungkirl Ahn, Robert J Lefkowitz 2023. Whole-cell radioligand binding for receptor internalization. protocols.io <https://dx.doi.org/10.17504/protocols.io.x54v9dz91g3e/v1>

Manuscript citation:

Li A, Liu S, Huang R, Ahn S, Lefkowitz RJ (2023) Loss of biased signaling at a G protein-coupled receptor in overexpressed systems. PLoS ONE 18(3): e0283477. doi: [10.1371/journal.pone.0283477](https://doi.org/10.1371/journal.pone.0283477)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 28, 2023

Last Modified: February 28, 2023

Protocol Integer ID: 77746

Keywords: receptor internalization this protocol detail, receptor internalization, cell radioligand, receptor, coupled receptor, overexpressed system

Funders Acknowledgements:

Howard Hughes Medical Institute

National Institutes of Health

Grant ID: R01 HL16037

Abstract

This protocol details an experimental procedure used to generate results described in the manuscript Li, A., Liu, S., Huang, R., Ahn, S., & Lefkowitz, R. J. (2023). Loss of biased signaling at a G protein-coupled receptor in overexpressed systems.



Day 1

- 1 Grow cells on 150 mm dish to ~70% confluency
- 2 Wash twice with 10 mL PBS
- 3 Detach cells with 1.5 mL trypsin-EDTA + 10.5 mL media. Trypsinize for 5 minutes and check under microscope for complete detachment; tap dish if necessary
- 4 Count cells with hemacytometer and dilute collected cells to 164,000 cells/mL
- 5 Pipet 1 mL into each well of poly-D-lysine coated 12-well plate

Day 2

- 6 Add 250 mL of 5x doxycycline for all desired concentrations to cells 14 hours before assay

Day 3

- 7 Prepare > 300 uL ligand-containing media per well (e.g. 10 uM AngII) for stimulation
- 8 Remove media from plate wells by vacuum in cell culture hood
- 9 Add 300 uL ligand-containing media per well and return to incubator for 1 hour
- 10 Prepare saturating concentration of membrane-impermeable radioligand in media on ice. A ligand with high specific activity is preferred. For example, 30 uL 3H-AngII in 12 mL media. Prepare at least 4 mL per plate. Vortex briefly to mix
- 11 Take 1/3 of media with radioligand and add high concentration (orders of magnitude above radioligand concentration) of nonlabelled antagonist (e.g. 200 uM candesartan



- final concentration) for nonspecific binding measurement. Vortex briefly to mix
- 12 After 1 hour ligand incubation, place plates on ice and remove media by vacuum. Wash twice with 1 mL ice-cold PBS with MgCl₂ and CaCl₂ per well.
 - 13 Strip ligand for 5 minutes with 2 mL acid wash buffer (0.2 M acetic acid, 0.5 M NaCl, pH 2.5) per well
 - 14 Remove acid wash buffer by vacuum and wash 2-3 times with 1 mL ice-cold PBS with MgCl₂ and CaCl₂ per well.
 - 15 Remove wash by vacuum. Add 300 uL media with radioligand per well on ice. Use top two rows of plate for total binding duplicates and bottom row for nonspecific binding. Wait ≥ 1.5 hours on ice
 - 16 Pre-chill 50 mL per plate PBS with MgCl₂ and CaCl₂ on ice and pre-warm 0.1% SDS + 0.5 N NaOH in 42 °C water bath
 - 17 Remove radioligand media into dedicated radioactive liquid waste vial using P1000 on ice. Remove nonspecific binding media last
 - 18 Wash wells with 1 mL cold PBS with MgCl₂ and CaCl₂ on ice and collect using P1000 into radioactive liquid waste vial. Wash 4 times total. Check plates under microscope after last wash
 - 19 Remove plates from ice and add 500 uL warm 0.1% SDS + 0.5 N NaOH per well. If there is precipitate in the solution, allow it to warm for longer. After adding, shake for ≥ 15 minutes
 - 20 Collect samples from wells using P1000 and add to labeled scintillation vials. Add 300 uL of each radioligand media mix to two scintillation vials each for total counts. Add 5 mL scintillation cocktail (Lefko-Fluor) to each vial. Cap vials
 - 21 Wipe vials with moist KimWipe and place in racks for scintillation counter. Place racks in counter and count using 3h_cpm_3_min protocol
 - 22 Perform second count using different protocol flag on next day

Analysis

- 23 Copy binding counts into Excel spreadsheet



- 24 Calculate nonspecific counts for each receptor concentration
- 25 Calculate specific counts for each condition by averaging duplicates and subtracting nonspecific count
- 26 Calculate percent internalization as percent reduction in specific counts for ligand-stimulated condition compared to control