

Apr 15, 2025

Whole genome amplification of human parechovirus type 3 (PEV-A3) utilizing tiling-PCR

DOI

<https://dx.doi.org/10.17504/protocols.io.q26g758q3lwz/v1>

Nora Deezsi-Magyar¹

¹National Center for Public Health and Pharmacy



Nora Deezsi-Magyar

National Center for Public Health and Pharmacy

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g758q3lwz/v1>

Protocol Citation: Nora Deezsi-Magyar 2025. Whole genome amplification of human parechovirus type 3 (PEV-A3) utilizing tiling-PCR. protocols.io <https://dx.doi.org/10.17504/protocols.io.q26g758q3lwz/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 15, 2025

Last Modified: April 15, 2025

Protocol Integer ID: 126707

Keywords: whole genome amplification of human parechovirus type, whole genome amplification of pev, human parechovirus type, pev, pcr the procol, whole genome amplification, pcr, genome, sequencing

Abstract

The procol allows whole genome amplification of PEV-A3, prior to next generation library preparation and sequencing. Suitable for Nanopore and Illumina sequencing as well.



DNase I treatment (to remove contaminating human-derived gDNA)

40m

1 Reagents

- TURBO DNaseTM (2 U/ μ L), Invitrogen
- Nuclease-free destillated water

1.1 Prepare the following master mix

	A	B
	TURBO DNase (2 U)	1.0 μ L
	TURBO DNase puffer	4.0 μ L
	Destillated water	6.0 μ L
	RNA	40.0 μL

1.2 Incubation in thermocycler

37 °C

00:30:00

30m

2 Post-DNase I treatment clean-up

Reagents

- AMPure XP Beads, Beckman Coulter
- 70% EtOH (prepared freshly)
- Nuclease-free destillated water

2.1 Add 1.8 μ L AMPure XP per 1.0 μ L of sample.

2.2 Bind DNA fragments to paramagnetic beads.

2.3 Separation of beads + DNA fragments from contaminants.

2.4 Wash beads + DNA fragments twice with 70% EtOH to remove contaminants.



2.5 Elute purified DNA fragments from beads.

2.6 Transfer to new tubes.

Reverse transcription

1h 20m

3 Reagents

- SuperScript[™] VILO[™] cDNA Synthesis Kit, Invitrogen
- Nuclease-free destillated water

3.1 Prepare the following master mix

	A	B
	5X VILO [™] Reaction Mix	4.0 uL
	10X SuperScript [™] Enzyme Mix	2.0 uL
	Destillated water	4.0 uL
	RNA	8.0 uL

3.2 Incubation in thermocycler

🌡️ 25 °C ⌚ 00:10:00

10m

3.3 Incubation in thermocycler

🌡️ 42 °C ⌚ 01:00:00

3.4 Incubation in thermocycler

🌡️ 85 °C ⌚ 00:05:00

5m



Polymerase chain reaction (tiling-PCR)

2h 30m

4 Reagents

- Phusion™ High-Fidelity DNA Polymerase, New England Biolabs
- Nuclease-free destillated water
- dNTP mix

4.1 Prepare the following master mixes

Mix_1

	A	B
	5x Phusion [®] High-Fidelity PCR Master Mix with GC Buffer	4.0 uL
	dNTP mix (10 mM)	0.4 uL
	Phusion [®] High-Fidelity DNA Polymerase	0.2 uL
	Primer pool_2	3.0 uL
	Destillated water	7.4 uL
	cDNA template	5.0 uL

The mix includes pool_1 primers (see section 4.2)

Mix_2

	A	B
	5x Phusion [®] High-Fidelity PCR Master Mix with GC Buffer	4.0 uL

A	B
dNTP mix (10 mM)	0.4 uL
Phusion [®] High-Fidelity DNA Polymerase	0.2 uL
Primer pool_2	3.0 uL
Destillated water	7.4 uL
cDNA template	5.0 uL

The mix includes pool_2 primers (see section 4.2)

4.2 Tiling-PCR primers

A	B	C	D	E
Primer	Pool	Sequence (5' - 3')	Size	Position (bp)
hPeV-3_1_LEFT	1	ATACCCCGATTGCTGAGC TTC	22	41 – 987
hPeV-3_1_RIGHT	1	GCCATTGAATGAAAGTTAT CCACATTTC	28	
hPeV-3_2_LEFT	2	GATGTAGTGCAAGCTACGA CCA	22	899 – 1895
hPeV-3_2_RIGHT	2	CTACTGTGGCTGCCCAATC AAA	22	
hPeV-3_3_LEFT	1	TGATCCTAGAACTGCAGGG AGT	22	1750 – 2745
hPeV-3_3_RIGHT	1	TCTTCAGTGTCATATGTATG AGCCAC	26	
hPeV-3_4_LEFT	2	AAGAGGGTCATGGCATGTT GTC	22	2625 – 3631

	A	B	C	D	E
	hPeV-3_4_RIG HT	2	TGTAGACAAACAAGCAGTG GTTAGA	25	
	hPeV-3_5_LEF T	1	TCCTCAGCAGCCACAGAA ATTG	22	3500 – 4487
	hPeV-3_5_RIG HT	1	TCTCTTCAAGATGTGCCAT TGGG	23	
	hPeV-3_6_LEF T	2	GCCAGTGAGTTCATGGATG GTT	22	4340 – 5339
	hPeV-3_6_RIG HT	2	GCTGGTATCCTGCACAATG TGT	22	
	hPeV-3_7_LEFT	1	GCCAAACCAAAGAGTGCT TTCC	22	5192 – 6188
	hPeV-3_7_RIGH T	1	CGGACTTAACAAAGCTGTA CCCT	23	
	hPeV-3_8_LEF T	2	GGCAAGGTGTTAAAGCATG TGTC	23	6039 – 7057
	hPeV-3_8_RIG HT	2	GCTGCTTGAATGTGCTGAA GTTT	23	
	hPeV-3_9_LEF T	1	ACCACCATCTTTAACACTT GTCTCA	25	6296 – 7323
	hPeV-3_9_RIG HT	1	TGGTATGTCCAATATTCCA AATTAGTGTTT	30	

4.3 Tiling-PCR thermocycling

6m 50s

- 1 98 °C 00:00:30
- 2 98 °C 00:00:10
- 3 54 °C 00:00:20
- 4 72 °C 00:00:50
- 5 Repeat steps 2-3-4 for 40 times.
- 6 72 °C 00:05:00



7  4 °C hold

4.4 Post PCR-clean-up

Reagents

- Ampure XP Beads, Beckman Coulter
- 70% EtOH (freshly prepared)
- Nuclease-free distilled water

4.5 Add 0.8 µl AMPure XP per 1.0 µl of sample.

4.6 Bind DNA fragments to paramagnetic beads.

4.7 Separation of beads + DNA fragments from contaminants.

4.8 Wash beads + DNA fragments twice with 70% EtOH to remove contaminants.

4.9 Elute purified DNA fragments from beads.

4.10 Transfer to new tubes.

Tiling-amplicon quantification using Qubit Flex fluorometer

10m

5 Reagents

- Qubit 1X dsDNA HS Assay Kit

5.1 Add 190 µL of working solution to all wells of two Qubit strips. Label the top of the tubes S1 and S2. Then add 10 µL of the respective HS Standard (1 or 2) to each tube of the two strips.

5.2 For each of your sample tubes, add 198 µL of working solution to labeled tubes. Then add 2 µL of your samples to each tube.

5.3 Vortex all strips and spin down if necessary. Incubate for 3-5 minutes at room temperature.



- 5.4 On the Qubit, select 1x dsDNA → dsDNA High Sensitivity → Read standards
- 5.5 Insert strip standard 1 and then select Read standard. Repeat for standard 2.
- 5.6 Select Run samples and read the respective wells of the strips. Use the 2 uL sample volume. When you are done.