

May 02, 2025

Whole-genome STARR-seq Protocol – Reddy Lab 2023

DOI

<https://dx.doi.org/10.17504/protocols.io.3byl4z7povo5/v1>

Graham Johnson¹, Kari Strouse¹, Keith Siklenka¹, Shauna Morrow¹, Tim Reddy¹

¹Duke University

IGVF



Alejandro Barrera

Duke University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl4z7povo5/v1>

Protocol Citation: Graham Johnson, Kari Strouse, Keith Siklenka, Shauna Morrow, Tim Reddy 2025. Whole-genome STARR-seq Protocol – Reddy Lab 2023. [protocols.io https://dx.doi.org/10.17504/protocols.io.3byl4z7povo5/v1](https://dx.doi.org/10.17504/protocols.io.3byl4z7povo5/v1)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 27, 2025

Last Modified: May 02, 2025

Protocol Integer ID: 126117

Keywords: qPCR, DNA library, Electrocompetent cells, Reverse transcription, STARR-seq, genome starr, whole genome starr, seq protocol, seq experiment from the reddy lab, genome, seq experiment, reddy lab 2023 reddy, output library generation, protocol product, reddy lab, seq, protocol, encode project, expected protocol product, transfection

Funders Acknowledgements:

IGVF

Grant ID: HG012053

CEGS

Grant ID: HG011123

Abstract

Reddy lab whole genome STARR-seq protocol. The document includes the input, transfection and output library generation, which should be the expected protocol products. Adapted from the document found in STARR-seq experiment from the Reddy lab in the ENCODE project (document attached with additional details).

Attachments



wgSTARRseq_ENCODE_

KP...

230KB



Materials

-  NEBNext DNA Library Prep Master Mix Set for Illumina - 60 rxns **New England Biolabs Catalog #E6040L**
- KAPA HiFi HotStart PCR Kit (Cat. # KK2502)
-  SYBR™ Safe DNA Gel Stain **Thermo Scientific Catalog #S33102**
-  GeneJET Gel Extraction Kit **Thermo Fisher Catalog #K0692**
-  NEBuilder HiFi DNA Assembly Master Mix **New England Biolabs Catalog #E2621**
- MN Nucleobond PC 10000 EF (Cat. #740548)
-  RNase Block **Agilent Technologies Catalog #300151**
- Dynabeads mRNA Purification Kit (Cat.# 61005)
-  TURBO DNA-free™ Kit **Thermo Scientific Catalog #AM1907**
- SuperScript III system (Invitrogen Cat.#18080)
-  Ambion™ Recombinant RNase A **Thermo Fisher Scientific Catalog #AM2269**

Part 1. INPUT LIBRARY GENERATION

18h 25m

1 Day 1: Prepare STARR-seq inserts

1.1 Shear DNA to ~500 bp.

1.2 Using 8 rxns with  5 µg of DNA per rxn (total  40 µg DNA per genome), follow the NEBNext DNA Library Prep Master Mix Set for Illumina (Cat. #E6040L) through the End - Repair and A-tailing steps.

1.3 Perform the ligation step with  15 micromolar (µM) STARR-seq adaptors instead of the supplied adaptor (final concentration  3 micromolar (µM)).

Note

STARR-seq adaptors:

Adaptor1: ACACTCTTTCCCTACACGACGCTCTTCCGATC*T

Adaptor2: [Phos]GATCGGAAGAGCACACGTCTGAACTCCGATC*T

1.4 Following the first SPRI cleanup after the ligation step, perform the 500 bp insert SPRI size selection. Elute with  30 µL of EB and collect  27 µL of sample.

1.5 Using the KAPA HiFi HotStart PCR Kit (Cat. # KK2502) amplify the STARR-seq library:

- Primers: TS-2-SS-F + TS-2-SS-R
- Reaction conditions per sample:

	A	B
		ul per sample
	HiFi polymerase	1
	GC 5x buffer	10
	10 uM primer mix	5 + 5

	A	B
	10 mM dNTPS	1.5
	DNA library	25
	H2O	2.5
	Sum	50

- Cycling conditions:

	A	B	C
	95°C	3 min	
	98°C	20 s	10 cycles
	63°C	15 s	
	72°C	1 min	
	72°C	5 min	
	4°C	Hold	

1.6 Pool enriched libraries and perform a 0.9x SPRI cleanup. Elute in  100 µL .

1.7 Assess size distribution on the Tape Station and concentration with Qubit Assay.

1.8 Store at  -20 °C .

2 Day 2: Linearize STARR-seq Vector

Note

If using previously linearized vector, run on a gel to check for degradation and double check concentration on Qubit.

2.1 Digest suitable amount of supercoiled STARR-seq plasmid with AgeI and SalI.

- 2.2 Run restriction products on a large 1% agarose gel and 1kb ladder.
- 2.3 Stain gel with SybrSafe (Thermo #S33102) on an orbital shaker, covered with foil.
- 2.4 Excise large band (~2500 bp) with clean razor blade and sterile tweezers.
- 2.5 Recover fragment with GeneJet Gel Extraction Kit (Thermo #K0692).
- 2.6 Assess size distribution on the Tape Station and concentration with Qubit Assay.

3 Day 3: Gibson Assembly

- 3.1 Pre-warm thermocycler to  50 °C .
- 3.2 Mix the following  On ice for a **[M]** 10 nanomolar (nM) Gibson reaction:
 - NEBuilder HiFi DNA Assembly Master Mix (Cat. #E2621)
 - For whole-genome STARR-seq, use 20 reactions per genome



A	B
	per reaction
Total DNA fragments	0.2 pmols
Insert:Vector molar ratio	3:1
NEBuilder HiFi DNA Assembly MM	10uL
H2O	up to 20uL total

- 3.3 Distribute up to  100 µL of the MM into each well.



3.4 Incubate for 00:15:00 - 00:30:00 at 50 °C .

30m



3.5 Immediately place the reactions On ice and add EDTA to 10 millimolar (mM) .



- Ex. 2 μL 0.5 Mass Percent EDTA to 100 μL volume

3.6 Pool sample reactions and begin an EtOH precipitation by adding 0.1X volume 3 Mass Percent NaOAc and 2.5X volume very cold 100% EtOH.

3.7 Incubate at -20 °C Overnight .

8h



4 Day 4: Transformation Gibson product

4.1 Centrifuge samples at 16000 x g, 4°C, 00:30:00 .

30m



4.2 Wash pellet with 500 μL cold 70% EtOH, centrifuge at 16000 x g, 4°C, 00:10:00 . Repeat once.

10m



4.3 Let DNA pellet air dry for a few minutes, then resuspend in 40 μL H₂O. Continue with transformation or store at -20 °C .

4.4 Split gibson product across 8 cuvettes per genome, each with 300 μL of in-house grown and frozen Endura electrocompetent cells.

4.5 Transform on the BioRad electroporator with the preset protocol: *E. coli* 2mm 3kV.

4.6 Immediately add 1 mL Lucigen Recovery Media warmed to 37 °C . Wash cuvettes 4X and collect in a tube for a total of 5 mL . Recover by shaking for

1h



 01:00:00 at  37 °C .

4.7 Plate diluted culture on carbenicillin plates to determine transformation efficiency (puc19) and estimate library diversity the following day.

4.8 Pool transformations into  1 L LB with carbenicillin (4 cuvettes per flask) and shake for  08:00:00 at  37 °C .

8h

4.9 Pellet transformed cells  5000 x g, 4°C, 00:15:00 and store at  -20 °C .

15m



5 Day 5: Harvest plasmid DNAs and prepare QC sequencing libraries

5.1 Harvest the STARR-seq plasmid library using 1 Giga-prep column per 2 flasks.

- MN Nucleobond PC 10000 EF (Cat. #740548)

5.2 Determine product concentration with Qubit assay.

5.3 Amplify STARR-seq input sequencing libraries off of plasmid DNAs using the Kapa HiFi HotStart PCR Kit (Cat.# KK2502).

- Primers: 208 TruSeq primer + i7 Index primer or i5 Index primer + i7 PCR primer
- Reaction conditions per sample:

	A	B	C
		ul per sample	x5.5
	20ng plasmid DNA	-	-
	GC buffer	5	27.5
	10 mM dNTPS	0.75	4.13
	10 uM primers	0.75 + 0.75	4.13 + 4.13
	Polymerase	0.5	2.75



A	B	C
H ₂ O	to 25uL	to 137.5uL

- Cycling conditions:

A	B	C
98°C	45 s	
98°C	15 s	9 cycles
64.5°C	30 s	
72°C	1 m. 10 s	
72°C	1 m	
4°C	Hold	

- 5.4 Perform a double-sided bead cleanup (0.5X, 0.9X). Elute with  30 µL of water.
- 5.5 Assess size distribution on the Tape Station and concentration with Qubit Assay.
- 5.6 If running multiple samples, pool and dilute to make  20 µL of a 4nM pool.
- 5.7 Run a MiSeq sequencing run to estimate library quality and diversity before sequencing deeper.

Part 2. TRANSFECTION & OUTPUT LIBRARY GENERATION

4h 16m

6 Day 6: Transfection

- 6.1 Transfect an appropriate amount of the STARR-seq library into the desired cell type using the Lonza 4D Nucleofector and a recommended kit for that cell line.



6.2 Harvest cells 6 hours post-transfection, add RLT buffer (from Qiagen RNeasy kit) with 1% beta mercaptoethanol, and vortex.

6.3 Flash freeze in 50mL conical tubes and store at  -80 °C .

7 Day 7: RNA Isolation

7.1 Pass thawed sample through a 18- to 20-guage needle 10 times.

7.2 Isolate RNAs through Qiagen RNeasy midi-columns, performing the on-column DNase step.

Note

can increase speed and decrease time - 1min instead of 5min

7.3 Elute twice using  125 µL and  75 µL RNase-Free H₂O.

7.4 Quantify RNA concentration with Qubit BR and assess R.I.N. values with RNA Tape.

7.5 Add  1 µL RNase Block (Agilent Cat.# 300151).



7.6 Store at  -80 °C .

8 Day 8: RNA-seq Library Construction Pt.1

8.1 **Prepare 50 mLs of each capture buffer fresh using nuclease-free stocks**

- 2x Binding buffer:



A	B
Tris-HCl pH 7.5	20 mM
LiCl	1.0 M
EDTA	2 mM

- Wash Buffer B:

A	B
Tris-HCl pH 7.5	10 mM
LiCl	0.15 M
EDTA	1 mM

- 10 mM Tris-HCl pH 7.5

8.2 Prepare RNA

2m

- Use the Dynabeads mRNA Purification Kit (Cat.# 61005).
- Adjust the volume of  75 µg total RNA to  100 µL with
 10 millimolar (mM) Tris-HCl, pH 7.5

Note

If using less than 75 µg RNA scale bead volume (below; 1.0 mg), but keep other volumes constant.

- Heat RNA to  65 °C for  00:02:00 to disrupt secondary structures. Place
 On ice .

8.3 Prepare Dynabeads



- Transfer  200 µL ( 1.0 mg) of resuspended Dynabeads to a micro centrifuge tube.
- Place the tube on the magnet and discard the supernatant.



- Remove tube from the magnet and add  100 μL 2x Binding Buffer.
- Place the tube on the magnet and discard the supernatant.
- Remove tube from the magnet and add  100 μL 2x Binding Buffer to the Dynabeads.

8.4 Isolate mRNA

7m



- Add the total RNA to the Dynabeads/2x Binding Buffer suspension.

Note

Optimal hybridization conditions are obtained in 2x Binding Buffer added in a 1:1 ratio relative to sample volume.

- Rotate for  00:05:00 at  Room temperature or mix occasionally by hand.
- Place the tube on the magnet and discard the supernatant.
- Remove the tube from magnet and wash the mRNA-bead complex twice with  200 μL Washing Buffer B.
- Elute with  100 μL of  10 millimolar (mM) Tris-HCl, pH 7.5. Heat to  75 $^{\circ}\text{C}$ for  00:02:00 and place the tube immediately on the magnet.
- Transfer the eluted mRNA to a new RNase-free tube  On ice . Save beads for the next step.

8.5 Wash Dynabeads for second capture



- Add  50 μL 2x binding buffer and  50 μL RNase-Free H₂O. Remove supernatant and repeat once.
- Remove supernatant and tube from the magnet.
- Add  100 μL 2x Binding Buffer to beads.

8.6 Perform a second mRNA capture

9m

- Heat the eluted RNA to  65 $^{\circ}\text{C}$ for  00:02:00 to disrupt secondary structures. Place  On ice .
- Add the  100 μL of RNA to the Dynabeads suspension.

- Rotate for 00:05:00 at Room temperature or mix occasionally by hand.
- Place the tube on the magnet and discard the supernatant.
- Remove the tube from magnet and wash the mRNA-bead complex twice with 200 μL Washing Buffer B.
- Elute with 40 μL of [M] 10 millimolar (mM) Tris-HCl, pH 7.5. Heat to 75 $^{\circ}\text{C}$ for 00:02:00 and place the tube immediately on the magnet.
- Transfer 40 μL of mRNA to a new RNase-free tube On ice .

8.7 Remove contaminating plasmid and genomic DNA

38m

- Prepare DNase digestion master mix using the Turbo DNase (Ambion # AM1907).



A	B
	μL per sample
10 $\times$ TURBO DNase Buffer	5.0
TURBO DNase	1.0
RNase Block	1.0
H ₂ O	3.0
Sum:	10 μL

- Add 10.0 μL master mix to each well (40 μL mRNA) for a final volume of 50 μL .
- Incubate at 37 $^{\circ}\text{C}$ for 00:30:00 .
- Add 5 μL resuspended DNase Inactivation Reagent and mix well.
- Incubate 00:05:00 at Room temperature , mixing occasionally.
- Centrifuge at 10000 $\times$ g, 00:01:30 .
- Transfer 40 μL RNA to a fresh tube, make sure bead slurry does not carry over.
- Add 35 μL of H₂O to each tube/well with beads to perform a back extraction with a final volume of $\sim$ 50 μL .
- Centrifuge at 10000 $\times$ g, 00:01:30 .
- Transfer 27 μL of supernatant to matched supernatant from first spin, make sure bead slurry does not carry over.

- Add  1 μL of Agilent RNase Block (40 units/ μL), vortex or pipette to mix for a final volume of  68 μL .
- **Optional:** Quantify RNA and DNA concentration (Qubit). Replace volume with H₂O for a final volume of  68 μL .

8.8 Reverse Transcription

2h 20m



- Use the SuperScript III system (Invitrogen Cat.#18080) and a reporter-specific primer (RT_UMI) to perform first-strand cDNA synthesis.
 1. Each sample will be primed in a single reaction using the first master mix.
 2. In the second reaction, the sample will be split into two unequal amounts; 5/6 ( 65 μL) will be reverse-transcribed and 1/6 ( 13 μL) will serve as a no-RT control.
- Prepare first master mix.

A	B
	μl per sample
4 uM RT_UMI Primer	3.0
10 mM dNTPs	6.0
RNase Block	1.0
Sum:	10 μl

- Add  10 μL master mix to each RNA sample (~  68 μL) for a final volume of  78 μL .
- Heat samples to  65 $^{\circ}\text{C}$ for  00:05:00 and incubate  On ice for at least 2 minutes.
- From each sample, remove  13 μL to a new well. This is your RT- control.
- Prepare the RT+ master mix.

A	B
	μl per sample
5X 1st Strand Buffer	20.0
0.1 M DTT	5.0



A	B
RNase Block	5.0
SuperScript III (200U/ul)	5.0
Sum:	35 µl

- Add  35 µL master mix to each primed RT+ sample ( 65 µL) for a final volume of  100 µL .
- Prepare the RT- master mix.

A	B
	µl per sample
5X 1st Strand Buffer	4.0
0.1 M DTT	1.0
RNase Block	1.0
H2O	1.0
Sum:	7 µl

- Add  7 µL master mix to each primed RT- sample ( 13 µL) for a final volume of  20 µL .
- Incubate at  50 °C for  02:00:00 and inactivate the reaction by heating at  70 °C for  00:15:00 .
- This is a good stopping point. Transfer samples to  -20 °C .

9 Day 9: RNA-seq Library Construction Pt.2

9.1 RNase digestion

- Add  1 µL DNase-free recombinant RNase A (Thermo #AM2269) and incubate at  37 °C for  01:00:00 .

1h



Note

For the RT- samples dilute the RNase 1:5 in H₂O and add 1 µl to each well.

- Perform a 1.5x SPRI bead cleanup
- 1. Add  150 µL of beads to the RT+ well
- 2. Add  30 µL of beads to the RT- well.
- Wash 2x with 80% EtOH.
- Elute the RT+ well with  145 µL H₂O and transfer  20 µL of each sample to the first seven wells of a column of a new plate. Each column in this new plate corresponds to a sample.
- Elute the RT- well with  22 µL H₂O and transfer  20 µL to the eighth well of the column containing the control's corresponding RT+ samples.

9.2 PCR enrichment

1. One master mix will be made for each sample and corresponding RT- reaction. The RT+ portion of the sample will be enriched in 7 individual PCRs and the RT- sample in an 8th reaction. Therefore, prepare enough master mix for ≥8.5 reactions:

Note

Reactions must be setup on ice.

A	B
	µl per sample (RT+ and RT-)
5x GC Buffer	85.0
10 uM i7PCR primer	42.5
10 mM dNTPs	12.75
1 U/µL Kapa HiFi HotStart Polymerase	8.5
H ₂ O	63.75
Sum:	212.5 µl



2. Add  25 μL master mix to each well in a column.
3. Add  5 μL of  10 micromolar (μM) i5 index primer to each well in a column.
Use a different primer for each column/sample.
4. Enrich cDNAs using the following protocol:

A	B	C
98°C	45 s	
98°C	15 s	5 cycles
64.5°C	30 s	
72°C	1 m 10 s	
4°C	Hold	

5. Setup qPCR survey reactions using  5 μL of pre-enriched sample (from first well of a sample column) as template to determine the optimal number of PCR cycles needed for a sample. The RT- well is not sampled.

Note

Keep the pre-enriched plate covered on ice while the survey reaction runs. The polymerase is active after the hot-start and if let at room temp it will degrade the primers.

A	B
	μl per sample
5x GC Buffer	2.0
10 μM i7PCR primer	1.0
10 mM dNTPs	0.3
1 U/ μL Kapa HiFi HotStart Polymerase	0.2
100x Sybr or 20X EvaGreen	0.09 or 0.75
H2O	5.41 or 4.75
	9 μl



6. For each sample, add  9 μL of mix,  5 μL of sample, and  1 μL of the corresponding i5 index primer to a well of a qPCR optical plate for a final volume of  15 μL .

7. qPCR cycling conditions:

	A	B	C
	98°C	45 s	
	98°C	15 s	45 cycles
	64.5°C	30 s	
	72°C	1 m 10 s	

8. Plot signal (R) vs. cycle number. Calculate the number of additional PCR cycles needed for each sample by determining the number of cycles needed to reach 1/3 of the maximum R. If running many samples and Ct's are 'reasonably' similar, calculate the median.

9. Continue PCR on partially-amplified samples for the appropriate number (N) of cycles:

	A	B	C
	98°C	45 s	
	98°C	15 s	N cycles
	64.5°C	30 s	
	72°C	1 m 10 s	
	72°C	1 m	
	4°C	hold	

10. Perform a 0.9x ( 45.0 μL) SPRI bead cleanup.

Note

can pool 7 PCR wells into a deepwell plate

11. Elute the seven RT+ wells with  10 μL EB and collect  8 μL per well, pooling elutions from the same column. The final volume should be 56 μL of library per sample.

12. Elute the No-RT wells with  14 µL EB. Keep this solution in the well with the beads bound to the magnet. Draw directly from the well when preparing to run a tape in the next step.

Note

Eluting the no-RT control in this volume preserves the initial ratio of the control being 1/6 of the total captured sample (14 µl / (14 + 70 µl)). Eluting in a greater amount could dilute unwanted product.

13. Run RT+ and RT- products on the tape station with a D1000 tape.
14. **Optional**, adjust volume to 0.05% Tween for longer term storage of library.
15. Quantify DNA concentration with the Qubit system.
16. Run a MiSeq sequencing run to QC library quality before sequencing deeper.

Primer sequences

10

A	B
Primer name	sequence
TS-2-SS-F	TAGATTGATCTAGAGCATGCACCGGACACACTCTTTCCCTACACGAC
TS-2-SS-R	CGAAGCGGCCGGCCGAATTCGTCGAGTGACTGGAGTTCAGACGTGT
208 TruSeq	AATGATACGGCGACCACCGAGATCTACACTCTTTCCCTACACGACGCTCTTCCGATC
i7 PCR	CAAGCAGAAGACGGCATACGA*G*A*T
i7 Index (#1)	CAAGCAGAAGACGGCATACGAGATACATCGGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
i5 Index (\$1)	AATGATACGGCGACCACCGAGATCTACACAGCGCTAGACACTCTTTCCCTACAC*G*A*C
RT_UMI	CAAGCAGAAGACGGCATACGAGATNNNNNNNNBGTGACTGGAGTTCAGACGTGTGCTCTTCCG*A*T*C



Protocol references

Arnold CD, Gerlach D, Stelzer C, Boryń ŁM, Rath M, Stark A. Genome-wide quantitative enhancer activity maps identified by STARR-seq. *Science*. 2013 Mar 1;339(6123):1074-7. doi: 10.1126/science.1232542. Epub 2013 Jan 17. PMID: 23328393. <https://pubmed.ncbi.nlm.nih.gov/23328393/>

Johnson GD, Barrera A, McDowell IC, D'Ippolito AM, Majoros WH, Vockley CM, Wang X, Allen AS, Reddy TE. Human genome-wide measurement of drug-responsive regulatory activity. *Nat Commun*. 2018 Dec 21;9(1):5317. doi: 10.1038/s41467-018-07607-x. PMID: 30575722; PMCID: PMC6303339. <https://pubmed.ncbi.nlm.nih.gov/30575722/>