



Nov 06, 2023

Wholemout Edu Staining (Zebrafish Larvae)

DOI

dx.doi.org/10.17504/protocols.io.n92ld9429g5b/v1

Henna Myllymaki¹, Abigail Elliot¹, Yi Feng¹

¹Centre for Inflammation Research, Institute for Tissue Repair and Regeneration, Edinburgh BioQuarter, University of Edinburgh



Yi Feng

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n92ld9429g5b/v1>

Protocol Citation: Henna Myllymaki, Abigail Elliot, Yi Feng 2023. Wholemout Edu Staining (Zebrafish Larvae). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.n92ld9429g5b/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working for larvae between 1-5 dpf. Our analysis is often focused on the skin tissue, but we saw cells in other part of the fish were also stained.



Created: October 22, 2021

Last Modified: November 06, 2023

Protocol Integer ID: 54435

Keywords: zebrafish embryo, evaluating cell proliferation, tissue cultured cell, assays by invitrogen, cultured cell, zebrafish larvae, cell biology, cell proliferation in vivo, wholemount edu staining, cell proliferation, cell, crucial edu injection step, assay, cancer biology, staining protocol, variable edu incorporation efficiency, inconsistent results due to variable edu incorporation efficiency, developmental biology, field of developmental biology, fast staining speed, vivo, standard kit from invitrogen

Funders Acknowledgements:

Wellcome Trust

Grant ID: 100104/Z/12/Z

Cancer Research UK

Grant ID: C38363/A26931

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](https://www.protocols.io) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](https://www.protocols.io), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

Methods for evaluating cell proliferation are important for research in the field of developmental biology, cancer biology and cell biology. Recent years, Click-iT EdU Assays by Invitrogen is gaining popularity due to its fast staining speed, mild tissue treatment condition, reproducibility and specificity. However, the standard kit from Invitrogen is designed for tissue cultured cells, hence the manufactory provided protocol is optimized for cells cultured on coverslips. EdU staining became a popular method to evaluate cell proliferation in vivo in zebrafish embryos. However, the protocols that we have tested in the past often yield inconsistent results due to variable EdU incorporation efficiency. Here we describe an EdU staining protocol that we optimized for zebrafish larvae, which incorporated a crucial EdU injection step. Our protocol was used successfully by many new students in training, without any issues.

Troubleshooting



EdU staining

4h 50m

- 1 Inject 2nl bolus of 500uM EdU into the yolk of embryos

30m

Note

Reagents for EdU staining can be purchased from -
<https://www.thermofisher.com/uk/en/home/references/protocols/cell-and-tissue-analysis/protocols/click-it-plus-edu-imaging-protocol.html>

embryos were immobilized in 3% methylcellulose in 0.3 x DANEAU's solution (or E3) and returned to induction solution after injection.

- 2 Incubate at **28.5°C** in the dark for **2.5 hours (this can be between 2 -3 hours, but be consistent between experiments; we have done 2 hours)**

2h 30m

- 3 Fix in 4% PFA (**30min**, Room temperature) 🌡️ 25 °C

30m

⌚ 00:30:00

- 4 Wash in 0.1% PBT (**5min**, RT shaking) 🌀 30 rpm, 25°C, 00:05:00 repeat 3 times

15m

- 5 Wash in 3% BSA, Dilute in PBT (**5min**, RT shaking) -grams/100ml, i.e. 3g/100ml, 300mg in 10ml. 🌀 30 rpm, 25°C

5m

- 6 Block in 3% BSA, Dilute in PBT (**1h**, shaking) 🌀 30 rpm, 25°C, 01:00:00

1h

- 7 Incubate in staining solution (as per kit protocol) – we use 250ul per sample (10-15 larvae) **30 min**, 🌀 30 rpm, 25°C, 00:30:00

30m

🧴 250 µL per 10-15 larvae



	1	2	4	5	10	25	50
1X Click-iT® reaction buffer	440 μL	880 μL	1.84 mL	2.25 mL	4.4 mL	10.9 mL	21.9 mL
Copper protectant	10 μL	20 μL	40 μL	50 μL	100 μL	250 μL	500 μL
Alexa Fluor® picolyl azide (Component B)	1.2 μL	2.5 μL	5 μL	6 μL	12.5 μL	31 μL	62 μL
1X Click-iT® EdU buffer additive	50 μL	100 μL	200 μL	250 μL	500 μL	1.25 mL	2.5 mL
Total volume	500 μL	1 mL	2 mL	2.5 mL	5 mL	12.5 mL	25 mL
*Note: Add the ingredients in the order listed in the table.							

Immunostaining in dark (EGFP or hRAS)

1d 17h 50m

- 8 Wash in PBT (**5min**, RT shaking **x3**) 30 rpm, 25°C, 00:05:00 X3 15m
- 9 Wash in 5% Goat Serum, Dilute in PBT (**5min**, RT shaking) 30 rpm, 25°C, 00:05:00 5m
- 10 Block in 5% Goat Serum, Dilute in PBT (**at least 2h**, RT shaking) 30 rpm, 25°C, 02:00:00 2h
- 11 Incubate with anti-GFP Primary Ab 1:200, dilute in 5% Goat Serum in PBT (O/N, +4 °C shaking) 30 rpm, 4°C, 16:00:00 Over night 16h
 250 μL 10-30 larvae

Note

GFP (D5.1)XP Rabbit mAB (Cell Signalling Technology 2956)



12 Wash in PBT (**25-30 min**, RT shaking) **x 8-10 à >4 hours in total!!!**

🌀 30 rpm, 25°C, 00:30:00 X10

5h

13 Incubate with Secondary Ab [] 1:250, Dilute in 3% Goat Serum in PBT (overnight at +4 °C -- preferred) 🌀 30 rpm, 4°C, 16:00:00 Or 25C 2 hours

16h

Note

Alexa Fluor 488 goat anti-rabbit IgG (H+L) (Invitrogen A11008)]
This step over night is preferred, however, 2 hours at RT also works.

14 Wash in PBT (**15 min**, RT shaking) **x4** 🌀 30 rpm, 25°C, 00:15:00 X 4

1h

15 Stain with Hoechst 33258 Dilute 1 in 1000 in PBT (**30mins**, RT shaking)

🌀 30 rpm, 25°C, 00:30:00 in dark

30m

Note

Or other nuclear stain such as DAPI

16 Wash in PBT (15mins, PBT) **x4** 🌀 30 rpm, 25°C, 00:15:00 X 4

1h

Embed and Image

10m

17 Incubate in AF1 solution.

10m

Note

<https://www.cntech-labsupplies.co.uk/products/sample-preparation/film-forming-polymer-solutions>

Sample can be stored in the dark for few weeks

18 Mounting with cover-slip and glass slides.

