


Oct 11, 2023

Version 2

Wnt-3a and R-spo1 conditioned media reporter assay V.2

DOI

<https://dx.doi.org/10.17504/protocols.io.5qpvorz99v4o/v2>



Natalia Martagón¹, Gurdrun Kliem²

¹Charité Universitätsmedizin Berlin; ²Robert Koch Institute



Natalia Martagón

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5qpvorz99v4o/v2>

Protocol Citation: Natalia Martagón, Gurdrun Kliem 2023. Wnt-3a and R-spo1 conditioned media reporter assay. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.5qpvorz99v4o/v2> Version created by [Natalia Martagón](#)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: October 11, 2023

Last Modified: October 11, 2023

Protocol Integer ID: 89164

Keywords: Wnt-3a, conditioned media, activity assay, luciferase reporter assay, R-spondin conditioned media, conditioned media reporter assay protocol, media reporter assay protocol, luciferase reporter assay, reporter hek cell line, conditioned media, rspo1, cell lysi, expressing luciferase

Funders Acknowledgements:

Deutsche Forschungsgemeinschaft (DFG)

Grant ID: 437638965

Abstract

Protocol designed to measure the activity of Wnt-3a or R-spondin-1 (Rspo1) conditioned media.

A reporter HEK cell line expressing luciferase under Wnt-3a stimulation is cultured with conditioned media followed by cell lysis and a luciferase reporter assay. Activity is compared to previous media batches or references.

Image Attribution

Created with BioRender.com



Materials

Assay-reagent solution:

- [M] 20 millimolar (mM) Tricine (MW 179,2)
- [M] 1.07 millimolar (mM) Mg Carbonate x 5 H₂O (MW 485,7)
- [M] 2.67 millimolar (mM) Mg Sulfate x 7 H₂O (MW 246,5)
- [M] 0.1 millimolar (mM) EDTA (MW 372,2)

Preparation (warm up for preparation):

Dissolve:

3.584 g Tricine

0.519 g Mg-CO₃

in ca. 850 mL water

Add 54 µL of the [M] 0.4 Molarity (M) MgSO₄ pH 8 solution

adjust pH to pH 07.8

bring volume to 1 L Store at RT in glass bottles

Lysis solution

luciferase solution stocks preparation

Co Enzym A: .

Sigma C3019-100 mg

Dissolve 100 mg in 1,27 ml in water. Make aliquotes of 100 µL

Luciferin

Dissolve 25 mg D-Luciferin-Na salt in 4,135 ml water. [Stock]=20 mM. Make aliquotes of 500 µL

0.5 M CDTA pH 8:

Dissolve 8,7 g in 50 ml. Buffer the solution first with NaOH-pearls and fine adjustment with liquid NaOH ca. 10 M

0.5 M MgSO₄-Stock, pH 8:

Dissolve 6,15 g in 50 ml in water

1M EDTA-Stock:

Dissolve 29,2 g in 100 ml water

0.5 M Tris-Stocksolution, pH 7,8:

Dissolve 15,14 g Tris-Base(MW 121,14) in 250 ml water. Adjust pH with phosphoric acid



1M DTT-Stock:

Dissolve 1,54 g (MW 154,3) in 10 ml water

Lysis solution stocks

125 mM Tris-Phosphate

10 mM DTT

10 mM CDTA

50% Glycerin

5 % Triton X 100

5x Lysis-Solution preparation

12,5 ml Tris-Phosphat 0,5 M, pH 7,8

500 µL DTT 1 M

1 ml CDTA 0.5 M, pH 8

25 g Glycerin

2,5 g Triton X 100

Bring to a total vol. of 50 ml with water. Make 10 ml aliquotes and store at -20°C

HEK-medium

DMEM F12 HEPES (Gibco Cat. No. 11330),

20% FCS

200 µg/ml G418

material for cell counting

Microplatte for luminescence readers

Greiner Bio-one Cat. no. 655094



Before start

Be sure to have access to:

- Luminometer with or without injection system.
- 96 well plates appropriate for luminescence signal measurements
- Reagents for culturing reporter cell line. See reference (1)

*reagents might take a while to gather and prepare but once prepared they last for a long time

*time to produce conditioned media around 3 weeks. See reference protocols (2), (3).



Day 1: Seeding of Hek 293 STF cells

- 1 Culture HEK 293 STF CRL-3249™ (from now on referred as HEK-STF) according to the company specifications (1) until confluent and not too many passages old.

Seed 3 wells of HEK-STF cells per sample and following controls:

- negative control: conditioned media from L-cells (not transfected with luciferase construct) (4)
- positive control: previous batch with known activity, HEK-STF cells with recombinant Wnt-3a or agonist stimulation
- control lysates as blank for luminescence: HEK-SFT cells with HEK-medium only

- 1.1 Start with one almost confluent T75 culture bottle of HEK-STF cells
- 1.2 1x wash with DPBS: take out medium, add 5 mL DPBS, turn gently the bottle, take out DPBS 
- 1.3 Detach cells with 1-2 mL Trypsin/EDTA ( 37 °C) and transfer to a conical tube with 8-9 mL HEK-medium 5m
- 1.4 Centrifuge at 100-200 g for  00:05:00 5m
count cells (Neubauer chamber or automatic cell counter) 
- 1.5 count cells (Newbauer chamber or automatic cell counter) 5m
- 1.6 Seed cells in a 24 well plate. HEK-STF : 15m
24 well Platte (0.05×10^6 cells/well). 1.3×10^6 cells / 13ml for the whole plate
cover each well with  0.5 mL of HEK medium.

Day 2: cell stimulation

- 2 Aspirate medium and discard. 30m



Add a total of  500 μL of HEK medium with the following amount of conditioned medium (CM) to test:

- For R-spondin CM:  12.5 % volume Wnt-3a CM +  2.5 % volume Rspo1 CM (3)
- For Wnt3a CM .  50 % volume

2.1 Incubate at stadard culture conditions for aprox.  24:00:00

1d



Day3 : luicferasse reaction

20m

3 Luiferase solution

Prepare on the day luciferase solution. Keep in the fridge until ready to use

When using a plate reader with automatic injection, calculate the extra volume needed for it (ca 3 mL)

fill to final volume with **Assay reagent** solution

	A	B	C	D
	Reagent	units	[stock]	[working]
	DTT	M	1	0,033
	CoA	mM	10	0,266
	Luciferin	mM	20	0,467
	ATP	mM	100	0,633

4 Prepare 1x Lysis medium from the stock with destiled water

20m

5 Cell lysis

Aspirate medium from HEK-STF cells

2m

5.1 Add  150 μL /well of 1x Lysis buffer.

2m



- 6 Leave the plate for  00:20:00 at  Room temperature . (Note: pippeting seems to make cell aggregates and bubbles) 20m 
- 7 Check under the microscope that most cells are lysed 
- 8 Shake gently for  00:05:00 on plate shaker 5m
- 9 Transfer  20 μ L of each sample to a 96 well plate for luminescence read. 10m
Avoid bubbles or cell clumps. Recommended pipetting scheme as in the 24 well plate of HEK-STF culture-
- 10 If using a plate reader without automatic injection: 20m
Take to the plate reader with luci. solution, multichanel pipette, and reagent reservoir
Once the plate reader is set, add  100 μ L of luciferse solution per well and read.
For accuracy, use a multichanel pipette and add the solution covering not one sample at a time but one of the triplicates of all samples at one time.
- 11 
(signal decreases ca. 20% in the first 10 minutes)
- 12 Check that the relative luminescence units from negative controls are orders of magnitude lower that test samples. Average values from triplicates.
Values below more than 50-60% lower than a working batch of conditioned media are considered of poor quality.
Examples of equipment and methods: Berthold Tristar-GAS ISRE Luciferase Assay, SpectraMax i3x-SpectraMax Glo Steady-Luc Reporter Assay



Protocol references

(1) HEK 293 STF CRL-3249TM

<https://www.atcc.org/products/crl-3249>

(2) Culture and establishment of self-renewing human and mouse adult liver and pancreas 3D organoids and their genetic manipulation

<https://www.nature.com/articles/nprot.2016.097>

(3) Establishment and Culture of Human Intestinal Organoids Derived from Adult Stem Cells

<https://currentprotocols.onlinelibrary.wiley.com/doi/10.1002/cpim.106>

(4) Notes on the use of conditioned medium

<https://web.stanford.edu/~rnusse/assays/actvwnt.html#:~:text=Comparing%20the%20Wnt%20conditioned%20medium,conditioned%20medium%20may%20therefore%20contain>

*"Comparing the Wnt conditioned medium to medium from control L cells is a reasonable way of controlling for specificity but it should be realized that Wnt expression can lead to the induction of genes encoding other secreted factors (such as FGFs, see **target genes**) and that the conditioned medium may therefore contain such factors."*

image created with BioRender.com