

Oct 25, 2019

## WTc11 iPSC Culture and Maintenance

DOI

[dx.doi.org/10.17504/protocols.io.8rh hv36](https://dx.doi.org/10.17504/protocols.io.8rh hv36)

WTc11 iPSC  
Culture and Maintenance

Connor Ludwig<sup>1</sup>, Ruilin Tian<sup>1</sup>, Martin Kampmann<sup>1</sup>

<sup>1</sup>University of California, San Francisco

Neurodegeneration Met...

KampmannLab



**Martin Kampmann**

University of California, San Francisco

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.8rh hv36](https://dx.doi.org/10.17504/protocols.io.8rh hv36)

**Protocol Citation:** Connor Ludwig, Ruilin Tian, Martin Kampmann 2019. WTc11 iPSC Culture and Maintenance. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.8rh hv36>

**Manuscript citation:**

Tian et al (2019). CRISPR Interference-Based Platform for Multimodal Genetic Screens in Human iPSC-Derived Neurons. Neuron pii: S0896-6273(19)30640-3. [Epub ahead of print] PubMed PMID: 31422865.

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** October 25, 2019

**Last Modified:** October 25, 2019

**Protocol Integer ID:** 29193

**Keywords:** CRISPR interference, CRISPRi, functional genomics, neuron, stem cell, essential genes, CROP-seq, Perturb-Seq, single-cell RNA sequencing, high-content microscopy

## Abstract

This protocol explains general culture and maintenance of the WTc11 iPSC cell line.

## Attachments



WTc11 iPSC CULTURE

S...

17KB

## Materials

### MATERIALS

⊗ StemPro™ Accutase™ Cell Dissociation Reagent **Thermo Fisher Scientific Catalog #A1110501**

⊗ DPBS, no calcium, no magnesium **Thermo Fisher Catalog #14190136**

⊗ Trypan Blue Solution, 0.4% **Thermo Fisher Catalog #15250061**

⊗ StemFlex™ Medium **Thermo Fisher Scientific Catalog #A3349401**

⊗ Matrigel **Corning Catalog #356231**

⊗ KnockOut™ DMEM **Thermo Fisher Scientific Catalog #10829018**



⊗ Rock Inhibitor Y-27632 Dihydrochloride **Tocris Catalog #1254**

## Safety warnings





⚠ Please refer to the Safety Data Sheets (SDS) for safety and health data.

## Before start

Before use, warm complete medium required for that day at room temperature until it is no longer cool to the touch.



Alternatively, an aliquot for use that day may be pre-warmed in a  37 °C waterbath until no longer cool to the touch. Avoid extended dwell times at  37 °C .

## Prepare complete StemFlex Medium





- 1 Thaw the frozen StemFlex Supplement 10X at  Room temperature for ~  02:00:00 or overnight at  2 °C to  8 °C .

### Note

IMPORTANT! Do not thaw the frozen supplement at 37°C.

- 2 Mix the thawed supplement by gently inverting 3–5 times. 
- 3 Aseptically transfer 50 mL of StemFlex Supplement 10X to the bottle of StemFlex™ Basal Medium (450 mL fill).
- 4 Gently invert the bottle several times to obtain 500 mL of homogenous complete medium. 

### Note

Following reconstitution, complete StemFlex™ Medium can be stored at 2°C to 8°C for up to 2 weeks or aliquoted and stored at  -5 °C to  -20 °C for up to 6 months. Alternatively, usage size aliquots of the supplement can be made and frozen at  -5 °C to  -20 °C for up to 6 months. Avoid multiple freeze-thaw cycles.

## Feeding WTc11 iPSC

- 5 Feed the PSCs the day after seeding followed by every-other-day thereafter.

### Note

If the cells are to be left without feeding for two days (for example, over a weekend), then double the feed volume (i.e., 4 mL added per well of 6-well plate).

## Passaging WTc11 iPSC



- 6 iPSCs should be split when cells are ~ 80% confluent.
- 7 Thaw Matrigel on ice and dilute in pre-chilled Knockout DMEM for a final volume of [M] 100 microgram per milliliter ( $\mu\text{g}/\text{mL}$ ) .
- 8 Coat desired wells/plates with diluted Matrigel and incubate at  $37^\circ\text{C}$  for 00:30:00 - 01:00:00 using the following table for volumes to add **per well**:

Per:	96-well	24-well	12-well	6-well	10-cm dish	15-cm dish
Volume to add:	40 $\mu\text{L}$	200 $\mu\text{L}$	0.5 mL	1 mL	5 mL	10 mL

#### Note

Matrigel may be re-used during this time only to coat additional plates. Original plates should have PBS or media to prevent the matrix from drying out. Matrigel coated plates must be used within 14 days of coating.


- 9 Tilt cell-containing plate towards you and aspirate existing media.
- 10 Wash wells once with ample PBS (about 2x amount of media).
- 11 Add accutase to well(s) using the following table for volumes per well and incubate at  $37^\circ\text{C}$  for 00:03:00 ; add another 01:00:00 - 00:02:00 if cells have not mostly lifted/dissociated.

Per:	96-well	24-well	12-well	6-well	10-cm dish	15-cm dish
Volume to add:	20 $\mu\text{L}$	100 $\mu\text{L}$	250 $\mu\text{L}$	0.5 mL	2 mL	4 mL




- 12 Add ample PBS to accutase-containing well(s) to dilute accutase using the following table for volumes per well:

Per:	96-well	24-well	12-well	6-well	10-cm dish	15-cm dish
Volume to add:	200 $\mu$ L	1 mL	2.5 mL	5 mL	10 mL	20 mL

- 13 Pipette up and down gently to mechanically release remaining cells, collect, and add to appropriately-sized conical tubes. 

- 14 Spin cells at  200 x g for  00:05:00 at  Room temperature . 


- 15 Carefully aspirate supernatant from pelleted conicals.


- 16 Add appropriate volume of StemFlex + Rock inhibitor at  10 micromolar ( $\mu$ M) to conicals according to pellet size for counting.

#### Note

Rock inhibitor should only be used when cells are individualized or in small colonies (typically the first two days after passaging); the presence of Ri at higher densities results in cell stress/death, and in general, Rock inhibitor greatly reduces proliferation.

#### Note

For first time use of Rock inhibitor, it is suggested to aliquot at 10mM [1000x], diluting in DPBS, and use on cells at concentration  10 Mass Percent .

- 17 Triturate to resuspend cells in StemFlex + Rock inhibitor and remove 10 $\mu$ L and add this volume to a 1.5mL Eppendorf tube. 

#### Note

Be careful to minimize contact of pipette with the side of the conical wall.



- 18 Count cells and calculate desired number of cells to seed, and dilute this volume with additional StemFlex + Rock inhibitor to plate using the following table for volume to add per well:

<b>Per:</b>	96-well	24-well	12-well	6-well	10-cm dish	15-cm dish
<b>Volume to add:</b>	50-100 $\mu$ L	250-500 $\mu$ L	0.75-1 mL	1.5-2 mL	8-12 mL	15-25 mL

#### Note

For general passaging where exact cell number seeded is not important, adding resuspended cells at 1:100, 1:50, and 1:20 the final well volume typically provides near-confluency in 7 days, 5 days, and 3 days, respectively.

- 19 iPSCs can be frozen in StemFlex + 10% DMSO.

