



Oct 03, 2024

X-region PolyU/UC mutagenesis

DOI

<https://dx.doi.org/10.17504/protocols.io.8epv5rxy6g1b/v1>

Carolina Lopez¹

¹Washington University



Cecilia Escudero

wustl

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8epv5rxy6g1b/v1>

Protocol Citation: Carolina Lopez 2024. X-region PolyU/UC mutagenesis. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.8epv5rxy6g1b/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 28, 2024

Last Modified: October 03, 2024

Protocol Integer ID: 100753

Keywords: uc mutagenesis protocol to pcr, uc mutagenesis protocol, pcr

Abstract

Protocol to PCR and in vitro transcribe the HCV X region w/out mutagenesis

Materials

STAR METHOD

REAGENT	SOURCE	IDENTIFIER
Platinum® <i>Taq</i> DNA Polymerase kit	Invitrogen	Cat#10966-018
Set of dATP, dCTP, dGTP, dTTP	Promega	Cat# U1420
MEGAscript® T7 Transcription Kit	Ambion Thermo Fisher Technologies	Cat#1334
HCV polyU/UC plasmid	Gift from Michael Gale lab	NA
HCV X-region plasmid	Gift from Michael Gale lab	NA
UltraPure Water	Invitrogen	Cat#10977-015

Material:

Primer: X-region F/ X-region R (10 micromolar (μ M))

PolyU/UC F and R (10 micromolar (μ M))

dNTPs (10 millimolar (mM)) Fermentas –For PCR

HCV PolyU/UC X-region plasmid



HCV X-region PolyU/UC PCR - T7 in vitro transcription

1 PCR:

- 1) Thaw component and cDNA, spin down before use
- 2) Assemble reaction: Master Mix: 10X (each reaction assemble 3 tubes)

Reagent	Volume (ul)
H2O	26
10X Buffer	5
MgCl ₂	5
dNTP (10uM)	2
Forward Primer (10uM)	5
Reverse Primer (10uM)	5
Taq pol. (5unit/ul)	0.25
Plasmid (10ng/ul)	1.5

Total:  50 µL

3) Run PCR

Program	Temperature (C)	Time	
Denaturation	95	3 min	Hold
Denaturation	95	30 sec	35 cycles
Annealing	55	30 sec	
Extension	72	90 sec	
Final extension	72	5 min	Hold
	4	∞	Hold

II. In vitro transcription:

15m

- 2
 - Spin-down all reagents
 - Keep the 10X Reaction Buffer at room temperature while assembling the reaction.
 - Vortex all solutions until they are completely in solution.
 - Assemble dNTP mix (2µL ATP+ 2µL CTP+ 2µL GTP+ 2µL TTP/reaction).
 - Store the ribonucleotides on ice!

15m



Assemble the reaction in an appropriate RNase free PCR tube.

- 1)  1 μg DNA plasmid (digested)
- 2)  8 μL NTP Mix
- 3)  2 μL 10X Reaction Buffer
- 4)  0.5-1 μL RNaseOUT (RNase Inhibitor)
- 5)  2 μL T7 enzyme mix
- 6) RNase Free water up to  20 μL

- Pipette the mixture up and down gently and microfuge
- Incubate 3-6hrs at  37 °C using PCR machine. The shorter the ivtRNA, the longer the reaction time needed.
- (Optional but do most times) Template plasmid Digestion: 1 μL Turbo DNase and mix well,  00:15:00 at  37 °C

IV. RNA precipitation: Transfer to an 1.5ml RNase free vial

35m

- 3 If you expect very low yields of RNA, do not dilute the transcription reaction with water prior to adding the LiCl.

35m

- Add  30 μL H₂O
 30 μL LiCl
- Mix well, incubate  00:30:00 at  -20 °C
- Centrifuge at 4°C/20 min/max
- Remove supernatant
- Add  1 mL 70% COLD ethanol
- Centrifuge at 4°C/5 min/max
- Remove supernatant and allow the pellet to dry
- Resuspend in  15-20 μL RNase-free H₂O
- Place the tube  00:05:00 at  65 °C water bath
- Dilute  1 μL in  19 μL H₂O to measure concentration

Protocol references

Jie Xu (Lopez Lab)