

Feb 18, 2025

XChem crystallographic fragment screening

 Forked from [XChem solvent test](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.n2bvjnzwxgk5/v1>



Peter Marples^{1,2,3}, Daren Fearon^{1,2,3}, Blake Balcomb^{1,2,3}

¹Diamond Light Source; ²Research Complex at Harwell; ³ASAP Discovery Consortium

ASAP Discovery

OpenBind Consortium



Peter Marples

Diamond Light Source

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n2bvjnzwxgk5/v1>

External link: <https://www.jove.com/v/62414/achieving-efficient-fragment-screening-at-xchem-facility-at-diamond>

Protocol Citation: Peter Marples, Daren Fearon, Blake Balcomb 2025. XChem crystallographic fragment screening. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.n2bvjnzwxgk5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 26, 2024

Last Modified: February 18, 2025

Protocol Integer ID: 104139

Keywords: XChem, FBDD, fragment, drug discovery, fragment screening, xchem crystallographic fragment, crystallographic fragment screening protocol, large scale crystal graphic fragment screening, xchem facility at diamond light source, crystallographic fragment, essential part to the fbdd toolkit, fbdd toolkit, large scale crystal, xchem facility, diamond light source, entire crystal to deposition process, entire crystal, graphic fragment screening, xchem, structural information about the interaction, structural information, weak binder, molecular level

Funders Acknowledgements:

National Institutes of Health/National Institute Of Allergy and Infectious Diseases (NIH/NIAID)

Grant ID: U19AI171399

Disclaimer

The content is solely the responsibility of the authors and does not necessarily represent the official views of the National Institutes of Health.

Acknowledgements:

Diamond Light Source Ltd, Harwell Science and Innovation Campus, Didcot OX11 0QX, UK
Research Complex at Harwell, Harwell Science and Innovation Campus, Didcot OX11 0FA, UK

Abstract

The XChem Facility at Diamond Light Source provides routine, large scale crystal graphic fragment screening, supporting the entire crystal to deposition process. X-ray crystallography is an essential part to the FBDD toolkit. It is sensitive enough to identify weak binders and directly yields structural information about the interactions on a molecular level. Here we summarise a generic layout of a crystallographic fragment screening protocol.

Guidelines

N/A



Materials

SwissCI 3 lens crystallization plates <https://swissci.com/product/3-lens-crystallisation-plate/> Codes: Midi: UVXPO-3LENS 3W96T-PS 3W96T-UVP

Uni Pucks <https://moleculardimensions.com/en/product/MD7-601> Codes: Catalog # MD7-601

Uni Wand <https://moleculardimensions.com/en/product/MD7-411> Codes: Catalog # MD7-411

Loops <https://moleculardimensions.com/en/category/mounted-litholoops> Codes: Catalog # MD7-131, Catalog # MD7-133, Catalog # MD7-135, Catalog # MD7-138

DMSO <https://www.sigmaaldrich.com/GB/en/substance/781367685> Codes: **D5879**

Glycerol <https://www.sigmaaldrich.com/GB/en/substance/920956815> Codes: **G5516**

Ethylene Glycol <https://www.sigmaaldrich.com/GB/en/substance/6207107211> Codes: **324558**

Safety warnings

 Follow all handling warning for the chemicals use.

Use of cryogen liquid nitrogen; cold burns can occur with exposure of skin to very low temperatures, contact with eyes can cause permanent damage, evaporating liquid nitrogen can cause oxygen depletion, asphyxia and death. Wear appropriate footwear and PPE when handling including cryo gloves and safety glasses. Use appropriate volumes of liquid nitrogen for the ventilation of the laboratory while in use.

Ensure centrifuge is properly balanced when used.

Before start

Grow any necessary crystal plates in sitting drop plates.



Imaging

- 1 Image all sitting drop crystal plates required for the experiment in a crystal plate imager.



Selecting crystals

- 2 Based on the images, select drops with suitable crystals for the experiment. Fragment dispensing is done in singletons. i.e. one compound per drop of crystals. Ensure enough drops for the total number of compounds required for the experiment.
- 3 If you have a targeting system, target the drops for dispensing compound, not directly targeting the crystals directly to avoid damaging crystals.

Fragment dispensing

3h 2m

- 4 Based on the results of your solvent test, calculate the % v/v required of compound to dispense into crystal drops. This percentage chosen will now include compound to dispense into crystal drops.

Note

When changing from solvent to compound, the tolerance of crystals often is reduced, so the optimal soaking time may need to be lowered by 5% v/v compared to the solvent test results.

Protocol



NAME

XChem solvent test

CREATED BY

Peter Marples[Preview](#)



- 5 Take the compound plate containing the fragments and spin the plate in a benchtop centrifuge for 00:02:00 at 200 x g .
- 6 Dispense the required volume of compound into the corresponding crystal drop.
- 7 Store the plate in the incubator for 01:00:00 to 03:00:00 depending on the results of your solvent test at either 4 °C or 20 °C depending on the crystal growth temperature.

2m



3h



Harvesting crystals

- 8 In preparation for crystal harvesting, if cryoprotectant is required add the necessary volume of cryo onto crystal drops just before full soaking time has ended.
- 9 Select a tray of the appropriately sized loops for the crystals (35 µm , 75 µm , 150 µm or 300 µm). Importantly, choose a loop size that matches the size of the crystal as closely as possible to improve the data quality.
- 10 Place the pucks in a foam dewar and cool them down with liquid nitrogen.

Safety information
Wear appropriate PPE when handling liquid nitrogen.
- 11 If the crystal has survived, mount the crystal in the loop and plunge into the liquid nitrogen, placing the sample into a puck.
- 12 If the crystal has been successfully mounted, record the success and continue to the next sample. Otherwise record the failure and continue to the next sample until all samples that can be harvested have been mounted.
- 13 Place lids on the pucks and store.
- 14 Send the pucks to be collected.





Protocol references

Douangamath, A., Powell, A., Fearon, D., Collins, P. M., Talon, R., Krojer, T., Skyner, R., Brandao-Neto, J., Dunnett, L., Dias, A., Aimon, A., Pearce, N. M., Wild, C., Gorrie-Stone, T., von Delft, F. Achieving Efficient Fragment Screening at XChem Facility at Diamond Light Source. *J. Vis. Exp.* (171), e62414, doi:10.3791/62414 (2021)

<https://www.jove.com/v/62414/achieving-efficient-fragment-screening-at-xchem-facility-at-diamond>