


Part of **SPRINGER NATURE**

Feb 18, 2025

DOI

<https://dx.doi.org/10.17504/protocols.io.5jyl82518l2w/v1>



Peter Marples^{1,2,3}, Daren Fearon^{1,2,3}, Blake Balcomb^{1,2,3}

¹Diamond Light Source; ²Research Complex at Harwell; ³ASAP Discovery Consortium

ASAP Discovery

OpenBind Consortium



Peter Marples

Diamond Light Source

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5jyl82518l2w/v1>

External link: <https://www.jove.com/v/62414/achieving-efficient-fragment-screening-at-xchem-facility-at-diamond>

Protocol Citation: Peter Marples, Daren Fearon, Blake Balcomb 2025. XChem solvent test. [protocols.io](https://www.protocols.io)

<https://dx.doi.org/10.17504/protocols.io.5jyl82518l2w/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 25, 2024

Last Modified: February 18, 2025

Protocol Integer ID: 98776

Keywords: XChem, solvent test, FBDD, fragment, drug discovery, solvent test the xchem facility, large scale crystal graphic fragment screening, soaking conditions for fragment screening, generic solvent test protocol, solvent test, fragment screening, xchem facility, solvent, entire crystal to deposition process, xchem, graphic fragment screening, large scale crystal, identifying soaking condition, fbdd toolkit, essential part to the fbdd toolkit, diamond light source, entire crystal, structural information about the interaction, molecular level, screening

Funders Acknowledgements:

National Institutes of Health/National Institute Of Allergy and Infectious Diseases (NIH/NIAID)

Grant ID: U19AI171399

Disclaimer

The content is solely the responsibility of the authors and does not necessarily represent the official views of the National Institutes of Health.

Acknowledgements:

Diamond Light Source Ltd, Harwell Science and Innovation Campus, Didcot OX11 0QX, UK
Research Complex at Harwell, Harwell Science and Innovation Campus, Didcot OX11 0FA, UK

Abstract

The XChem Facility at Diamond Light Source provides routine, large scale crystal graphic fragment screening, supporting the entire crystal to deposition process. X-ray crystallography is an essential part to the FBDD toolkit. It is sensitive enough to identify weak binders and directly yields structural information about the interactions on a molecular level. Here we summarise a generic solvent test protocol necessary for identifying soaking conditions for fragment screening.

Guidelines

N/A



Materials

SwissCI 3 lens crystallization plates <https://swissci.com/product/3-lens-crystallisation-plate/> Codes: Midi: UVXPO-3LENS 3W96T-PS 3W96T-UVP

Uni Pucks <https://moleculardimensions.com/en/product/MD7-601> Codes: Catalog # MD7-601

Uni Wand <https://moleculardimensions.com/en/product/MD7-411> Codes: Catalog # MD7-411

Loops <https://moleculardimensions.com/en/category/mounted-litholoops> Codes: Catalog # MD7-131, Catalog # MD7-133, Catalog # MD7-135, Catalog # MD7-138

DMSO <https://www.sigmaaldrich.com/GB/en/substance/781367685> Codes: **D5879**

Glycerol <https://www.sigmaaldrich.com/GB/en/substance/920956815> Codes: **G5516**

Ethylene Glycol <https://www.sigmaaldrich.com/GB/en/substance/6207107211> Codes: **324558**

Safety warnings

 Follow all handling warning for the chemicals use.

Use of cryogen liquid nitrogen; cold burns can occur with exposure of skin to very low temperatures, contact with eyes can cause permanent damage, evaporating liquid nitrogen can cause oxygen depletion, asphyxia and death. Wear appropriate footwear and PPE when handling including cryo gloves and safety glasses. Use appropriate volumes of liquid nitrogen for the ventilation of the laboratory while in use..

Before start

Grow any necessary crystal plates in sitting drop plates.



Imaging

- 1 Image all sitting drop crystal plates required for the experiment in a crystal plate imager.



Selecting crystals

- 2 Based on the images, select drops with suitable crystals for the experiment. In a standard solvent test, 38 drops are required; 2 control drops and 12 drops for varying solvent concentrations at three time points
- 3 If you have a targeting system, target the drops for dispensing solvent, not directly targeting the crystals to avoid damaging crystals.

Solvent dispensing

3h

- 4 Use the drop volume to calculate the required solvent needed to have drops with a final % of solvent equal to 5, 10, 15, 20, 25 and 30 % v/v. For control crystals, leave two drops without the addition of solvent. i.e. 0% v/v.
- 5 Dispense the 5, 10, 15, 20, 25 and 30% of solvent into individual drops in duplicate.
- 6 Store the plate in the incubator for 01:00:00 at either 4 °C or 20 °C depending on the crystal growth temperature.
- 7 Repeat steps 5 & 6 [go to step #5](#) two times, resulting in a 01:00:00 , 02:00:00 and 03:00:00 solvent soaking duration for the entire solvent test.

1h



2h

Harvesting crystals

- 8 In preparation for crystal harvesting, if cryoprotectant is required add the necessary volume of cryo onto crystal drops just before full soaking time has ended.
- 9 Select a tray of the appropriately sized loops for the crystals (for example 35 μm , 75 μm , 150 μm or 300 μm). Importantly, choose a loop size that matches the size of the crystal as closely as possible to improve the data quality.



- 10 Place the pucks in a foam dewar and cool them down with liquid nitrogen.

Safety information

Wear appropriate PPE when handling liquid nitrogen.

- 11 If the crystal has survived, mount the crystal in the loop and plunge into the liquid nitrogen, placing the sample into a puck.
- 12 If the crystal has been successfully mounted, record the success and continue to the next sample. Otherwise record the failure and continue to the next sample until all samples that can be harvested have been mounted.
- 13 Place lids on the pucks and store.
- 14 Send the pucks to be collected. Based on the results, select the highest % v/v of solvent for the longest soaking time that results in viable data for your needs to decide your solvent % in future soaking experiments.

Protocol references

Douangamath, A., Powell, A., Fearon, D., Collins, P. M., Talon, R., Krojer, T., Skyner, R., Brandao-Neto, J., Dunnett, L., Dias, A., Aimon, A., Pearce, N. M., Wild, C., Gorrie-Stone, T., von Delft, F. Achieving Efficient Fragment Screening at XChem Facility at Diamond Light Source. *J. Vis. Exp.* (171), e62414, doi:10.3791/62414 (2021)

<https://www.jove.com/v/62414/achieving-efficient-fragment-screening-at-xchem-facility-at-diamond>