

Oct 13, 2019

Xylem inoculation with flg22

DOI

dx.doi.org/10.17504/protocols.io.752hq8e

Cleo Bagchus¹

¹Wageningen University

iGEM Wageningen 2019



Cleo B.

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account





DOI: https://dx.doi.org/10.17504/protocols.io.752hq8e

Protocol Citation: Cleo Bagchus 2019. Xylem inoculation with flg22 . protocols.io

https://dx.doi.org/10.17504/protocols.io.752hq8e

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol



Created: October 13, 2019

Last Modified: October 13, 2019

Protocol Integer ID: 28570

Keywords: plant infection with xanthomonas campestris campestri, xylem inoculation, xanthomonas campestris campestri,

plant infection

Abstract

Protocol can be used after: Plant Infection with Xanthomonas Campestris Campestris

Guidelines

Depending on the country, X. campestris may be an quarantine organism. Please make sure you follow the specific guidelines of your country.

Materials

sthree-week-old Brassica oleracea plants in climate chamber Needle 1 uM flg22 at room temperature

Troubleshooting



Day 1

- 1 Place \perp 15 μ L of flg22 3 cm above the lesion on the central vein of the leaf.
- 2 With the needle puncture through the droplet into the central vein three times. Do not puncture the leaf completely. Due to a difference in pressure, the droplet should be sucked into the xylem.
- 3 Place a second \perp 15 μ L of flg22 3.5 cm above the lesion on the central vein of the leaf.
- 4 Repeat step 2.