


May 17, 2019

## Yale - Alkaline Phosphatase

DOI

[dx.doi.org/10.17504/protocols.io.yz7fx9n](https://dx.doi.org/10.17504/protocols.io.yz7fx9n)



Gary Cline<sup>1</sup>, John Stack<sup>1</sup>

<sup>1</sup>Yale University

Mouse Metabolic Phenotyping Centers  
Tech. support email: [info@mmpc.org](mailto:info@mmpc.org)



Lili Liang

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.yz7fx9n](https://dx.doi.org/10.17504/protocols.io.yz7fx9n)

External link: <https://mmpc.org/shared/document.aspx?id=208&docType=Protocol>

**Protocol Citation:** Gary Cline, John Stack 2019. Yale - Alkaline Phosphatase. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.yz7fx9n>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** March 12, 2019

**Last Modified:** May 17, 2019

**Protocol Integer ID:** 21279

**Keywords:** Alkaline Phosphatase activity



## Abstract

### Summary:

Procedure to measure the amount of Alkaline Phosphatase activity. Alkaline Phosphatase (ALP) activity is measured from the hydrolysis of 4-nitrophenylphosphate to 4-nitrophenoxide ion (monitored at 405 nm) and phosphate.

## Materials

### MATERIALS

⊗ Alkaline Phosphatase Reagent **Prolabs(cliniqa) Catalog #R85120**

⊗ Assayed Control Serum 1 **Prolabs(cliniqa) Catalog #R83082**

⊗ Assayed Control Serum 2 **Prolabs(cliniqa) Catalog #R83083**

### Reagent Preparation:

***Alkaline Phosphatase Reagent:*** Add the appropriate amount of water (6.5mL) to the reagent bottle. Invert to mix, allowing 15 minutes for the reagent to settle.

***Assayed Control Serum 1:*** Add the appropriate amount of water (6.5mL) to the chemical control bottle. Invert to mix, allowing 15 minutes for the reagent to settle.

***Assayed Control Serum 2:*** Add the appropriate amount of water (6.5mL) to the chemical control bottle. Invert to mix, allowing 15 minutes for the reagent to settle.

## Before start

Analysis by automated system Cobas Mira Plus.



- 1 Calibrate Cobas for Alkaline Phosphatase Activity analysis by running two assayed control serum.
- 2 Sample handling as performed by the Cobas Mira Plus.
  - a) Pipette 3  $\mu\text{L}$  of sample into a cuvette slot.
  - b) Add 150  $\mu\text{L}$  of Alkaline Phosphatase Reagent.
  - c) Mixture is incubated at 37°C and spun for 10 minutes.
  - d) Absorbance is measured at 405 nm.