



May 16, 2019

Yale - Aspartate Amino Transferase

DOI

dx.doi.org/10.17504/protocols.io.yz9fx96



Gary Cline¹, John Stack¹

¹Yale University

Mouse Metabolic Phenotyping Centers
Tech. support email: info@mmpc.org



Lili Liang

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yz9fx96>

External link: <https://mmpc.org/shared/document.aspx?id=207&docType=Protocol>

Protocol Citation: Gary Cline, John Stack 2019. Yale - Aspartate Amino Transferase. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.yz9fx96>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 11, 2019

Last Modified: May 16, 2019

Protocol Integer ID: 21281

Keywords: Aspartate Amino Transferase (AST), activity of aspartate amino transferase, aspartate amino transferase summary, aspartate amino transferase, reactions of ast, oxaloacetate from aspartate, ast activity, oxidation of nadh, ast, aspartate, rate of nad formation, dehydrogenase, conversion of oxaloacetate, nad formation, change in absorbance, nad, oxaloacetate, absorbance


Abstract

Summary:

Procedure used to measure the activity of Aspartate Amino Transferase (AST). AST activity is measured by the enzymatically coupled reactions of AST (to form oxaloacetate from aspartate and α -ketoglutarate) and malate dehydrogenase (conversion of oxaloacetate to malate with oxidation of NADH to NAD). The rate of NAD formation is monitored by the change in absorbance at 340 nm.

Materials

MATERIALS

 Aspartate Amino Transferase Reagent **Prolabs(cliniqa) Catalog #R85121**

 Assayed Control Serum 1 **Prolabs(cliniqa) Catalog #R83082**

 Assayed Control Serum 2 **Prolabs(cliniqa) Catalog #R83083**

Reagent Preparation:

Aspartate Amino Transferase Reagent: Add the appropriate amount of water (6.5mL) to the reagent bottle. Invert to mix, allowing 15 minutes for the reagent to settle.

Assayed Control Serum 1: Add the appropriate amount of water (6.5mL) to the chemical control bottle. Invert to mix, allowing 15 minutes for the reagent to settle.

Assayed Control Serum 2: Add the appropriate amount of water (6.5mL) to the chemical control bottle. Invert to mix, allowing 15 minutes for the reagent to settle.

Troubleshooting

Before start

Analysis by automated system Cobas Mira Plus.



- 1 Calibrate Cobas for Aspartate Amine Transferase Activity analysis by running twoassayed control serum.
- 2 Sample handling as performed by the Cobas Mira Plus.
 - a) Pipette 16 μ L of sample into a cuvette slot.
 - b) Add 145 μ L of Aspartate Amino Transferase Reagent.
 - c) Mixture is incubated at 37°C and spun for 10 minutes.
 - d) Absorbance is measured at 340 nm.