

May 17, 2019

Yale - Blood Glucose

DOI

dx.doi.org/10.17504/protocols.io.y2kfycw



John Stack¹, Gary Cline¹

¹Yale University

Mouse Metabolic Phenotyping Centers
Tech. support email: info@mmpc.org



Lili Liang

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.y2kfycw

External link: <https://mmpc.org/shared/document.aspx?id=202&docType=Protocol>

Protocol Citation: John Stack, Gary Cline 2019. Yale - Blood Glucose. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.y2kfycw>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 11, 2019

Last Modified: May 17, 2019

Protocol Integer ID: 21292

Keywords: Blood Glucose



Abstract

Summary:

Procedure used to measure glucose concentrations. Glucose is measured by the enzymatically coupled reactions of hexokinase and glucose-6-P dehydrogenase. The rate of NADH formation is monitored by the change in absorbance at 340 nm.

Materials

MATERIALS

⊗ Glucose Reagent 1 **Prolabs(cliniqa) Catalog #R84682**

⊗ Glucose Reagent 2 **Prolabs(cliniqa) Catalog #R84682**

⊗ Multi Analyte Calibrator **Prolabs(cliniqa) Catalog #R60010**

⊗ Assayed Control Serum 1 **Prolabs(cliniqa) Catalog #R83082**

⊗ Assayed Control Serum 2 **Prolabs(cliniqa) Catalog #R83083**

Reagent Preparation:

Glucose Reagent 1: As supplied by Vendor

Glucose Reagent 2: As supplied by Vendor

Multi Analyte Calibrator: Add the appropriate amount of water (6.5mL) to the chemical control bottle. Invert to mix, allowing 15 minutes for the reagent to settle.

Assayed Control Serum 1: Add the appropriate amount of water (6.5mL) to the chemical control bottle. Invert to mix, allowing 15 minutes for the reagent to settle.

Assayed Control Serum 2: Add the appropriate amount of water (6.5mL) to the chemical control bottle. Invert to mix, allowing 15 minutes for the reagent to settle.

Note: Author switched vendors and are now using a **Sekisui 235-60 reagent** for glucose assay. (There is now only 1 reagent instead of 2).

Before start

Analysis by automated system Cobas Mira Plus



- 1 Calibrate Cobas for Glucose analysis by running a multi analyte standard and two control serum.
- 2 Sample Handling as performed by the Cobas Mira Plus.
 - a) Pipette 3 μL of sample into cuvette.
 - b) Absorbance is measured at 340 nm.
 - c) Add 100 μL of Glucose liquid reagent.
 - d) Mixture is incubated at 37°C for 10 minutes.
 - e) Absorbance is measured at 340 nm. Change in absorbance is calculated.