



May 17, 2019

Yale - Blood Urea Nitrogen

DOI

dx.doi.org/10.17504/protocols.io.y3qfymw



John Stack¹, Gary Cline¹

¹Yale University

Mouse Metabolic Phenotyping Centers
Tech. support email: info@mmpc.org



Lili Liang

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.y3qfymw>

External link: <https://mmpc.org/shared/document.aspx?id=203&docType=Protocol>

Protocol Citation: John Stack, Gary Cline 2019. Yale - Blood Urea Nitrogen. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.y3qfymw>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 12, 2019

Last Modified: May 17, 2019

Protocol Integer ID: 21328

Keywords: Blood Urea Nitrogen, blood urea nitrogen summary, conversion of ammonia, ammonia, urea, reactions of urease, rate of nad formation, glutamate dehydrogenase, oxidation of nadh, serum, blood

Abstract

Summary:

Procedure used to measure the concentration of Blood Urea Nitrogen(BUN) in blood, plasma, and serum. Urea is determined by the enzymatically coupled reactions of urease (to form ammonia) and glutamate dehydrogenase (conversion of ammonia and glutamate to glutamine with oxidation of NADH to NAD). The rate of NAD formation is monitored by the change in absorbance at 340 nm.

Materials

MATERIALS

⊗ BUN liquid Reagent **Prolabs(cliniqa) Catalog #R84533**

⊗ Multi Analyte Calibrator **Prolabs(cliniqa) Catalog #R60010**

⊗ Assayed Control Serum 1 **Prolabs(cliniqa) Catalog #R83082**

⊗ Assayed Control Serum 2 **Prolabs(cliniqa) Catalog #R83083**

Reagent Preparation:

BUN liquid Reagent: As supplied by vendor

Multi Analyte Calibrator: Add the appropriate amount of water (6.5mL) to the chemical control bottle. Invert to mix, allowing 15 minutes for the reagent to settle.

Assayed Control Serum 1: Add the appropriate amount of water (6.5mL) to the chemical control bottle. Invert to mix, allowing 15 minutes for the reagent to settle.

Assayed Control Serum 2: Add the appropriate amount of water (6.5mL) to the chemical control bottle. Invert to mix, allowing 15 minutes for the reagent to settle.

Troubleshooting

Before start

Analysis by automated system Cobas Mira Plus

- 1 Calibrate Cobas for BUN analysis by running a multi analyte standard and two assayed control serums.

- 2 Sample handling as performed by the Cobas Mira Plus.
 - a) Cobas pipettes 2 μ L of sample into a cuvette slot.
 - b) Absorbance is measured at 340 nm.
 - c) Add 200 μ L of BUN liquid reagent.
 - d) Mixture is incubated at 37°C for 10 minutes.
 - e) Absorbance is measured at 340 nm. Change in absorbance is calculated.