



May 17, 2019

## Yale - Magnesium

DOI

[dx.doi.org/10.17504/protocols.io.y3zfyp6](https://dx.doi.org/10.17504/protocols.io.y3zfyp6)



John Stack<sup>1</sup>, Gary Cline<sup>1</sup>

<sup>1</sup>Yale University

Mouse Metabolic Phenotyping Centers  
Tech. support email: [info@mmpc.org](mailto:info@mmpc.org)



Lili Liang

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.y3zfyp6>

External link: <https://mmpc.org/shared/document.aspx?id=218&docType=Protocol>

**Protocol Citation:** John Stack, Gary Cline 2019. Yale - Magnesium. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.y3zfyp6>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** March 12, 2019

**Last Modified:** May 17, 2019

**Protocol Integer ID:** 21337

**Keywords:** Magnesium, magnesium summary, concentration of magnesium, magnesium, chlorophosphonazo iii complex, calcium, calcium interference, concentration, addition of egta

## Abstract

### Summary:

Procedure used to determine the concentration of magnesium in blood, serum, and plasma. Magnesium is measured as the chlorophosphonazo III complex (addition of EGTA prevents calcium interference) and monitored at 600 nm.

## Materials

### MATERIALS

 Magnesium Reagent **Prolabs(cliniqa) Catalog #125-12**

 Magnesium Calibrator **Prolabs(cliniqa) Catalog #125-12**

### Reagent Preparation:

**Magnesium Reagent:** As supplied by vendor

**Magnesium Calibrator:** As supplied by vendor

## Troubleshooting

## Before start

*Analysis by automated system Cobas Mira Plus*

- 1 Calibrate Cobas for Magnesium analysis by running the Magnesium Calibrator and a Water Blank.
- 2 Sample handling as performed by the Cobas Mira Plus
  - a) Pipette 4  $\mu\text{L}$  of sample into a cuvette slot.
  - b) Add 350  $\mu\text{L}$  of Magnesium Reagent.
  - c) Mixture is incubated at 37°C for 10 minutes
  - d) Absorbance is measured at 600 nm.