


Apr 13, 2020

## Yeast growth profile analysis

DOI

[dx.doi.org/10.17504/protocols.io.33jgqkn](https://dx.doi.org/10.17504/protocols.io.33jgqkn)



Monica Rieth<sup>1</sup>

<sup>1</sup>Southern Illinois University-Edwardsville

Labyrieth

Yeast Protocols, Tools, a...



Monica Rieth

Southern Illinois University-Edwardsville

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.33jgqkn>

**Document Citation:** Monica Rieth 2020. Yeast growth profile analysis. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.33jgqkn>



**License:** This is an open access document distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Created:** June 11, 2019

**Last Modified:** April 13, 2020

**Document Integer ID:** 24395

**Keywords:** yeast, growth, profile, growth curve of baker, assessment of growth media, growth curve, saccharomyce, growth media, incorporation of heterologous gene, heterologous gene, addition of metabolite, strain

## Abstract

This method outlines a procedure for analyzing the growth curve of Baker's yeast, *Saccharomyces cerevisiae*, for the assessment of growth medias, the addition of metabolites to media, and the incorporation of heterologous genes to our yeast strain.

## Troubleshooting



## PROCEDURE

### Day 1:

1. Streak a YPD-agar plate from a glycerol stock.
2. Incubate for 48-72 hours at 25-30 °C or until colonies are large enough to pick (1-2 mm diam.)

### Day 2:

3. Check plates. There likely won't be any colonies that are large enough to pick

### Day 3:

4. Colonies should be large enough to pick (1-2 mm diameter)
5. Propagate colony (ies) in 5 mL YPD media
6. Allow colonies to grow overnight at 30 °C shaking at 240-250 rpm

### Day 4:

7. Propagate triplicates in 5 mL fresh YPD using 50 µL of the overnight growth added directly to the 5 mL of fresh media
8. Begin collecting data for the growth curve. This marks the t=0 time point.
9. To begin collecting time points, 50 µL of culture can be added to 750 µL sterile media (or ddH<sub>2</sub>O) using a polystyrene semi-microcuvette. Measure the absorbance at 600 nm  
NOTE: Yeast grow more slowly than bacteria and time points can safely be collected every 2 hours. Additionally, after the fresh media is inoculated and time point 0 is collected, it can take up to 4 hours for any appreciable absorbance to be reached.
10. Time points data can be collected up to and including 48 hours. After 12 hours data was collected every 4 hours up to 24 hours and every 12 hours thereafter.
11. After an OD of 1.0-1.2 is achieved, yeast can be gently harvested at RT and fresh YPG added to induce protein expression
  - Fix the cap of the culture tube so it is sealed (2nd stop)
  - Using the Sartorius benchtop centrifuge (or other), spin at 500 x g for 3-5 minutes or until the supernatant is clear
  - Gently remove the supernatant with a micropipette (save if interested in metabolite studies at - 20 °C in sterile 15-mL conical tubes)
  - Replace with 5 mL of fresh YPG media (warmed to 30 °C) and place back in the incubator at 30 °C. Continue collecting time point data
12. Pellets can be retained and saved at the end of the study for additional analysis

[Yeast Growth Curve Analysis\\_unlin...](#)