



Jan 06, 2023

Yeast Peptone Dextrose (YPD) medium

DOI

dx.doi.org/10.17504/protocols.io.x54v9d7x1g3e/v1



Andreas Sagen¹

¹University of Oslo



Andreas Sagen

University of Oslo

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.x54v9d7x1g3e/v1>

Protocol Citation: Andreas Sagen 2023. Yeast Peptone Dextrose (YPD) medium. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.x54v9d7x1g3e/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: January 05, 2023

Last Modified: January 06, 2023

Protocol Integer ID: 74794

Keywords: YPD, Yeasts, *Saccharomyces cerevisiae*, *S. cerevisiae*, yeast peptone dextrose, medium yeast peptone, yeast extract, yeast, others *saccharomyces cerevisiae*, many common yeast, glucose, peptone, protein, addition of protein, primary carbon

Abstract

Yeast peptone dextrose is a medium composition used to grow many common yeasts, among others *Saccharomyces cerevisiae*. Glucose is the primary carbon and energy source, while nitrogen and essential amino acids are provided by the yeast extract and peptone. It is nonselective and the addition of protein and yeast extract allows faster growth so that during exponential or log-phase growth, the cells divide every 90 minutes in *S. cerevisiae*.

Guidelines

Follow step by step, unless stated otherwise. Equipment needed should be standard to a microbiology lab.

Materials

Analytical scale, autoclave, bottle(s), weight vessel, LAF bench

Protocol materials

☒ Dextrose Merck MilliporeSigma (Sigma-Aldrich) Catalog #G7021

☒ Yeast Extract Merck MilliporeSigma (Sigma-Aldrich) Catalog #Y0875

☒ Peptone Merck MilliporeSigma (Sigma-Aldrich) Catalog #P5905

☒ Peptone Merck MilliporeSigma (Sigma-Aldrich) Catalog #P5905

☒ Yeast Extract Merck MilliporeSigma (Sigma-Aldrich) Catalog #Y0875

☒ Agar Merck MilliporeSigma (Sigma-Aldrich) Catalog #A1296

Troubleshooting



Safety warnings

! You can mix Dextrose from the beginning with the other compounds, and autoclave together. While this is more time efficient and easier, it is important to take into account the possibility of toxic byproducts produced by the Millard reaction when autoclaving, producing Acrylamide, a probable human carcinogen (IARC Group 2A).

Furthermore, when removing autoclaved components, be sure to take care as these can be very hot. If using antibiotics, use sufficient PPE to protect yourself, as some can be toxic to humans.

Before start

Prepare glassware by cleaning it, and ensure that scale is sufficiently calibrated




100 mL Dextrose solution


1 All compounds are measured using a high precision analytical scale from powdered compounds. Each compound is measured to within 1% of the target weight. All compounds are mixed in a Duran bottle

1.1 Fill the bottle with  60 mL double-distilled water

1.2 Measure  10000 mg Dextrose


Powdered compounds:

 Dextrose **Merck MilliporeSigma (Sigma-Aldrich) Catalog #G7021**

1.3 Add powdered solids into bottle, and use a magnetic mixer with a stir bar to mix for  00:05:00

5m

1.4 Adjust pH while mixing to  5.8 using concentrated sodium hydroxide

1.5 Add distilled water to a total of  100 mL

1.6 Autoclave liquid at  121 °C for  00:15:00

15m


500 mL YPD (broth) medium

2 All compounds are measured using a high precision analytical scale from powdered compounds. Each compound is measured to within 1% of the target weight. All compounds are mixed in a Duran bottle

2.1 Fill the bottle with  300 mL double-distilled water

2.2 Measure  10000 mg Peptone and  5000 mg Yeast extract

Powdered compounds:

 Peptone **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P5905**



Yeast Extract **Merck MilliporeSigma (Sigma-Aldrich) Catalog #Y0875**

2.3 Add powdered solids into bottle, and use a magnetic mixer with a stir bar to mix for

00:05:00

5m

2.4 Adjust pH while mixing to 5.8 using concentrated sodium hydroxide

2.5 Autoclave liquid at 121 °C for 00:15:00

15m

2.6 In a LAF bench, add 100 mL sterile Dextrose solution

2.7 Add sterile water to a total of 500 mL

Note

Cool to 50°C and supplement with antibiotics as appropriate

500 mL YPD (agar) medium

15m

3 All compounds are measured using a high precision analytical scale from powdered compounds. Each compound is measured to within 1% of the target weight. All compounds are mixed in a Duran bottle

3.1 Fill the bottle with 300 mL double-distilled water

3.2 Measure 10000 mg Peptone, 5000 mg Yeast extract and 7500 mg agar

Powdered compounds:

Peptone **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P5905**

Yeast Extract **Merck MilliporeSigma (Sigma-Aldrich) Catalog #Y0875**

Agar **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A1296**

3.3 Add powdered solids into bottle, and use a magnetic mixer with a stir bar to mix for

00:05:00

5m




3.4 Adjust pH while mixing to  5.8 using concentrated sodium hydroxide

3.5 Autoclave liquid at  121 °C for  00:15:00

15m

3.6 In a LAF bench, add  100 mL sterile Dextrose solution

3.7 Add sterile water to a total of  500 mL

Note

Cool to 50°C and supplement with antibiotics as appropriate

Agar can be stored, then reheated to 50°C to be poured