

Nov 26, 2024

Version 1

Yeast Sample Preparation for Fluorescence Live Cell Imaging V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.bp2l629xdgqe/v1>



Mathias Hammer¹, Ammeret Rossouw¹, Brain Tran¹, Anyee Esther Li¹, Graham Latino¹, Pieter Fop van Velde¹, Azra Lari², Ben Montpetit³, David Grunwald¹

¹UMass Chan Medical School, RNA Therapeutics Institute, Worcester, MA, USA;

²University of Alberta, Department of Cell Biology, Edmonton, AB, Canada;

³University of California, Department of Viticulture and Enology, Davis, CA, USA



Mathias Hammer

UMASS Chan Medical School/RTI

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bp2l629xdgqe/v1>

Protocol Citation: Mathias Hammer, Ammeret Rossouw, Brain Tran, Anyee Esther Li, Graham Latino, Pieter Fop van Velde, Azra Lari, Ben Montpetit, David Grunwald 2024. Yeast Sample Preparation for Fluorescence Live Cell Imaging. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bp2l629xdgqe/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 14, 2024

Last Modified: November 26, 2024

Protocol Integer ID: 99800

Keywords: yeast live cell imaging, yeast sample preparation, fluorescence live cell imaging, preparation of yeast cell, live cell imaging, yeast cell, creating single layer yeast suspension, single layer yeast suspensions on the cover glass, yeast, fluorescence, molecule tracking, cell, media preparation, cover glass

Funders Acknowledgements:

NSF

Grant ID: 1917206

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to **protocols.io** is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with **protocols.io**, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

This protocol describes the preparation of yeast cells for fluorescence single molecule tracking live cell imaging. This protocol utilizes linked protocols for specifics on media preparation. Special attention is given to creating single layer yeast suspensions on the cover glass.



Materials

Concanavalin A (ConA) from Canavalia ensiformis:

Sigma-Aldrich

Cat#: C2272-10MG

Lot#: 129H0322

UltraPure Distilled Water

Cover slides ø 25mm, No. 1.5 or 1.5H

Attofluor Cell Chamber

Thermo Fisher Scientific

Cat#: A7816

Equipment:

10 ml centrifugation tube

10 ml pipette

Spectrometer: Molecular Devices SpectraMax m5

Plasma cleaner: Quorum Emitech K100X glow discharger
mixer

bench top centrifuge



Stock sample preparation

2d 12h

- 1 Plate the yeast strain of interest onto pre-imaging solid agar plate.

Pre-Imaging Solid Growth Medium - Yeast protocol:

<https://dx.doi.org/10.17504/protocols.io.8epv5r9wjg1b/v1>

- 1.1 Incubate at  25 °C until colonies are visible. (~2-3 days  60:00:00).

2d 12h

Note

Temperature and duration depending on the specific yeast strain.

- 1.2 The agar plate with grown yeast can be stored, sealed in  4 °C fridge, and can be used as stock for 2-3 months.

Preparation for imaging the next day

12h

- 2 Pick a colony from your stock sample and dilute it in  5 mL liquid pre-imaging medium in a test tube.

Note

Choose the pre-imaging medium according to your labeling system. If it includes labeled pp7-stem loop (pp7-SL), its expression can be suppressed through higher Methionine concentrations in the medium [2]. A recipe for such a medium provides the protocol below that note.

For otherwise labeled yeast it is recommended to use media with higher concentrations of Adenine to reduce autofluorescence [1]. You can follow the "Pre-Imaging Liquid Growth Medium":

<https://dx.doi.org/10.17504/protocols.io.81wgbzxe1gpk/v1>

Pre-Imaging Liquid Growth Medium for Methionine promoted Labeling Systems - Yeast protocol:

<https://dx.doi.org/10.17504/protocols.io.kqdg329yqv25/v1>

- 2.1 Incubate at  200 rpm, 25°C  Overnight .

12h



Note

The incubation temperature depends on the yeast strain!

Yeast conditioning

4h 5m

- 3 Dilute the cultured yeast to $OD_{600} = 0.25 - 0.3$ and incubate  200 rpm, 25°C, 03:00:00 -  04:00:00 . Let the sample grow to an OD_{600} of 0.8 -1.0 to establish exponential growth phase.

4h

Sample preparation

4h 5m

- 4 Pipette  1.5 mL of the incubated yeast into an Eppendorf-tube.
- 4.1 Centrifuge at  560 x g, 23°C, 00:05:00 60
- 4.2 Remove supernatant and re-dilute the pellet in  750 μ L pre-imaging medium.
- 4.3 Plate the re-diluted cells into a cell chamber with a ConA prepared cover slide.
- Protocol for the ConA cover slide preparation:
<https://dx.doi.org/10.17504/protocols.io.kxygxy3rzi8j/v1>
- 4.4 Cover the cell chamber with Kin wipe and let the cells settle in an incubator at  25 °C for  00:10:00 .

5m**10m**

Note

Incubator temperature depends on the Yeast strain and can differ in respect to the experiment.

Sample


10m

5 Remove the supernatant with a 1 ml pipette.

5.1 Tilt the chamber slightly and add  1 mL of medium from the top corner. Remove the medium at the bottom corner to remove loose cells.

5.2

Repeat the last step to remove the second layer of cells. Continue washing the sample as described, until you expect a single layer of cells in the central area of microscope slide.



Note

The pressure to remove the cells from the cover slide respectively the second layer of cells differs between yeast strains and the age (~use cycles) of the ConA.

Note

The removal of the second layer of cells reduces background fluorescence and increases the probability to have the nuclei in the same plan. Space around the cell groups allows the yeast to grow in xy direction and less into z direction creating a new second layer, which allows for longer observation duration.

6 Cautiously add  1 mL refractive index adjusted imaging medium (RI=1.41).

Refractive index adjusted imaging medium: Iodixanol (RI ~ 1.41 and 1.42) - Yeast:



<https://dx.doi.org/10.17504/protocols.io.5qpvek369l4o/v1>

Note

The ideal refractive index for the imaging medium depends on your experimental setup, especially on the immersion medium.

The RI of the medium can be adjusted by varying the RI gradient substance (here Iodixanol). Examples:

- RI = 1.4 <https://dx.doi.org/10.17504/protocols.io.dm6gpz3kjlzp/v1>
- RI > 1.42 <https://dx.doi.org/10.17504/protocols.io.3byl4933ogo5/v1>

- 7 Place the the chamber in the stage incubator. Wait until the temperature reached equilibrium before you start imaging.

Protocol references

[1] Kokina, Agnese et al. "Adenine auxotrophy—be aware: some effects of adenine auxotrophy in *Saccharomyces cerevisiae* strain W303-1A." *FEMS yeast research* 14.5 (2014): 697-707.

doi:10.1111/1567-1364.12154

[2] Lari, Azra, et al. "Live-Cell Imaging of mRNP–NPC Interactions in Budding Yeast." *Imaging Gene Expression: Methods and Protocols* (2019): 131-150.

doi.org/10.1007/978-1-4939-9674-2_9