

Sep 16, 2025

🌐 Yeast suspension for SEM (with HMDS)

DOI

<https://dx.doi.org/10.17504/protocols.io.n2bvje2epgk5/v1>



Alexander Mironov¹

¹University of Manchester



Alexander Mironov

University of Manchester

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n2bvje2epgk5/v1>

Protocol Citation: Alexander Mironov 2025. Yeast suspension for SEM (with HMDS). [protocols.io](https://dx.doi.org/10.17504/protocols.io.n2bvje2epgk5/v1)

<https://dx.doi.org/10.17504/protocols.io.n2bvje2epgk5/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: September 16, 2025

Last Modified: September 16, 2025

Protocol Integer ID: 227339

Keywords: SEM, HMDS, yeast, yeast suspension for sem, suspension of yeast, yeast suspension, rt sem imaging, using hmd, sem

Abstract

Suspension of yeast prepared for RT SEM imaging using HMDS drying

Materials

- 1) Stock solution of glutaraldehyde (could be from 10 to 25%)
- 2) Ethanol
- 3) Hexamethyldisilazane (HMDS)
- 4) OsO4 stock solution in water (not less than 1%)
- 5) 12mm metal SEM stubs
- 6) Sticky carbon adhesive tabs
- 7) Silver DAG
- 8) Round 10mm glass coverslips
- 9) Glow discharger
- 10) Sputter coater



Substrate preparation

1h 2m

- 1 Use glass coverslips that fit to your regular SEM metal stubs. If a stub is 12mm in diameter use 10mm round coverslips.
- 2 Glow discharge coverslips for about 1-2 min with standard settings (20-25mA).
- 3 Briefly soak coverslips in Poly-L-lysine commercial solution (any), remove the solution excess by gentle blotting and leave them to dry in the open air for 30 min.

1m

1m

Sample fixation

1h

- 4 Put coverslips on parafilm flat and apply about 30-50ul of yeast suspension on top of them taking care not to overspill on the parafilm. Wait for about 20-30 min, so yeast cells will drawn to the coverslip.
- 5 Apply 2-5 ul of 10%-12% glutaraldehyde directly into yeast suspension and wait until the fixative will reach equilibrium concentration throughout.

30m

30m

Osmium postfixation

1h 20m

- 6 Wash the coverslips with distilled water buffer 2 times for 5 min each.
- 7 Put coverslips into appropriate glass container and apply 1%OsO₄ solution in water.
- 8 Wash coverslips with water 2 times for 5 min each.

10m

1h

10m

Dehydration

1h 15m

- 9 Incubate coverslips in 50% ethanol
- 10 Incubate coverslips in 75% ethanol

15m

15m



11 Incubate coverslips in 90% ethanol

15m

12 Incubate coverslips in 100% ethanol 2 times for 15 min each

30m

Drying

1h

13 Incubate coverslips in 50% HMDS in ethanol

30m

14 Incubate coverslips in 100% HMDS

15m

15 Incubate coverslips in 100% HMDS

15m

16 Remove most of HMDS from the containers with coverslips. Leave containers overnight opened in the fume cupboard.

17 Remove coverslips from the containers and glue them to SEM stubs with carbon sticky tape taking care about the side with attached cells.

Coating

18 To get better conductivity paint the conductive bridge between top of a coverslip and metal SEM stub with silver DAG. Wait until it is completely dry.

19 Use your sputter coater to create 7-9nm of metal coating with you favourite target (Au/Pd will do good).

20 Sample can be imaged with most room temperature SEMs using standard settings for moderately conductive samples.