

Sep 22, 2019

Yeast transformation

DOI

dx.doi.org/10.17504/protocols.io.7jkhkkw

Marijn Ceelen¹

¹Wageningen University

iGEM Wageningen 2019



Santi Castanedo

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.7jkhkkw

Protocol Citation: Marijn Ceelen 2019. Yeast transformation. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.7jkhkkw>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 22, 2019

Last Modified: September 22, 2019

Protocol Integer ID: 27980

Keywords: Transformation, DNA, yeast, *Saccharomyces cerevisiae*

Abstract

This is a protocol for the transformation of yeast cells with linear or circular DNA.



Materials

MATERIALS

⊗ Polyethylene Glycol 3350

⊗ Salmon Sperm DNA Carrier **Thermofisher Catalog #15632011**

⊗ 100ml Lithium acetate [1M] **G-Biosciences Catalog #R039**

⊗ YPD Broth **Thermo Fisher Catalog #A1374501**

Synthetic complete (SC) medium + agar

Before start

Have all the DNA fragments that will be transformed in high enough concentrations to reach from 0.5-4 µg DNA per transformation.

Have SC plates without the right amino acid that will be used as a selective marker.

- 1 Make a transformation mix consisting of:

DNA	34 μ L (0.5- 4 μ g per fragm ent)
Salmon sperm DNA 2 mg/mL	50 μ L
PEG- 3350	240 μ L
1M LiOAc	36 μ L
MQ water	up to a final volume of 360 μ L

- 2 Centrifuge the 1.5 mL tubes with 100 μ L of yeast competent cells for 2 minutes at 5000x g and remove the supernatant.
- 3 Add 360 μ L of the transformation mix into the 1.5 mL with the cells and vortex thoroughly.
- 4 Heat-shock the solution at 42 °C for 40 minutes.
- 5 Centrifuge the cells for 2 minutes at 5000x g and remove the supernatant.
- 6 Resuspend the cells in 1 mL of YPD medium. Divide the volume into two different 1.5 mL tubes.
- 7 Leave cells to recover for 1.5 hours at 30 °C and 300 rpm.



- 8 Centrifuge cells for 2 minutes at 5000x g and decant supernatant, leaving a bit of the medium inside the tubes.
- 9 Resuspend cells in the resting medium and plate 50 μ L onto a SC plate without the amino acid used as a selection marker.
- 10 Incubate the cells for 3 days at 30 °C.