



Sep 11, 2025

2nd generation Lentiviral production

DOI

<https://dx.doi.org/10.17504/protocols.io.ewov11z4ovr2/v1>

Hannah Lucas-Clarke¹

¹UCL/The Francis Crick Institute



Hannah Lucas-Clarke

UCL/The Francis Crick Institute

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ewov11z4ovr2/v1>

Protocol Citation: Hannah Lucas-Clarke 2025. 2nd generation Lentiviral production. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.ewov11z4ovr2/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 11, 2025



Last Modified: September 11, 2025

Protocol Integer ID: 227021

Keywords: lentivirus for future transduction, 2nd generation lentiviral production, lentiviral production, lentivirus, plasmid, grna, future transduction, grna of interest

Abstract

This is a protocol to amplify and purify plasmids (with gRNA of interest) and package into lentivirus for future transduction.

Materials

Equipment

Centrifuges/Ultracentrifuge (RT and 4 degree)
Tissue Culture Hood appropriate for viral work
Bunsen burner
37 degree incubators one static and one shaking

Materials

Single-gRNA Expression Lentiviral Vector
pxPAX2 (Addgene, #12260)
pMD2.G (Addgene, #12259)
Plastic loops
Agar plates with the appropriate antibiotic
Ampicillin/Or other appropriate antibiotic
Liquid Broth
Large plastic centrifuge tubes
PureLink Maxi-Prep kits
Isopropanol
Ethanol
TE bufffer
HEK 293T cells
FuGene HD (Promega, cat. no. E2311)
DMEM (Gibco, cat. no. 31966021)
FBS
Pen-Step (cat. no. 15070063)
Glutamax (cat. no. 35050038)
10cm cell culture dishes
Falcon and Eppendorf tubes
0.45um filters
PBS



Vector Design

- 1 Design Mammalian Single-gRNA Expression Lentiviral Vector. Typically we use Vector Builder which deliver E.Coli stocks.

Bacteria growth on agar plates

1d

- 2 For this you need plastic loops, bunsen burner and agar plates with the **appropriate** antibiotic based on vector design.
- 3 Sterilise plastic loop over bunsen flame right before dipping into bacteria stock
- 4 Dip into stock and streak onto agar plate. Make sure the agar plate is labelled!
- 5 Incubate agar plate overnight at 37 degrees for bacteria to grow in colonies.

Liquid Broth

1d

- 6 Add 500uL of ampicillin (**or appropriate antibiotic**) to bottle of LB broth, place LB broth into large glass flask (+1 extra for control)
- 7 Remove agar plates from 37 degrees and check for successful growth of colonies
- 8 Take a pipette tip and sterilise over bunsen burner flame, pick colony and place into flask - cover flask with foil (loosely)
- 9 As a control take a tip, sterilise and place into LB flask without picking a colony
- 10 Incubate flasks shaking at 37 degrees overnight

Maxi-Prep

2h

- 11 Large plastic tubes that fit into the centrifuge are needed for this step. Remove the flasks from the incubator - there should be visible growth of bacteria in all but control flask. Transfer contents of flask into plastic centrifuge tube (**without the tip**)
- 12 Weigh each tube - add PBS to make up difference (this is to ensure the centrifuge is balanced. Centrifuge at 4000g 4 degrees for 10 minutes.
- 13 We use the PureLink Maxi-Prep kits from Thermo Fisher Scientific. Place the red plastic top onto a glass beaker and then the column on top. Add 30mL of EQ1 buffer to each column.
- 14 Remove plastic tubes from centrifuge - discard supernatant. Resuspend sediment in 10mL of R3 buffer and pipette up and down until cloudy.
- 15 Add 10mL of L7 lysis buffer to the tube - invert to mix 4-5 times. Incubate for 5 minutes
- 16 Add 10mL of N3 buffer and shake. Transfer contents to a labelled falcon tube and centrifuge 12000g for 10 minutes (or 4000g for 30 minutes)
- 17 Remove falcons from centrifuge - there should be clear layers. Pipette liquid into column, avoid top layer and bottom layer by tilting tube. Allow to filter through column
- 18 Add 60mL wash buffer to each column. Flow through can be discarded.
- 19 Take washed column and place on top of a labelled falcon tube. Add 15mL elution buffer to each column. Once eluted - throw column.
- 20 Add 10.5mL isopropanol to each falcon - shake and centrifuge for 30 minutes at 12000g at 4 degrees.
- 21 DNA should have pelleted at base of falcon. Discard supernatant. Add 5mL 70% ethanol to each DNA pellet. Further centrifuge for 10 minutes to remove any remaining bacteria.
- 22 Remove all liquid and allow pellet to air dry for 10 minutes.
- 23 Add 200-500uL TE buffer to resuspend DNA pellet, transfer to eppendorf. Nanodrop DNA to get the concentration. Store in freezer.



Lentiviral Production

1w

- 24 Prior to this make sure you have HEK 293T cells in culture. Media: 500mL DMEM, 10% FBS, 1% Pen-Step, 1% Glutamax. When confluent, plate into 10cm cell culture dish. Good to do this on a Friday so they are 80% confluent on a Monday ready for transfection. Do multiple dishes per plasmid - i.e. 3 dishes of HEK. Remember to scale the rest of the protocol depending on the number of dishes you are doing.
- 25 3-4 hours before transfection change media on HEKs (10mL) - aim for transfection to be in the afternoon (4pm)
- 26 Add 18uL FuGene HD to 300uL DMEM and incubate for 5 minutes.
- 27 Add 3ug transfer plasmid with 2ug psPAX2 and 1ug pMD2.G in 300uL of DMEM
- 28 Mix both reactions gently and incubate for 30 minutes at room temp
- 29 Add mixture (600uL) to HEK cells
- 30 18 hours after the transfection replace the media
- 31 Turn on Ultracentrifuge - set to 4 degrees.
Harvest the media after 48 hours into 50mL falcons (curve bottom) and centrifuge 1500rpm for 3 minutes at room temperature
- 32 Filter supernatant through 0.45um filter
- 33 Transfer to ultracentrifuge tubes. Weigh to make sure this is balanced.
- 34 Centrifuge for 3 hours at 48000g
- 35 Pellet will be hard to see - circle when you remove from centrifuge. Remove supernatant in the hood and resuspended pellet in ice cold sterile PBS. 3 dishes - 850uL of PBS. Leave these on ice in fridge/cold room overnight.



36 Next day aliquot virus - snap freeze and transfer to -80.