

Jun 13, 2016

## CviJI Purification From IL-3A Virus Infected NC64A Chlorella

DOI

[dx.doi.org/10.17504/protocols.io.er3bd8n](https://dx.doi.org/10.17504/protocols.io.er3bd8n)

David Dunigan and Irina Agarkova<sup>1</sup>

<sup>1</sup>The University of Nebraska-Lincoln

VERVE Net



Irina Agarkova

UNL

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.er3bd8n>

**Protocol Citation:** David Dunigan and Irina Agarkova 2016. CviJI Purification From IL-3A Virus Infected NC64A Chlorella. [protocols.io https://dx.doi.org/10.17504/protocols.io.er3bd8n](https://dx.doi.org/10.17504/protocols.io.er3bd8n)

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working



**Created:** March 25, 2016

**Last Modified:** March 17, 2018

**Protocol Integer ID:** 2587

**Keywords:** cviji purification, nc64a chlorella, cviji, chlorella, virus, nc64a, purification

## Guidelines

### **MATERIALS AND ASSAY CONDITIONS:**

1) 7 hour IL-3A virus infected NC64A, pellets frozen at -80°C

2) Buffer A:

- 10 mM Tris-Acetate, pH 8.0
- 10 mM 2-M
- 50 µg/mL PMSF

3) Buffer B:

- 20 mM Tris-Acetate, pH 8.0
- 0.5 mM EDTA
- 7 mM 2-ME, 10% Glycerol

4) Buffer B, pH 8.5:

- 20 mM Tris-Acetate, pH 8.5
- 0.5 mM EDTA
- 7 mM 2-ME

5) Storage buffer:

- 20 mM Tris-Acetate, pH 8.0
- 0.5 mM EDTA
- 0.1 mM DTT
- 50 mM KOAc
- 5 mM MgAc, 50% Glycerol
- 100 µg/mL BSA (the BSA is added after dialysis into the storage buffer from a 10 mg/mL stock)

6) 4 M NaCl

7) 28% Polyethyleneglycol (PEG) 8000

8) Assay buffer and conditions:

1X CviJI assay buffer:

- 20 mM GlycylGlycine, pH 8.5 (with KOH)
- 10 mM MgAc
- 0.1 mM DTT
- 50 mM KOAc
- 100 µg/mL casein (optional)

All assays are carried out in 20.0 µL volumes with 1 µg of pUC19 DNA as substrate for 60 to 120 min at 25°C. The assays are electrophoresed on 2.0% agarose gels (100 mL gels with double 20 lane combs) for 1 hour at 100 volts in 1X TPE buffer. Gels are stained with 0.5 µg/mL ethidium bromide for 30 min and photographed on a UV light box.

## Troubleshooting

### Before start

7 hour IL-3A virus infected NC64A, pellets frozen at -80°C



## 1 Prepare Buffer A:

### Protocol

NAME

**CviJI Buffer A**

CREATED BY

Irina Agarkova

[Preview](#)

1.1 10 mM Tris-Acetate, pH 8.0

1.2 10 mM 2-ME

1.3 50 µg/mL PMSF

## 2 Prepare Buffer B:

### Protocol

NAME

**CviJI Buffer B**

CREATED BY

Irina Agarkova

[Preview](#)

2.1 20 mM Tris-Acetate, pH 8.0

2.2 0.5 mM EDTA

2.3 7 mM 2-ME, 10% Glycerol



### 3 Prepare Buffer B, pH 8.5:

#### Protocol

NAME

**CviJI Buffer B, pH 8.5**

CREATED BY

Irina Agarkova

[Preview](#)

3.1 20 mM Tris-Acetate, pH 8.5

3.2 0.5 mM EDTA

3.3 7 mM 2-ME

### 4 Prepare Storage Buffer:

#### Protocol

NAME

**CviJI Storage Buffer**

CREATED BY

Irina Agarkova

[Preview](#)

4.1 20 mM Tris-Acetate, pH 8.0

4.2 0.5 mM EDTA

4.3 0.1 mM DTT



- 4.4 50 mM KOAc
- 4.5 5 mM MgAc, 50% Glycerol
- 4.6 100 µg/mL BSA (the BSA is added after dialysis into the storage buffer from a 10 mg/mL stock)

## 5 Prepare 1X CviJI Assay Buffer:

### Protocol

NAME

**1X CviJI Assay Buffer**

CREATED BY

Irina Agarkova


Preview

- 5.1 20 mM GlycylGlycine, pH 8.5 (with KOH)
- 5.2 10 mM MgAc
- 5.3 0.1 mM DTT
- 5.4 50 mM KOAc
- 5.5 100 µg/mL casein (optional)
- 6 Thaw 7 hour IL-3A virus infected NC64A *chlorella* and suspend in MSK flasks with Buffer A.

**Note**

Suspend with 20 mL per flask per  $1.0\text{--}1.5 \times 10^{11}$  infected cells.


- 7 Homogenize the cells in the MSK mechanical homogenizer with 15 gm of 0.3 mm glass beads at 4,000 rpm for 90 sec (2 X 45 sec) with CO<sub>2</sub> cooling.

 00:01:30

- 8 Recover the homogenate to clean tubes.

- 9 Wash the glass beads 3X with 5 mL of Buffer A and combine with the homogenate.

- 10 Centrifuge the homogenate in the Sorvall SS34 rotor at 10,000 rpm, 20 min, 4°C.

 00:20:00


- 11 Save the supernatant.

- 12 Adjust the homogenate supernatant to 70% saturation with (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub> at 4°C with gentle stirring.


**Note**

Add the (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub> gradually.

- 13 Incubate at 4°C for 60-90 min without stirring.

 01:30:00

- 14 Centrifuge the material in the Sorvall SS34 rotor at 10,000 rpm 10 min, 4°C.

 00:10:00

- 15 Save the pellet.

- 16 Suspend the pellets with Buffer A.



17 Per mL of suspension add: 0.45 mL of 4 M NaCl and 0.45 mL of 28% PEG 8000 (heated to 65°C).

18 Mix gently by inversion for 5-10 min.

00:10:00

19 Centrifuge the material in the Sorvall SS34 rotor at 10,000 rpm, 10 min, 4°C.

00:10:00

20 Save the supernatant.

#### Note

If there is small particulate matter in the supernatant, centrifuge a second time as before and save the supernatant.

21 Dilute the supernatant with 10-15 volumes of Buffer B to reduce the NaCl concentration.

22 Load the material overnight onto a Heparin-Sepharose column equilibrated with Buffer B in the cold room.

18:00:00

23 Elute the Heparin-Sepharose column with Buffer B using a 0-2.0 M KOAc gradient.

24 Assay the column fractions and pool the active fractions.

#### Protocol

NAME

**CviJI Assay Conditions**

CREATED BY

Irina Agarkova

Preview

24.1 All assays are carried out in 20.0  $\mu$ L volumes with 1  $\mu$ g of pUC19 DNA as substrate for 60 to 120 min at 25°C.

02:00:00





- 24.2 The assays are electrophoresed on 2.0% agarose gels (100 mL gels with double 20 lane combs) for 1 hour at 100 volts in 1X TPE buffer.

01:00:00

- 24.3 Gels are stained with 0.5 µg/mL ethidium bromide for 30 min and photographed on a UV light box.

00:30:00

- 25 Dilute the pooled fractions with 10-15 volumes of Buffer B to reduce the salt concentration.

- 26 Load the material overnight onto a Blue-Sepharose column equilibrated with Buffer B in the cold room.

18:00:00

- 27 Elute the Blue-Sepharose column with Buffer B using a 0-2.0 M KOAc gradient.

- 28 Assay the column fractions and pool the active fractions.

#### Protocol

NAME

**CviJI Assay Conditions**

CREATED BY

Irina Agarkova

Preview

- 28.1 All assays are carried out in 20.0 µL volumes with 1 µg of pUC19 DNA as substrate for 60 to 120 min at 25°C.

02:00:00


- 28.2 The assays are electrophoresed on 2.0% agarose gels (100 mL gels with double 20 lane combs) for 1 hour at 100 volts in 1X TPE buffer.

01:00:00

- 28.3 Gels are stained with 0.5 µg/mL ethidium bromide for 30 min and photographed on a UV light box.

00:30:00



- 29 Dilute the pooled fractions with 10-15 volumes of Buffer B, pH 8.5 to reduce the salt concentration.
- 30 Load the material overnight onto a Q-Sepharose column equilibrated with Buffer B, pH 8.5 in the cold room.  
 18:00:00
- 31 Elute the Q-Sepharose column with Buffer B, pH 8.5 using a 0-2.0 M KOAc gradient.
- 32 Assay the column fractions and pool the active fractions.

#### Protocol





NAME

### CviJI Assay Conditions

CREATED BY

Irina Agarkova

Preview

- 32.1 All assays are carried out in 20.0  $\mu$ L volumes with 1  $\mu$ g of pUC19 DNA as substrate for 60 to 120 min at 25°C.  
 02:00:00
- 32.2 The assays are electrophoresed on 2.0% agarose gels (100 mL gels with double 20 lane combs) for 1 hour at 100 volts in 1X TPE buffer.  
 01:00:00
- 32.3 Gels are stained with 0.5  $\mu$ g/mL ethidium bromide for 30 min and photographed on a UV light box.  
 00:30:00
- 33 Concentrate the pooled enzyme by dialysis overnight into storage buffer at 4°C.  
 18:00:00
- 34 Add BSA (10 mg/mL) to a final concentration of 100  $\mu$ g/mL.
- 35 Store the enzyme at -20°C.

