

Feb 09, 2016

Obtaining pure cyanophage stocks (plaque purification)

 In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.dqq5vv

Mathias Middelboe¹, Amy M. Chan¹, and Sif K. Bertelsen¹

¹Manual of Aquatic Viral Ecology

VERVE Net

Suttle Laboratory of Mar...



Amy Chan

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dqq5vv>

External link: http://www.aslo.org/books/mave/MAVE_118.pdf

Protocol Citation: Mathias Middelboe, Amy M. Chan, and Sif K. Bertelsen 2016. Obtaining pure cyanophage stocks (plaque purification). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.dqq5vv>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: September 03, 2015

Last Modified: November 09, 2017

Protocol Integer ID: 1520

Keywords: obtaining pure cyanophage stock, pure cyanophage stock, plaque purification, plaque, purification

Troubleshooting



- 1 Make a dilution series of the lysate (assume 10^4 to 10^5 PFU per mL in the plaque lysate)

Note

Use this to perform a second round of plaque assays to purify the phage.

- 2 Repeat the plaque purification procedure 2 more times to ensure that the cyanophage isolated is clonal.
- 3 Repeat the plaque purification procedure again.
- 4 Finally, prepare a primary cyanophage stock using lysate from the final purification via method A "**Liquid Amplification**" **OR** method B "**Plate Amplification**"
- 5 Method A: liquid amplification

Protocol

NAME

Liquid Amplification

CREATED BY

Amy Chan

Preview

- 5.1 Add some of the lysate to target host in liquid culture.
- 5.2 After the culture has lysed, remove cell debris via centrifugation.
- 5.3 Filter sterilize the stock.
- 5.4 Store at 4°C until further analysis.
- 6 Method B: plate amplification



Protocol

NAME

Plate Amplification

CREATED BY

Amy Chan

[Preview](#)

- 6.1 Prepare plaque assays with a dilution series of lysate from the final purification.


Note

Plates with confluent lysis of the host lawn (typically ca. 10^4 PFUs) can then be used to obtain cyanophage stocks by elution of phages from the plates.

- 6.2 Add 5 mL sterile seawater to the plate.

- 6.3 Scrape off the top agar layer into the seawater.

- 6.4 Leave at 4°C overnight.

 18:00:00

- 6.5 Remove agar and cell debris by centrifugation.

- 6.6 Filter sterilize the stock.

- 6.7 Store at 4°C until further analysis.

- 7 Titer the final stock via plaque assay.

- 8 Cyanophage stocks stored at 4°C in the dark are stable for at least a year.