

Feb 02, 2016

## Roseobacter Screening Of Surface Waters For Viruses

DOI

[dx.doi.org/10.17504/protocols.io.dbt2nm](https://dx.doi.org/10.17504/protocols.io.dbt2nm)

Matthew Sullivan<sup>1</sup>

<sup>1</sup>Matthew Sullivan Lab, University of Arizona, Ohio State University

VERVE Net

Sullivan Lab



Verve Team

University of Arizona

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dbt2nm>

**Protocol Citation:** Matthew Sullivan 2016. Roseobacter Screening Of Surface Waters For Viruses. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.dbt2nm>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working



**Created:** July 08, 2015

**Last Modified:** November 09, 2017

**Protocol Integer ID:** 1107

**Keywords:** roseobacter screening of surface water, roseobacter screening, surface water, virus

## Guidelines

### Notes:

Goal initially is to screen several seawater samples on several roseobacter hosts to determine breadth of rosephage present at each location and then to focus on a single site.

Nomenclature for viruses should include a designation for the site of seawater collection followed by a number and the host it was isolated on.

Put information into an Excel worksheet. Create a separate Excel file for each source of virus. Use a different tab in each worksheet for different hosts.

Turbid plaques indicate possible lysogenic phage that may be of interest to Feng Chen's laboratory in Maryland. Keep the plaque plug at 4°C. Take 50-100 µl from a few and grow up small 1 ml cultures to cryopreserve (glycerol) them. Here trying to get the putative lysogenic cells (prophage containing confers resistance to colonies of cells within the plaque = turbid).

## Troubleshooting



## Growth

- 1 Pull *Roseobacter* cultures from freeze and grow in 5 ml broth cultures

### Note

Matt has 14 isolates from Wendy Ye in Mary Ann Moran's lab in Georgia. Two of them are from open ocean (CCS1 and CCS2) and they are grown at 20-22°C in 1/10 YTSS. All others are from coastal waters and are grown at 30°C in 1/2 YTSS. For virus screening, grow them to log phase ( $OD_{600} = 0.6 - 0.8$ ). Some will grow quickly (overnight) while others will take 1-4 days to reach a log phase.

## Screening

- 2 Perform virus screening in sterile 96 well flat-bottom tissue culture treated microtiter plates

### Note

Viruses will be screened from seawater collected from various sources.

- 3 Store seawater filtrates at 4°C

### Note

Ice crystal formation at -20°C **bad**.

## Absorption

- 4 Absorb the viruses to the bacteria
- 5 Place 15 µl of log phase cells and 15 µl of seawater filtrate (containing virus) into wells of a 96 well plate

### Note

Uninfected controls should surround infected wells on each plate in a checkerboard pattern.

- 6 Allow the viruses to absorb to the bacteria for 1 hour



01:00:00



- 7 Add 200  $\mu$ l growth media (1/2 or 1/10 YTSS) to the wells
- 8 Seal the plate with parafilm and incubate with shaking
- 9 Look for infected wells that show lysis with respect to uninoculated controls

## Enrichment


- 10 Enrich for viruses by transferring from initial 96 well plate to another 96-well plate containing fresh liquid media ("frogging" into new wells)

### Note

Does not require fresh log-phase host cells, the cells come with the transfer. Ideally frog **before** cells into stationary phase.

### Note

Enrichment for viruses may take several rounds of culturing because we're assuming initial concentrations of phages in the seawater are low.

- 11 Once enriched sufficiently, filter out cells with 0.2 $\mu$ m filtration into 1.5ml microfuge tubes
- 12 Spin for 5 mins. at max to pellet cells  
 00:05:00
- 13 Transfer supernatant to fresh tube for storage
- 14 Store at 4°C

## Plaque Purification

- 15 Plaque-purify wells containing virus by growing cells plus virus on solid media



- 16 Prepare solid agar media and 0.5% overlay agar
- 17 Absorb virus to cells using as a starting point, 25  $\mu$ l virus from 96 well plate enrichments and 200  $\mu$ l log-phase host bacteria for 1hr

01:00:00

#### Note

Try growing lawns first, maybe needs to be up to 1 ml cells.

- 18 Mix the virus-absorbed cells in the 5 ml overlay agar
- 19 Pour over the solid media

- 20 Incubate and look for plaque formation

#### Note

Will have several possible results: no plaques (not likely if using virus-enriched cultures from 95 well plates), too many plaques to purify them (in which case, a dilution of virus is made and plaque assay repeated), or well-resolved plaques. If well-resolved plaques, then cored from agar using Pasteur pipet and dispensed into 100 $\mu$ l YTSS broth in 1.5ml tubes.

- 21 Pick representatives of all plaque types present and record appearance of the plaques

#### Note

Clear, well-lysed plaques are of the most interest to us. Turbid plaques may indicate lysogenic phage and some should be picked, but these will be sent to another laboratory for further study (see guidelines).


- 22 Store the broth that contains the agar plugs at 4°C and allow the viruses to diffuse out of the agar into the broth

## DNA Purification Determination

- 23 Once the viruses are plaque-purified they can be scaled-up in 5 ml broth cultures for DNA purification



- 24 Take 25  $\mu$ l of the plaque-purified virus (agar plug in  $\sim$ 100  $\mu$ l broth) and absorb to 200  $\mu$ l log-phase bacteria for 1 hr

 01:00:00

- 25 Add to 5 ml broth and grow with shaking

- 26 Centrifuge out cells and filter culture through 0.2  $\mu$ m filter

- 27 SYBR stain a portion of the filtrate to determine how much virus is present

#### Note

If a lot of virus is present, then may be enough to purify DNA.

## DNA Purification

- 28 Purify DNA using the Promega Wizard Lambda DNA kit

#### Note

Goal is to obtain about 1  $\mu$ g DNA for future work. If 5 ml does not yield that amount, will need to grow up more (e.g. 225 ml volume) and repeat DNA purification.