

Jan 12, 2023

Version 1

ARTIC-like Bacillus anthracis MLVA amplicon sequencing protocol for MinION V.1

DOI

dx.doi.org/10.17504/protocols.io.3byl4jzozlo5/v1

Ágnes Nagy¹, Gábor Tóth²

¹Hungarian Defence Forces Medical Centre; ²National Laboratory of Virology (Hungary), University of Pecs



Ágnes Nagy

Hungarian Defence Forces Medical Centre

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl4jzozlo5/v1>

Protocol Citation: Ágnes Nagy, Gábor Tóth 2023. ARTIC-like Bacillus anthracis MLVA amplicon sequencing protocol for MinION. protocols.io <https://dx.doi.org/10.17504/protocols.io.3byl4jzozlo5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: January 05, 2023

Last Modified: January 12, 2023

Protocol Integer ID: 74817

Keywords: Bacillus anthracis MLVA, environmental samples, amplicon sequencing, bacillus anthracis vntr loci, protocol for bacillus anthracis mlva, like bacillus anthracis mlva amplicon, related bacillus anthraci, bacillus anthracis mlva, bacillus anthraci, protocol on hungarian virulent bacillus anthraci, clear cultures of bacillus anthraci, hungarian virulent bacillus anthraci, 34f2 bacillus anthracis vaccine strain, genome sequencing, generating pcr amplicon, silico analysis of whole genome sequencing, whole genome sequencing, clear isolates of typeable strain, pcr amplicon, capillary electrophoresis of pcr amplicon, sequencing protocol, generating amplicon, typeable strain, equipped high biosafety level laboratory, high biosafety level laboratory, applicable amplicon, locus variable number of tandem repeat, tandem repeat, longer amplicon, genome, level subtyping methods for outbreak

Funders Acknowledgements:

Hungarian Defence Forces Medical Centre

Abstract

Multiple-Locus Variable Number of Tandem Repeats (VNTR) Analysis (MLVA) is one of the gold standard strain-level subtyping methods for outbreak-related *Bacillus anthracis* strains. The repeat numbers of 31 VNTR loci can be determined by capillary electrophoresis of PCR amplicons spanning repeat regions or in silico analysis of whole genome sequencing (WGS) data. However these methods require clear isolates of typeable strains and can be performed in fixed well-equipped high biosafety level laboratories.

We developed field-applicable amplicon sequencing protocol for *Bacillus anthracis* MLVA typing directly from environmental samples without isolating clear cultures of *Bacillus anthracis* strains. 62 primers were used for generating PCR amplicons for 31 *Bacillus anthracis* VNTR loci, according to MLVA31 typing scheme described by Beyer et al. 2012. The primers generating amplicons longer than 200 bp were used from MLVA31 typing scheme (44 primers for 22 VNTR loci). For amplicons shorter than 200 bp, primers were redesigned to generate longer amplicons (between 300-700 bp) suitable for MinION sequencing.

We optimized and tested this protocol on hungarian virulent *Bacillus anthracis* strains, a 34F2 *Bacillus anthracis* vaccine strain, and on spiked environmental samples in Hungarian Defence Forces field-deployable laboratory.

Citation

Beyer W, Bellan S, Eberle G, Ganz HH, Getz WM, Haumacher R, Hilss KA, Kilian W, Lazak J, Turner WC, Turnbull PC (2011). Distribution and molecular evolution of *bacillus anthracis* genotypes in Namibia.. PLoS neglected tropical diseases.

<https://doi.org/10.1371/journal.pntd.0001534>

LINK



Materials

Primers 25nm, desalted, ideally LabReady formulation from IDT



BaMLVA_primers.xlsx

Qiagen DNeasy Blood&Tissue kit Qiagen 69504

Q5 Hot Start HF Polymerase NEB M0493S

NEBNext Ultra II End Repair/dA-Tailing Module NEB E7546S

NEBNext Ultra II Ligation Module NEB E7595S

NEBNext Quick Ligation Module NEB E6056S

QuantiFluor ONE dsDNA System, 100rxn Promega E4871

Agencourt AMPure XP Beckman Coulter A63880

Native Barcoding Expansion 1-12 Nanopore EXP-NBD104

Native Barcoding Expansion 13-24 Nanopore EXP-NBD114

Native Barcoding Expansion Kit 1-96 Nanopore EXP-NBD196

Ligation Sequencing Kit Nanopore SQK-LSK109

Sequencing Auxiliary Vials Nanopore EXP-AUX001

Adapter Mix II Expansion Nanopore EXP-AMII001

Short Fragment Buffer Expansion kit Nanopore EXP-SFB001

Flow Cell Priming Kit Nanopore EXP-FLP002

R9.4.1 flow cells Nanopore FLO-MIN106

Troubleshooting

Safety warnings

- ! All procedures and manipulation with samples containing virulent *B. anthracis* spores should be performed in a biosafety level 3 laboratory.

Before start

Isolate DNA from environmental samples suspected to contain *Bacillus anthracis* spores with Qiagen DNeasy Blood&Tissue kit or similar suitable for DNA isolation from Gram positive bacteria.

It is recommended to apply an extra mechanical lysis step (for ex. bead beating) before DNA isolation to increase the effectiveness of spore disruption.






Before MLVA analysis check the isolated DNA with *Bacillus anthracis*-specific real-time PCR assay for *Bacillus anthracis* DNA content.



Primer pool preparation





- 1 If required resuspend lyophilised primers at a concentration of 100µM each. Primer names, characteristics, concentrations and volumes required for primer stocks are listed in the table below.

 BaMLVA_primers.xlsx

- 2 Generate  500 µL primer Pool 1 stock by adding  7 µL ,  13.5 µL or  15.5 µL of each primer to a  1.5 mL Eppendorf labelled "Pool 1 (stock)", following the table above.

Note

Primers should be diluted and pooled in the **mastermix cabinet** which should be cleaned with decontamination wipes and UV sterilised before and after use.

- 3 Dilute primer Pool 1 stock 1:10 in molecular grade water, to generate Pool 1 working solution.
- 4 Generate  100 µL primer Pool 2 stock by adding  5 µL of each odd region primer to a  1.5 mL Eppendorf labelled "Pool 2", and adjust final volume to  100 µL with molecular grade water.

Note













It is recommend that multiple aliquots of each primer pool are made to in case of degradation or contamination.

Multiplex PCR

12m 30s


- 5 In the mastermix hood set up the multiplex PCR reactions as follows in 0.2 mL 8-strip PCR tubes:

Component	Pool 1	Pool 2
-----------	--------	--------

5X Q5 Reaction Buffer	 5 µL	 5 µL
10 mM dNTPs	 0.5 µL	 0.5 µL
Q5 Hot Start DNA Polymerase	 0.25 µL	 0.25 µL
BaMLVA Primer Pool 1 working solution or Pool 2	 4.3 µL	 1.0 µL
Nuclease-free water	 12.45 µL	 15.75 µL
Total	 22.5 µL	 22.5 µL

Note

A PCR mastermix for each pool should be made up in the **mastermix cabinet** and aliquoted into PCR strip tubes. Tubes should be wiped down when entering and leaving the mastermix cabinet.

- 6 In the extraction and sample addition cabinet add  2.5 µL DNA to each tube and mix well by pipetting.

Note








The **extraction and sample addition cabinet** should be cleaned with decontamination wipes and UV sterilised before and after use.

- 7 Pulse centrifuge the tubes to collect the contents at the bottom of the tube.

- 8 Set-up the following programs on a gradient thermal cycler, or a thermal cycler suitable for running 2 or more different PCR cycles in one time, or 2 thermal cyclers:

12m 30s

Program for Pool 1 PCR:

Step	Temperature	Time	Cycles
Heat Activation	 98 °C	 00:00:30	1
Denaturation	 98 °C	 00:00:15	45
Annealing	 65 °C	 00:05:30	45
Hold	 4 °C	Indefinite	1

Program for Pool 2 PCR:

Step	Temperature	Time	Cycles
Heat Activation	🌡️ 98 °C	🕒 00:00:30	1
Denaturation	🌡️ 98 °C	🕒 00:00:15	45
Annealing	🌡️ 63 °C	🕒 00:05:30	45
Hold	🌡️ 4 °C	Indefinite	1

Equipment

Veriti 96-Well Thermal Cycler

NAME

Applied Biosystems


BRAND

4375786

SKU

<https://www.thermofisher.com/order/catalog/product/4375786#/4375786>^{LINK}

Quantification and normalisation

- Quantify  1 µL PCR product using the Quantus Fluorometer using the ONE dsDNA assay.

Protocol

NAME

DNA quantification using the Quantus fluorometer

CREATED BY

Josh Quick

Preview

Equipment

Quantus

NAME

Fluorometer

TYPE

Promega

BRAND


E6150

SKU

<https://www.promega.co.uk/products/microplate-readers-fluorometers-luminometers/fluorometers/quantus-fluorometer>






LINK

- 9.1 Remove Lambda DNA 400 ng/μL standard from the freezer and leave on ice to thaw. Remove ONE dsDNA dye solution from the fridge and allow to come to room temperature.

 QuantiFluor(R) ONE dsDNA System, 500rxn **Promega Catalog #E4870**



- 9.2 Set up two  0.5 mL tubes for the calibration and label them 'Blank' and 'Standard'



- 9.3 Add  200 μL ONE dsDNA Dye solution to each tube.
- 9.4 Mix the Lambda DNA standard 400 ng/ μL standard by pipetting then add  1 μL to one of the standard tube.
- 9.5 Mix each sample vigorously by vortexing for  00:00:05 and pulse centrifuge to collect the liquid.
- 9.6 Allow both tubes to incubate at room temperature for  00:02:00 before proceeding.
- 9.7 Selection 'Calibrate' then 'ONE DNA' then place the blank sample in the reader then select 'Read Blank'. Now place the standard in the reader and select 'Read Std'.
- 9.8 Set up the required number of  0.5 mL tubes for the number of DNA samples to be quantified.

Note


Use only thin-wall, clear, 0.5mL PCR tubes such as Axygen #PCR-05-C

- 9.9 Label the tubes on the lids, avoid marking the sides of the tube as this could interfere with the sample reading.
- 9.10 Add  199 μL ONE dsDNA dye solution to each tube.
- 9.11 Add  1 μL of each user sample to the appropriate tube.

Note

Use a P2 pipette for highest accuracy.




9.12 Mix each sample vigorously by vortexing for  00:00:05 and pulse centrifuge to collect the liquid.

9.13 Allow all tubes to incubate at room temperature for  00:02:00 before proceeding.

9.14 On the Home screen of the Quantus Fluorometer, select 'Protocol', then select 'ONE DNA' as the assay type.

Note



If you have already performed a calibration for the selected assay you can continue, there is no need to perform repeat calibrations when using ONE DNA pre diluted dye solution. If you want to use the previous calibration, skip to step 11. Otherwise, continue with step 9.



9.15 On the home screen navigate to 'Sample Volume' and set it to  1 μL then 'Units' and set it to ng/ μL .

9.16 Load the first sample into the reader and close the lid. The sample concentration is automatically read when you close the lid.


9.17 Repeat step 16 until all samples have been read.

9.18 The value displayed on the screen is the dsDNA concentration in ng/ μL , carefully record all results in a spreadsheet or laboratory notebook.

10 Label a  1.5 mL Eppendorf tube for each sample and assemble the following PCR dilution for each sample for final volume of  10 μL :




Pool 1 PCR reaction volume ( x μL) containing  196 ng PCR amplicon




Pool 2 PCR reaction volume ( x μL) containing  21 ng PCR amplicon

Nuclease-free water volume ( x μL) to a final volume of  10 μL


Total amount of PCR amplicons per sample  217 ng in  10 μL

**Note**

Input from Pool 1 and Pool 2 PCR reactions will vary depending on the starting amount of target DNA. If the Ct value of the target DNA is <30, it is possible to put the total  217 ng amount of PCR amplicons to  10 μL according to our experiences. If the Ct value of target DNA is >30, put as much PCR amplicons to  10 μL final volume as you can, keeping the 1:10 ratio of Pool 1 : Pool 2 PCR amplicons.

- 11 Dilute PCR amplicon pool of each sample 1:10 adding  90 μL molecular grade water, and mix well by pipetting.
- 12 Label a 0.2 mL PCR tube for each sample.
- 13 Use  10 μL input for the One-pot native barcoding reaction to give a total of  21.7 ng per sample.

Note

Input to the one-pot native barcoding reaction will vary depending on the amplicon length but we have determined  21.7 ng is the correct input for efficient barcoding of this amplicon length. Process at least 6 samples per native barcoded library in order to have sufficient material at the end.

Native barcoding

1h 49m

- 14 Barcode the amplicon pools using the one-pot native barcoding approach.

Note

We developed native barcoding protocol with modifications of Josh Quick 2020. One-pot native barcoding of amplicons v2. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bdp8i5rw>

- 14.1 Set up the following reaction for each sample:

**Component****Volume**

PCR dilution from previous step

10 µL

Nuclease-free water

2.5 µL

Ultra II End Prep Reaction Buffer

1.75 µL

Ultra II End Prep Enzyme Mix

0.75 µL

Total

15 µL

14.2 Incubate at room temperature for 00:10:00

21m

Incubate at 65 °C for 00:10:00

Incubate on ice for 00:01:00

14.3 In a new 0.2 mL PCR tube set up the following reaction:

Component**Volume**

Previous reaction mixture

3.5 µL

Nuclease-free water

3.7 µL

NBXX barcode

2.5 µL

Ultra II Ligation Master Mix

10 µL

Ligation Enhancer

0.3 µL

Total

20 µL

Note

Use one native barcode from the EXP-NBD104 (1-12) or EXP-NBD114 (13-24) or EXP-NBD196 (1-96) per sample. The minimum use of 6 barcodes is sufficient for effective application of R9 flow cells. Under these sample numbers the cost effectiveness of this method is highly decreased due to the low yield of overall data.

14.4 Incubate at room temperature for 00:20:00

36m

Incubate at 65 °C for 00:15:00

Incubate on ice for 00:01:00



**Note**

The 65°C incubation is to inactivate the DNA ligase to prevent barcode cross-ligation when reactions are pooled in the next step.

- 14.5 In a new 1.5 ml Eppendorf tube pool all  20 µL one-pot barcoding reactions together.

Note

It is recommended to pool maximum 6 one-pot barcoding reactions together in one tube. If more than 6 reactions are pooled into one tube, the next amplicon-cleaning step will take very long time due to the slow drying of high amount of SPRI beads.

- 14.6 Add 1.8x volume of SPRI beads to the sample tube and mix gently by either flicking or pipetting. For example add  216 µL SPRI beads to  120 µL 6-plex pooled one-pot native barcoding reactions.


Note

1.8x volume of SPRI will bind the shortest 200 bp amplicons in the presence of ligation buffer as in a one-pot reaction. It is recommended to use 1.8x volume of SPRI beads to not lose short amplicons even though this will result in excessive native barcode carryover.

Note





Vortex SPRI beads thoroughly before use to ensure they are well resuspended, the solution should be a homogenous brown colour.

- 14.7 Pulse centrifuge to collect all liquid at the bottom of the tube.



- 14.8 Incubate for  00:05:00 at room temperature.

5m




- 14.9 Place on magnetic rack and incubate for  00:02:00 or until the beads have pelleted and the supernatant is completely clear. 2m
- 14.10 Carefully remove and discard the supernatant, being careful not to touch the bead pellet.
- 14.11 Add  500 μL SFB and resuspend beads completely by pipette mixing.
- 14.12 Pulse centrifuge to collect all liquid at the bottom of the tube.
- 14.13 Place on magnetic rack and incubate until the beads have pelleted and the supernatant is completely clear.
- 14.14 Remove supernatant and discard.
- 14.15 Pulse centrifuge and remove any residual SFB.
- Note**
- You do not need to allow to air dry with SFB washes.
- 14.16 Bath the pellet in  500 μL of room-temperature 75% volume ethanol without resuspending the beads.
- 14.17 Carefully remove and discard ethanol, being careful not to touch the bead pellet.
- 14.18 Pulse centrifuge to collect all liquid at the bottom of the tube and carefully remove as much residual ethanol as possible using a P10 pipette.
- 14.19 With the tube lid open incubate for  00:05:00 or until the pellet loses it's shine (if the pellet dries completely it will crack and become difficult to resuspend). 5m



14.20 Resuspend pellet in  30 μL nuclease-free water, mix gently by either flicking or pipetting and incubate for  00:05:00 .

5m

14.21 Place on magnetic rack and transfer sample to a clean  1.5 mL Eppendorf tube ensuring no beads are transferred into this tube.

Note

If the barcoding reactions were pooled in 2 tubes (more than 6-plex), resuspend each pellet in 30-30 μL nuclease-free water, and pool into one tube the cleaned barcoded amplicon pools after incubation.

15 Set up the following AMII adapter ligation and clean-up with SFB.

Note

We developed adapter ligation protocol with modifications of Josh Quick 2020. Adapter ligation with AMII. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bdp9i5r6>

15.1 Set up the following AMII adapter ligation reaction:

Component**Volume**

Barcoded amplicon pool

 30 μL

NEBNext Quick Ligation Reaction Buffer (5X)



 10 μL

Adapter Mix (AMII)

 5 μL

Quick T4 DNA Ligase



 5 μL **Total** 50 μL **Note**

If the volume of barcoded amplicon pool is  60 μL , double each component of the adapter ligation reaction, thus final volume will be  100 μL .



15.2 Incubate at room temperature for  00:20:00 .


20m

15.3 Add 1x volume of SPRI beads to the sample tube (1:1 ratio of beads to sample volume) and mix gently by either flicking or pipetting. For example add  50 μ L SPRI beads to  50 μ L adapter ligation reaction.


Note

Vortex SPRI beads thoroughly before use to ensure they are well resuspended, the solution should be a homogenous brown colour.

15.4 Pulse centrifuge to collect all liquid at the bottom of the tube.

15.5 Incubate for  00:05:00 at room temperature.

5m

15.6 Place on magnetic rack and incubate for  00:02:00 or until the beads have pelleted and the supernatant is completely clear.

2m

15.7 Carefully remove and discard the supernatant, being careful not to touch the bead pellet.

15.8 Add  200 μ L SFB and resuspend beads completely by pipette mixing.

Note

SFB will remove excess adapter without damaging the adapter-protein complexes. Do not use 70% ethanol as in early clean-ups.

15.9 Pulse centrifuge to collect all liquid at the bottom of the tube.



15.10 Place on magnetic rack and incubate until the beads have pelleted and the supernatant is completely clear.


15.11 Remove supernatant and discard.

15.12 Repeat step 15.8-15.11. to perform a second SFB wash.

15.13 Pulse centrifuge and remove any residual SFB.

Note

You do not need to allow to air dry with SFB washes.

15.14 Add  15 µL EB and resuspend beads by pipette mixing.


15.15 Incubate at  37 °C for  00:08:00 .

8m

Note

The longer incubation at 37°C helps to eluate shorter amplicons.

15.16 Place on magnetic rack and transfer final library to a clean 1.5mL Eppendorf tube ensuring no beads are transferred into this tube.

16 Quantify  1 µL of the final library using the Quantus Fluorometer using the ONE dsDNA assay.



Protocol

NAME

DNA quantification using the Quantus fluorometer

CREATED BY

Josh Quick

[Preview](#)

Note

Final library can be now be stored at 4°C for up to a week if needed otherwise proceed directly to MinION sequencing.

MinION sequencing

- 17 Prime the flowcell and load 20-25 ng sequencing library onto the flowcell.

Protocol

NAME

Priming and loading a MinION flowcell

CREATED BY

Josh Quick

[Preview](#)

- 17.1 Thaw the following reagents at room temperature before placing on ice:



Sequencing buffer (SQB)

Loading beads (LB)

Flush buffer (FLB)




Flush tether (FLT)



- 17.2 Add  30 μL FLT to the FLB tube and mix well by vortexing.
- 17.3 If required place a new MinION flowcell onto the MinION by flipping open the lip and pushing one end of the flowcell under the clip and pushing down gently.
- 17.4 Rotate the inlet port cover clockwise by 90° so that the priming port is visible.
- 17.5 Take a P1000 pipette and tip and set the volume to  800 μL . Place the tip in the inlet port and holding perpendicularly to the plane of the flowcell remove any air from the inlet port by turning the volume dial anti-clockwise.


Note

Be careful not to remove so much volume that air is introduced onto the rectangular array via the outlet.


- 17.6 Load  800 μL of FLB (plus FLT) into the flow cell via the inlet port, dispense slowly and smoothly trying to avoid the introduction of any air bubbles.
- 17.7 Wait for  00:05:00 .
- 17.8 Gently lift the SpotON cover to open the SpotON port.
- 17.9 Load another  200 μL of FLB (plus FLT) into the flow cell via the inlet port, this will initiate a siphon at the SpotON port to allow you to load the library dilution.
- 17.10 In a new tube prepare the library dilution for sequencing:

Component**Volume**

SQB

 37.5 μL

LB

 25.5 μL


Final library

 12 μL

Total 75 µL**Note**

Mix LB immediately before use as they settle quickly.

Dilute library in EB if required.

- 17.11 Mix the prepared library gently by pipetting up and down just prior to loading.
- 17.12 Add the  75 µL library dilution to the flow cell via the SpotON sample port in a dropwise fashion. Ensure each drop siphons into the port before adding the next.
- 17.13 Gently replace the SpotON sample port cover, making sure the bung enters the SpotON port, close the inlet port and close the MinION lid.
- 18 Start the sequencing run using MinKNOW.
- 18.1 If required plug the MinION into the computer and wait for the MinION and flowcell to be detected.
- 18.2 Choose flow cell 'FLO-MIN106' from the drop-down menu.
- 18.3 Then select the flowcell so a tick appears.
- 18.4 Click the 'New Experiment' button in the bottom left of the screen.
- 18.5 On the New experiment popup screen, select the running parameters for your experiment from the individual tabs:

Experiment: Name the run in the experiment field, leave the sample field blank.



Kit: Selection: Select LSK109 and Native barcoding kit (EXP-NBD104 or EXP-NBD114 or EXP-NBD196).

Run Options: Set the run length to minimum 24 hours (you can stop the run once sufficient data has been collected as determined using RAMPART).

Basecalling: Leave basecalling turned but select 'superaccurate basecalling'.

Barcoding: Leave barcoding turned, turn on trim barcodes, but turn off Barcode both ends option.

Output: The number of files that MinKNOW will write to a single folder. By default this is set to 4000 but can be reduced to make RAMPART update more frequently.

Click 'Start run'.

Note

In case of using GridION or MinION with "high performance" computer, superaccurate basecalling and barcoding are recommended.

18.6 Monitor the progress of the run using the MinKNOW interface.