

May 29, 2025

Astrocyte Viral Transduction

 In 2 collections

DOI

<https://dx.doi.org/10.17504/protocols.io.n2bvjbwxpgk5/v1>

Camille Goldman¹

¹Icahn School of Medicine at Mount Sinai



Camille Goldman

Icahn School of Medicine at Mount Sinai, University of Oxford...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n2bvjbwxpgk5/v1>

Protocol Citation: Camille Goldman 2025. Astrocyte Viral Transduction. protocols.io

<https://dx.doi.org/10.17504/protocols.io.n2bvjbwxpgk5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 15, 2025

Last Modified: May 29, 2025

Protocol Integer ID: 218387

Keywords: ASAPCRN, astrocyte viral transduction protocol for transducing astrocyte, astrocyte viral transduction protocol, transducing astrocyte

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-024297

NASA

Grant ID: 80ARC022CA004

NIH/NINDS

Grant ID: R01NS114239

NIH/NINDS

Grant ID: UH3NS115064

NIH/NIA

Grant ID: T32AG04968

NIH/NINDS

Grant ID: F31NS13090

New York State Department of Health

Grant ID: NYSTEM-C32561GG

Abstract

Protocol for transducing astrocytes with lentivirus

- 1 Coate desired number of wells of a 6 well plate in 0.1% gelatin
- 2 Wash astrocytes with 1x PBS
- 3 Dissociate astrocytes in a trypsin enzyme solution (TrypLE Select, Gibco 12563011). For a 10cm plate of astrocytes, use 3mL per well
- 4 Transfer dissociated cells to a falcon tube with 10mL astrocyte media (AM; ScienCell #1801)
- 5 Centrifuge at 300rcf for 3 minutes at room temperature
- 6 Resuspend cells in 1mL of AM and count
- 7 Thaw lentivirus on ice
See <https://www.protocols.io/edit/lentiviral-production-dkdi4s4e> for lentivirus protocol
- 8 Transfer between 1×10^5 and 3×10^5 cells into a 1.5mL microfuge tube, 1 tube per viral transduction
- 9 Add lentivirus to the microfuge tube, so that the final concentration once cells are plated with the total volume of media will be 1:50 (40uL of lentivirus when final volume will be 2mL in a well of a 6-well plate)
NOTE: Viral titration may need to be determined on a batch to batch basis, depending on the lentivirus plasmid and the cells being transduced.
- 10 Incubate cells in the microfuge tube with the higher lentivirus concentration for 5-10 minutes at room temperature
- 11 Add the cell/lentivirus mixture to the prepared well and bring the final volume of AM to 2mL
- 12 After 24hrs, change to fresh AM media



- 13 If antibiotic selection is necessary, start adding at least 3 days after transduction, and continue selecting for at least 2 days