

Dec 28, 2023

Baermann Funnel Extraction

DOI

dx.doi.org/10.17504/protocols.io.n2bvj31rxlk5/v1

Casey A Schlenker¹

¹The Ohio State University

OSU Soybean Pathology...



Casey A Schlenker

The Ohio State University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n2bvj31rxlk5/v1>

Protocol Citation: Casey A Schlenker 2023. Baermann Funnel Extraction. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.n2bvj31rxlk5/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: December 27, 2023

Last Modified: February 27, 2024

Protocol Integer ID: 92763

Keywords: baermann funnel extraction, extraction

Abstract

Baermann Funnel Extraction

Materials

Materials for sample preparation:





- Glass funnel
- Funnel stand
- Filter paper
- Tubing
- 2mm mesh sieve
- DI water bottles

Troubleshooting



Sample Preparation

6d

- 1 Sieve soil to 2mm
- 2 Weigh out  50 g soil in a weigh boat
- 3 Place tissue paper on top of Baermann funnel apparatus
- 4 Place weighed and sieved soil on top of filter paper
- 5 Fill Baermann funnel with DI water using a DI water bottle so that soil sample is saturated but not submerged in water
- 6 Cover all Baermann funnels with a plastic Uline
- 7 At  24:00:00 and  48:00:00 mark, check soil and refill any funnels with decreased water levels to prevent the soil from drying out. 2d
- 8 At  72:00:00 empty funnel over a 25um mesh sieve to collect nematodes 1d
- 9 Transfer the nematodes that are on the sieve into a labelled centrifuge tube using a **hot fixative**

Protocol

NAME

Fixative preparation for Baermann funnel extraction

CREATED BY

Casey A Schlenker

[Preview](#)



- 10 Collect soil from top of Baermann funnel; dry and weigh to get gravimetric soil moisture content
- 11 Count the entire sample at 50x magnification using a dissecting scope
- 12 Identify 100 nematodes