



Aug 17, 2023

## BARseq - BARseq-styled in situ sequencing for barcoded rabies virus

DOI

<https://dx.doi.org/10.17504/protocols.io.n2bvj82q5gk5/v1>

Xiaoyin Chen<sup>1</sup>, Mara CP Rue<sup>1</sup>

<sup>1</sup>Allen Institute for Brain Science



Xiaoyin Chen

Allen Institute for Brain Science

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.n2bvj82q5gk5/v1>

**Protocol Citation:** Xiaoyin Chen, Mara CP Rue 2023. BARseq - BARseq-styled in situ sequencing for barcoded rabies virus. protocols.io <https://dx.doi.org/10.17504/protocols.io.n2bvj82q5gk5/v1>

**Manuscript citation:**

<https://doi.org/10.7554/eLife.87866.1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



**Protocol status:** Working

**We use this protocol and it's working**

**Created:** January 13, 2023

**Last Modified:** August 17, 2023

**Protocol Integer ID:** 75215

**Keywords:** BARseq, in situ sequencing, spatial transcriptomics, barcoding-based neuroanatomy, barcoded rabies virus, rabies virus, rabies virus this protocol, application of barseq, barseq, situ sequencing, style in situ sequencing, tracing experiment

**Funders Acknowledgements:**

NIH

Grant ID: 1DP2MH132940

## Abstract

This protocol describes the application of BARseq-style in situ sequencing adapted for barcoded rabies virus. Similar procedures for both trans-synaptic tracing and retrograde tracing experiments.

## Guidelines

Standard precautions with RNA samples should be taken to reduce RNA degradation during tissue processing and library preparation. Pipetting and suctioning should be gentle throughout the whole procedure, and sample should not be left dried.

## Materials

### DNA oligos

|  | A          | B                                                                                                         | C                                                         |
|--|------------|-----------------------------------------------------------------------------------------------------------|-----------------------------------------------------------|
|  | YS220      | GATCGTCGGACTGTAGAACTCTGAACCTGTCTG                                                                         | sequencing primer                                         |
|  | YS221      | /5Alex594N/GATCGTCGGACTGTAGAACTCTGAACCTGTCTG                                                              | Hybridization probe to visualize all genes                |
|  | XC1417     | GTTTCAGAGTTCTACAGTCCGACGATC                                                                               | RCA primer for SNAP25 RNA padlock                         |
|  | XC275<br>7 | /5AmMC12/NNNNNNNNNNNNNNNNNNNNNNNN                                                                         | N20 for random priming                                    |
|  | XC275<br>8 | /5Alex488N/AGTCAGCGTCGAGCACGCGGCACTTATTGCA                                                                | Hybridization probe to visualize Slc17a7                  |
|  | XC275<br>9 | /5Alex532N/TGAGTAGAGTTGACTAAGAGCCGTTAGATGCC                                                               | Hybridization probe to visualize Gad1                     |
|  | XC276<br>0 | /5Alex647N/TCGCTGTACTAATAGTTGTCGACAGATCGTCA                                                               | Hybridization probe to visualize B19G                     |
|  | XCAI5      | /5phos/agctccggcattttgttattcaTCCTCTATGATTACTGACTGCGTCTATTTAGTGGAGCCATTGCTATCTTCTTggatatacacatccgtagattgct | pRV-4mCherry-NheI-N20 BC gapfilling padlock               |
|  | XCAI6      | /5AmMC6/g+ac+at+at+tc+ga+gt+gactcataagaagt                                                                | primer 1 for pRV-4mCherry-NheI-N20 and pRV-4mCherry-CSCS2 |
|  | XCAI7      | /5AmMC6/g+aa+gt+tg+aa+ta+ac+aaaatgccggagc                                                                 | primer 2 for pRV-4mCherry-NheI-N20 and pRV-4mCherry-CSCS2 |
|  | XCAI6<br>3 | /5AmMC6/gggagtgactgacacctccctccctg                                                                        | RV B19G RT primer, Hyb XC2760 for detection               |
|  | XCAI6<br>4 | /5AmMC6/tatcaacatcaaggcagtcagggccc                                                                        | RV B19G RT primer, Hyb XC2760 for detection               |
|  | XCAI6<br>5 | /5AmMC6/tcctgagacctgattgtgcacatcgg                                                                        | RV B19G RT primer, Hyb XC2760 for detection               |
|  | XCAI6<br>6 | /5AmMC6/aactccatatgttgctggaggaggga                                                                        | RV B19G RT primer, Hyb XC2760 for detection               |
|  | XCAI6<br>7 | /5AmMC6/tgcttttccaaacccagggacaagtt                                                                        | RV B19G RT primer, Hyb XC2760 for detection               |
|  | XCAI6<br>8 | /5AmMC6/ttcgtctgagcgaaagtcgtgcagg                                                                         | RV B19G RT primer, Hyb XC2760 for detection               |
|  | XCAI6<br>9 | /5AmMC6/ttgcacgagacccatgttccatcca                                                                         | RV B19G RT primer, Hyb XC2760 for detection               |



|  | A      | B                                                                                              | C                                                        |
|--|--------|------------------------------------------------------------------------------------------------|----------------------------------------------------------|
|  | XCAI70 | /5AmMC6/aggcctctttcatctacaaagccgca                                                             | RV B19G RT primer, Hyb XC2760 for detection              |
|  | XCAI71 | /5AmMC6/acaaccctagctcgcgattctcgggc                                                             | RV B19G RT primer, Hyb XC2760 for detection              |
|  | XCAI72 | /5AmMC6/aagacaccgctactcctgagcacttc                                                             | RV B19G RT primer, Hyb XC2760 for detection              |
|  | XCAI73 | /5AmMC6/tggtttttacagttcgaagccagcgg                                                             | RV B19G RT primer, Hyb XC2760 for detection              |
|  | XCAI74 | /5AmMC6/caccggccatcttcagttgtacgcg                                                              | RV B19G RT primer, Hyb XC2760 for detection              |
|  | XCAI75 | /5AmMC6/agtgtaggttcagcctccgtcaciaa                                                             | RV B19G RT primer, Hyb XC2760 for detection              |
|  | XCAI76 | /5AmMC6/atgtaggagaaccctgacaggttggt                                                             | RV B19G RT primer, Hyb XC2760 for detection              |
|  | XCAI77 | /5phos/acacaatctcagagggacaggTCGCTGTACTAATAGTTGTCGACAGATCGTCACTTCGTTCCCTcaatcgatcagaa cctacgca  | RV B19G padlock, Hyb XC2760 for detection, use with XC63 |
|  | XCAI78 | /5phos/gggaagtatgtattactgagtgTCGCTGTACTAATAGTTGTCGACAGATCGTCACTTCGTTCCCTgacttgggtctcc cgaactgg | RV B19G padlock, Hyb XC2760 for detection, use with XC64 |
|  | XCAI79 | /5phos/tggtgaagttcaccttcccgaTCGCTGTACTAATAGTTGTCGACAGATCGTCACTTCGTTCCCTcggtgacgaggctga ggattt  | RV B19G padlock, Hyb XC2760 for detection, use with XC65 |
|  | XCAI80 | /5phos/tcccagagatgcaatcatcccTCGCTGTACTAATAGTTGTCGACAGATCGTCACTTCGTTCCCTgacctgacggcaat gtcttaa  | RV B19G padlock, Hyb XC2760 for detection, use with XC66 |
|  | XCAI81 | /5phos/gtttcagacgtctcagtcattTCGCTGTACTAATAGTTGTCGACAGATCGTCACTTCGTTCCCTtcatgacaaccaagtc agtga  | RV B19G padlock, Hyb XC2760 for detection, use with XC67 |
|  | XCAI82 | /5phos/cccgataagttggtgaacctgTCGCTGTACTAATAGTTGTCGACAGATCGTCACTTCGTTCCCTaatgaaaccaaagtg tgcctt  | RV B19G padlock, Hyb XC2760 for detection, use with XC68 |
|  | XCAI83 | /5phos/tggagttctaggacttagactTCGCTGTACTAATAGTTGTCGACAGATCGTCACTTCGTTCCCTgagcatgcaaactc aagttatg | RV B19G padlock, Hyb XC2760 for detection, use with XC69 |
|  | XCAI84 | /5phos/ccaaagggagtgagacttgcgTCGCTGTACTAATAGTTGTCGACAGATCGTCACTTCGTTCCCTatagtagaggga gagagcat   | RV B19G padlock, Hyb XC2760 for detection, use with XC70 |
|  | XCAI85 | /5phos/gattacaccatttgatgcccTCGCTGTACTAATAGTTGTCGACAGATCGTCACTTCGTTCCCTacctactgtccacta accac    | RV B19G padlock, Hyb XC2760 for detection, use with XC71 |

|  | A           | B                                                                                                       | C                                                              |
|--|-------------|---------------------------------------------------------------------------------------------------------|----------------------------------------------------------------|
|  | XCAI8<br>6  | /5phos/aggggtcttccttagcggggaagtTCGCTGTACTAATAG<br>TTGTCGACAGATCGTCACTTCGTTCC Ttatgacagatccct<br>tcactcg | RV B19G padlock, Hyb<br>XC2760 for detection, use<br>with XC72 |
|  | XCAI8<br>7  | /5phos/caatccgtaccctgactaccgTCGCTGTACTAATAGT<br>TGTCGACAGATCGTCACTTCGTTCC Tcccagatatgaaga<br>gtctctaca  | RV B19G padlock, Hyb<br>XC2760 for detection, use<br>with XC73 |
|  | XCAI8<br>8  | /5phos/ccagatgcatgtagagccgcgtTCGCTGTACTAATAG<br>TTGTCGACAGATCGTCACTTCGTTCC Tagaaagcatttcc<br>gccaaca    | RV B19G padlock, Hyb<br>XC2760 for detection, use<br>with XC74 |
|  | XCAI8<br>9  | /5phos/gttcacttgcacaggcggttgtTCGCTGTACTAATAGTT<br>GTCGACAGATCGTCACTTCGTTCC Ttcttagccataaaaagtg<br>aacgg | RV B19G padlock, Hyb<br>XC2760 for detection, use<br>with XC75 |
|  | XCAI9<br>0  | /5phos/tggaggacgaaggatgcaccaTCGCTGTACTAATAG<br>TTGTCGACAGATCGTCACTTCGTTCC Tgctgccaacaaca<br>tttgtag     | RV B19G padlock, Hyb<br>XC2760 for detection, use<br>with XC76 |
|  | XCAI13<br>1 | TGGAGCCATTGCTATCTTCTTggatatacacaatccgtagatt<br>gct                                                      | sequencing primer for<br>XCAI5                                 |

## Reagents

☒ Phusion high-fidelity PCR kit **Thermo Scientific Catalog #F553S**

☒ Tween-20 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P-7949**

BS(PEG)9, 100 mg (Note: BS(PEG)9 loses its effectiveness 1 month after reconstitution in DMSO. Prepare a  
☒ fresh batch every month, especially if it has been frozen and thawed repeatedly. **Thermo Scientific Catalog #21582**

☒ Formamide **Thermo Fisher Scientific Catalog #AM9342**

☒ 10x PBS **Thermo Fisher Scientific Catalog #AM9624**

☒ RNase-free water

☒ dNTP Mix (dATP, dCTP, dGTP, and dTTP, each at 10mM) **Thermo Fisher Scientific Catalog #R0192**

☒ phi29 DNA Polymerase (10 U/μL) **Thermo Scientific Catalog #EP0091**

☒ RNase H **Enzymatics Catalog #Y9220L**

☒ Glycerol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #G5516**

⊗ Ethanol **Merck Millipore (EMD Millipore) Catalog #100983**

⊗ Pierce&trade; MMTS (methyl methanethiosulfonate) **Thermo Fisher Catalog #23011**

⊗ SSC (20X), RNase-free **Thermo Fisher Catalog #AM9770**

⊗ RiboLock RNase Inhibitor (40 U/μL) **Thermo Fisher Catalog #EO0381**

⊗ RevertAid H Minus Reverse Transcriptase (200 U/μL) **Thermo Fisher Catalog #EP0452**

⊗ Paraformaldehyde 20% **Electron Microscopy Sciences Catalog #15713**

⊗ BSA Molecular grade **New England Biolabs Catalog #B9000S**

⊗ Ampligase DNA Ligase Kit **Lucigen Catalog #A8101**

⊗ KCl (2 M) RNase-free **Thermo Fisher Scientific Catalog #AM9640G**

⊗ Aminoallyl-dUTP Solution (50 mM) **Thermo Fisher Scientific Catalog #R1101**

⊗ Tris (1 M) pH 8.0 RNase-free **Thermo Fisher Scientific Catalog #AM9855G**

⊗ HiSeq SBS Kit v4 **Illumina, Inc. Catalog #FC-401-4003**

⊗ Grace Bio-Labs HybriWell-FL™ sealing system Fluor-friendly adhesive chamber **Merck MilliporeSigma (Sigma-Aldrich) Catalog #GBL612204**

Other equipment required include incubators set at 37 °C, 45 °C, and 60 °C. All tubes should be RNase-free. RNase-free filter tips should be used. A Crest Xlight v3 spinning disk confocal on a Nikon Ti2E with Photometrics Kinetix, and Lumencor Celesta was used for imaging the sequencing steps. The filters and lasers used are indicated in Table 1.

|  | A            | B     | C                                  | D                        |
|--|--------------|-------|------------------------------------|--------------------------|
|  | Channel<br>s | Laser | Dichroic                           | Emission filter          |
|  | G/YFP        | 514   | Zt405/514/6<br>35rpc               | FF01-565/24              |
|  | T/RFP        | 561   | FF421/491/5<br>67/659/776<br>-Di01 | FF01-441/511/593/684/817 |
|  | A            | 640   | Zt405/514/6<br>35rpc               | FF01-676/29              |
|  | C            | 640   | Zt405/514/6<br>35rpc               | FF01-775/140             |
|  | GFP          | 488   | FF421/491/5<br>72-Di01             | 69401m                   |
|  | DAPI         | 405   | FF421/491/5<br>72-Di01             | 69401m                   |
|  | TexasRe<br>d | 561   | FF421/491/5<br>72-Di01             | 69401m                   |

|  | A   | B   | C                    | D           |
|--|-----|-----|----------------------|-------------|
|  | Cy5 | 640 | Zt405/514/6<br>35rpc | ZET532/640m |

Table 1. Laser and filter settings for sequencing imaging.

## Safety warnings

 Use caution when handling liquids containing formaldehyde and formamide.

## Library preparation

- 1 Tissues with barcoded neurons should be cryo-sectioned to 20  $\mu\text{m}$  and mounted on slides. Slides can be stored at  $-80\text{ }^{\circ}\text{C}$  for up to a month.
- 2 **DAY 1**  
  
Take slide(s) out of  $-80\text{ }^{\circ}\text{C}$  and immerse immediately in 4% paraformaldehyde in 1x PBS (2 slides per 50mL falcon tube, back-to-back)
- 3 Incubate for 1 hour at room temperature on slow shaker
- 4 Wash the slides by immersing in 1x PBS (2 slides per 50ml falcon tube, back to back)
- 5 Wipe excess PBS off the surface of the chamber, then stick on the Hybriwell-FL chambers. Note that the ports on the chamber should be placed as far away from the tissue slices as possible.
- 6 Wash twice in PBST (1x PBS + 0.5% Tween-20)
- 7 Wash in 70% Ethanol for 5 mins
- 8 Wash in 85% Ethanol for 5 mins
- 9 Wash in 100% Ethanol for 5 mins
- 10 Replace with new 100% Ethanol, drop extra 100% Ethanol on top of slides and cover with ParaFilm to avoid evaporation. Incubate for at least 1.5 hrs at  $4\text{ }^{\circ}\text{C}$  (up to 3 hours)
- 11 Wash in PBST for 4-6 times, until all bubbles are cleared in the chamber and PBST flows in and out of the chamber smoothly.
- 12 Make reverse transcription mix: 50  $\mu\text{M}$  N20 primer (XC2757), 2  $\mu\text{M}$  XCAI6, 2  $\mu\text{M}$  XCAI7, 20 U/ $\mu\text{L}$  RevertAid H Minus M-MuLV reverse transcriptase, 500  $\mu\text{M}$  dNTP, 0.2  $\mu\text{g}/\mu\text{L}$  BSA, 1 U/ $\mu\text{L}$  RiboLock RNase Inhibitor, 1x RevertAid RT buffer.

For the monosynaptic tracing experiments, the RT primers additionally included 2  $\mu$ M of primers for the rabies glycoprotein (XCAI63 through XCAI76).

- 13 Incubate in reverse transcription mix overnight at 37 °C. Create a humidity chamber to avoid the slides drying out using kim-wipes and DI water.

- 14 **DAY 2:**

Wash with PBST once

- 15 Incubate in a mixture of 1 $\mu$ L BS(PEG)9 per 4  $\mu$ L PBST (e.g. 200ul BS(PEG)9 and 800ul PBST) for one hour at room temperature

- 16 Wash with 1M Tris pH 8.0, then incubate in new 1M Tris pH 8.0 for 30 mins

- 17 Wash twice in PBST

- 18 Make non-gap-filling ligation mix: 1x Ampligase buffer, 20 nM padlock probe each, 0.5 U/ $\mu$ L Ampligase, 0.4 U/ $\mu$ L RNase H, 1 U/ $\mu$ L RiboLock RNase Inhibitor, 50 mM KCl (extra of those already provided by the ampligase buffer), 20% formamide.

In the retrograde tracing experiments, the non-gap-filling padlock probe mix included all padlock probes for endogenous genes. In the monosynaptic tracing experiments, the non-gap-filling padlock probe mix included all padlock probes for endogenous genes except *Slc30a3*, and additionally included padlocks for the rabies glycoprotein (XCAI77 – XCAI90).

- 19 Incubate in ligation mix for at least 30 mins at 37 °C (can go longer but not shorter), then at least 45 mins at 45 °C (can go longer but not shorter).

- 20 Make the gap-filling ligation mix [same as the non-gap-filling mix with the rabies barcode padlock probe (XCAI5) as the only padlock probe, and with 50  $\mu$ M dNTP, 0.2 U/ $\mu$ L Phusion DNA polymerase, and 5% glycerol]

1x Ampligase buffer, 20 nM padlock probe each, 0.5 U/ $\mu$ L Ampligase, 0.4 U/ $\mu$ L RNase H, 1 U/ $\mu$ L RiboLock RNase Inhibitor, 50 mM KCl (extra of those already provided by the ampligase buffer), 20% formamide, 50  $\mu$ M dNTP, 0.2 U/ $\mu$ L Phusion DNA polymerase, and 5% glycerol.

- 21 Wash twice in PBST



- 22 Incubate in ligation mix for 5 mins at 37 °C, then 45 mins at 45 °C. **\*\*exact timing on this step!\*\***
- 23 Wash twice in PBST, then once in FISH Wash (2x SSC with 10% formamide)
- 24 Hybridize with 1 µM RCA primer (XC1417) in FISH wash for 10 mins at room temperature
- 25 Wash twice in FISH wash, then twice in PBST
- 26 Make RCA mix: 1 U/µL phi29 DNA polymerase, 1x phi29 polymerase buffer, 0.25 mM dNTP, 0.2 µg/µL BSA, 5% glycerol (extra of those from the enzymes), 125 µM aminoallyl dUTP
- 27 Incubate in RCA mix overnight at room temperature
- 28 **DAY 3:**  
  
Wash with PBST once
- 29 Incubate in a mixture of 1µL BS(PEG)9 per 4 µL PBST (e.g. 200ul BS(PEG)9 and 800ul PBST) for one hour at room temperature
- 30 Wash with 1M Tris pH 8.0, then incubate in new 1M Tris pH 8.0 for 30 mins
- 31 Wash twice in PBST

## Sequencing

- 32 **Hybridization of Gene sequencing primer:**  
  
Wash with FISH wash (2x SSC with 10% formamide)
- 33 Hybridize sequencing primer (YS220) with a primer concentration of 1 µM in FISH wash for 10 mins at room temperature



34 Wash with FISH wash three times, 2 mins each

35 Wash with PBST twice

36 **Sequence first cycle (genes and barcodes):**

Do the following incubations. Unless noted with incubation temperature, each step is performed at room temperature. For steps without incubation time, treat these as quick washes. large flat metal blocks can be used to place sample slides to quickly cycle through high and low temperatures. This version uses MiSeq Nano v2 kit:

<https://www.illumina.com/products/by-type/sequencing-kits/cluster-gen-sequencing-reagents/miseq-reagent-kit-v2.html>

36.1 Incorporation Buffer 60 °C 3 mins x1

36.2 2% PBST x1

36.3 Idoacetamide blocker: For 9.3mg vial dilute pellet in between 2.5-3.5mL 2% PBST. Make fresh tube daily and store out of light.

Idoacetamide blocker 60 °C 3 mins x1

36.4 2% PBST x1

36.5 Incorporation Buffer x2

36.6 IRM 60 °C 3 mins x2

36.7 2% PBST x1

36.8 2% PBST 60 °C 3 mins x4

36.9 Replace PBST with USM and Image \*\*if slides are dirty, clean with 70% Ethanol before adding USM\*\*

**37 Sequence subsequent cycles (genes and barcodes):**

Do the following incubations. Unless noted with incubation temperature, each step is performed at room temperature. For steps without incubation time, treat these as quick washes. large flat metal blocks can be used to place sample slides to quickly cycle through high and low temperatures.

37.1 Incorporation buffer x2

37.2 CRM 60 °C 3 mins x2

37.3 Incorporation buffer x1 - Wipe ports after adding the incorporation buffer, to ensure that no CRM is left on the slide's surface

37.4 2% PBST x1

37.5 Idoacetamide blocker: For 9.3mg vial dilute pellet in between 2.5-3.5mL 2% PBST. Make fresh tube daily and store out of light.

Idoacetamide blocker 60 °C 3 mins x1

37.6 2% PBST x1

37.7 Incorporation buffer x2

37.8 IRM 60 °C 3 mins x2

37.9 2% PBST x1

37.10 2% PBST 60 °C 3 mins x4

37.11 Replace 2% PBST with USM and image \*\*if slides are dirty, clean with 70% Ethanol before adding USM\*\*

## 38 After completing all gene sequencing imaging cycles:

### Hybridization cycle

#### 38.1 Hybridize probes:

Make strip buffer: 60% formamide 2xSSC 0.01% Tween20

Strip buffer 60 °C 5 mins x3

Cool down quickly on metal plates between washes, place on metal plates in 60 °C oven to heat up quickly.

#### 38.2 FISH wash (2x SSC with 10% formamide) 1x

#### 38.3 Hybridize probes (YS221, XC2758, XC2759, XC2760) with a primer concentration of 1 µM in FISH wash at 60 °C for 2 minutes, then for 10 mins at room temperature. Rotate plates in holder to ensure they cool down slowly.

#### 38.4 FISH wash x1

#### 38.5 0.002 mg/ML DAPI in PBST, room temperature for 5 mins

#### 38.6 Replace PBST with USM and image \*\*if slides are dirty, clean with 70% Ethanol before adding USM\*\*

## 39 After completing all gene sequencing imaging cycles and hybridization cycle:

### Hybridize barcode sequencing primers

#### 39.1 Strip buffer (60% formamide 2xSSC 0.01% Tween20)

60 °C 5 mins x3

Cool down quickly on metal plates between washes, place on metal plates in 60 °C oven to heat up quickly.

#### 39.2 FISH wash (2x SSC with 10% formamide) 1x

#### 39.3 Hybridize probe XCA131 with a primer concentration of 1 µM in 2x SSC and FISH wash at 60 °C for 2 minutes, then for 10 mins at room temperature. Rotate plates in holder to ensure they cool down slowly.

#### 39.4 FISH wash x1



- 40 Sequencing barcodes: The sequencing procedures are the same as those for gene sequencing (Steps 36 and 37)

After hybridizing barcode primers, go directly back to step 36, sequence first barcode imaging cycle. Then repeat step 37 for desired length of barcode sequence.