

Nov 16, 2023

## Cell isolation from dorsal mouse skin for single-cell RNA-seq

DOI

[dx.doi.org/10.17504/protocols.io.yxmvmn8m9g3p/v1](https://dx.doi.org/10.17504/protocols.io.yxmvmn8m9g3p/v1)

Cenk Celik<sup>1</sup>, Stella Yue Ting Lee<sup>1</sup>, Guillaume Thibault<sup>1</sup>

<sup>1</sup>Nanyang Technological University



Cenk Celik

Nanyang Technological University

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.yxmvmn8m9g3p/v1>

**Protocol Citation:** Cenk Celik, Stella Yue Ting Lee, Guillaume Thibault 2023. Cell isolation from dorsal mouse skin for single-cell RNA-seq. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.yxmvmn8m9g3p/v1>

**Manuscript citation:**

Cenk Celik, Stella Yue Ting Lee, Frederick Reinhart Tanoto, Mark Veleba, Kimberly Kline, Guillaume Thibault (2024) Decoding the complexity of delayed wound healing following *Enterococcus faecalis* infection eLife 13:RP95113

<https://doi.org/10.7554/eLife.95113.3>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** August 08, 2022

**Last Modified:** November 16, 2023

**Protocol Integer ID:** 68336

**Keywords:** enterococcus faecalis infection, cell isolation from dorsal mouse skin, heterogeneity of the skin cell population, enterococcus, skin cell population, faecalis infection, dorsal mouse skin, cell isolation, cell rna, rna

**Funders Acknowledgements:**

NMRC

Grant ID: MOH-000566

## Abstract

We identified heterogeneity of the skin cell populations upon *Enterococcus faecalis* infection.



## Materials

### Consumables

- Hank's Balanced Salt Solution  $\text{Ca}^{2+}/\text{Mg}^{2+}$ -free
- Dulbecco's Phosphate Buffered Saline  $\text{Ca}^{2+}/\text{Mg}^{2+}$ -free
- Dispase
- Collagenase type I
- Liberase TM Research Grade
- 0.25% Trypsin/EDTA
- Bovine Serum Albumin
- Trypan Blue
- Ice

### Instruments

- $\text{CO}_2$  supplemented 37°C incubator
- Chromium Instrument/X, 10X Genomics
- Countess 3 Automated Cell Counter, Invitrogen
- 15ml and 50ml tubes-compatible centrifuge
- 1.5ml tube-compatible centrifuge

### Surgical tools

- 6-mm punch
- Scissors
- Forceps
- Bistouri
- Blades G23
- 10mL syringes
- 21G needles

### Other labware

- 50ml centrifuge tubes
- 15ml centrifuge tubes
- 1.5ml Protein LoBind tubes
- Disposable haemocytometer (4 per sample)
- 10mm dishes
- 6-well plates
- Pasteur pipettes



## Protocol materials

- ✕ Dispace II, powder **Thermo Fisher Catalog #17105041**
- ✕ 1X Dulbecco's Phosphate Buffered Saline (DPBS) **Thermo Fisher Scientific Catalog #14190094**
- ✕ Collagenase Type I powder **Thermo Fisher Scientific Catalog #17100017**
- ✕ Hanks Balanced Salt Solution (HBSS) without Ca<sup>2</sup> Mg<sup>2</sup> **Thermo Fisher Scientific Catalog #88284**
- ✕ Hanks Balanced Salt Solution (HBSS) without Ca<sup>2</sup> Mg<sup>2</sup> **Thermo Fisher Scientific Catalog #88284**
- ✕ Liberase TM **Merck MilliporeSigma (Sigma-Aldrich) Catalog #000000005401119001**
- ✕ Hanks Balanced Salt Solution (HBSS) without Ca<sup>2</sup> Mg<sup>2</sup> **Thermo Fisher Scientific Catalog #88284**
- ✕ Bovine Serum Albumin **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A7030**
- ✕ Hanks Balanced Salt Solution (HBSS) without Ca<sup>2</sup> Mg<sup>2</sup> **Thermo Fisher Scientific Catalog #88284**
- ✕ 1X Dulbecco's Phosphate Buffered Saline (DPBS) **Thermo Fisher Scientific Catalog #14190094**
- ✕ Hanks Balanced Salt Solution (HBSS) without Ca<sup>2</sup> Mg<sup>2</sup> **Thermo Fisher Scientific Catalog #88284**
- ✕ Trypsin EDTA **Gibco - Thermo Fisher Scientific Catalog #25-051-Cl.**
- ✕ Hanks Balanced Salt Solution (HBSS) without Ca<sup>2</sup> Mg<sup>2</sup> **Thermo Fisher Scientific Catalog #88284**
- ✕ Standard Biopsy Punches, Disposable Standard biopsy punch; 6mm **Thermo Fisher Catalog #12460412**
- ✕ scalpel blades
- ✕ 10 mL syringes **Becton Dickinson (BD) Catalog #BD 309695**
- ✕ Falcon® 70 µm cell strainer **Corning Catalog #352350**
- ✕ 10 mL syringes **Becton Dickinson (BD) Catalog #BD 309695**
- ✕ Falcon® 70 µm cell strainer **Corning Catalog #352350**
- ✕ Falcon® 40 µm Cell Strainer **Corning Catalog #352340**

## Troubleshooting




## Enzyme solutions

20m

### 1 Dispace Solution

2m

Dissolve 0.5g of  Dispace II, powder **Thermo Fisher Catalog #17105041** in a 50ml of

 1X Dulbecco's Phosphate Buffered Saline (DPBS) **Thermo Fisher Scientific Catalog #14190094**

to make [M] 10 mg/mL of Dispace II solution

0.2 Filter the solution by using a 0.2µm filter

2m

0.3 Aliquot 1.25mL in tubes and store at  -20 °C


5m

### 1 Collagenase type I Solution

2m

Dissolve 1g of

 Collagenase Type I powder **Thermo Fisher Scientific Catalog #17100017** in 4mL of

 Hanks Balanced Salt Solution (HBSS) without Ca<sup>2+</sup> Mg<sup>2+</sup> **Thermo Fisher Scientific Catalog #88284**

to make a [M] 250 mg/mL collagenase type I solution

1.1 Make a stock solution

2m

i.e.  $67500\text{U/mL} \times 4\text{mL} = 50\text{U}/\mu\text{L} \times V_{\text{max}}$

$V_{\text{max}} = 5.4\text{mL}$

$5.4 - 4.0 = 1.4\text{mL}$  to add the stock to make the stock solution

1.2 Filter the solution by using a 0.2µm filter

2m

1.3 Aliquot 20µL in tubes and store at  -20 °C

5m

### 2 Liberase TM Solution

2m



Punch through the cap with a syringe to add 2mL of



Hanks Balanced Salt Solution (HBSS) without Ca<sup>2+</sup> Mg<sup>2+</sup> **Thermo Fisher Scientific Catalog #88284**

into 5mg



Liberase TM **Merck MilliporeSigma (Sigma-Aldrich) Catalog #00000000540119001**

bottle to make 1mL 2.5 mg/mL liberase solution

2.1 Dilute into 1mL 0.2 mg/mL in

2m



Hanks Balanced Salt Solution (HBSS) without Ca<sup>2+</sup> Mg<sup>2+</sup> **Thermo Fisher Scientific Catalog #88284**

(1.04 Wünsch unit/mL)

2.2 Filter the solution by using a 0.2µm filter

2m

2.3 Aliquot 750µL in tubes and store at -20 °C

5m

## Working Solutions

2m

### 3 Working Solution 1 (WS1)

2m

Make a 1mL 0.5 Mass / % volume of



Bovine Serum Albumin **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A7030** in

50mL of



Hanks Balanced Salt Solution (HBSS) without Ca<sup>2+</sup> Mg<sup>2+</sup> **Thermo Fisher Scientific Catalog #88284**

3.1 Filter the solution by using a 0.2µm filter and keep it at 4 °C

2m

### 4 Working Solution 2 (WS2)

2m

Make a 1mL 0.04 Mass / % volume of Bovine serum albumin solution in 50mL of



1X Dulbecco's Phosphate Buffered Saline (DPBS) **Thermo Fisher Scientific Catalog #14190094**



4.1 Filter the solution by using a 0.2µm filter and keep it at 4 °C

2m

## 5 Enzyme Cocktail 1 (EC1)–4 samples

3m

From the enzyme stocks that have been prepared before, mix:

- 1.25mL of Dispase
- 0.75mL of Liberase
- 20µL of type I Collagenase

in a 15mL tube

5.1 Add 8.0 mL of

1m

Hanks Balanced Salt Solution (HBSS) without Ca<sup>2+</sup> Mg<sup>2+</sup> **Thermo Fisher Scientific Catalog #88284**

5.2 Warm up in a 37 °C water bath before using

30m

## 6 Enzyme Cocktail 2 (EC2)–4 samples

2m

Dilute the Trypsin EDTA **Gibco - Thermo Fisher Scientific Catalog #25-051-Cl.** in

Hanks Balanced Salt Solution (HBSS) without Ca<sup>2+</sup> Mg<sup>2+</sup> **Thermo Fisher Scientific Catalog #88284**

to make a 0.05 Mass / % volume of trypsin/EDTA solution

6.1 Warm up in a 37 °C water bath before using

30m

## Tissue Dissociation and Cell Isolation Procedure

7m

7 Punch out the skins from mice using a sterile

5m

Standard Biopsy Punches, Disposable Standard biopsy punch; 6mm **Thermo Fisher Catalog #12460412**

and immediately float on WS1 on ice

8 Peel off any fat tissue beneath the skin using a blade scalpel blades










2m



- 9 Mince the fat tissue-cleaned skin samples into smaller pieces using sterile scissors while holding the tissue with sterile forceps on ice 2m
- 10 Incubate the minced tissues in 2 ml of EC1 in a 6-well plate at a 37 °C incubator with a 5% CO<sub>2</sub> supplementation for two hours 2h
- 02:00:00
- 10.1 Shake the plate orbitally in the incubator every 00:15:00 for better digestion 15m
- 11 Remove a piston from a sterile, 10 mL syringes **Becton Dickinson (BD) Catalog #BD 309695** to smash the partly-digested tissues in the plate 1m
- 12 Place a Falcon® 70 µm cell strainer **Corning Catalog #352350** into a 50mL tube and sift through the cell/tissue mixture 1m
- 13 Wash the EC1-digested tissue on the cell strainer with WS1 thrice 3m
- 14 Save the flow-through on ice
- 15 Invert the cell strainer with tissue remnants and place it on a sterile 6-well plate 1m
- 16 From the sieve-through side of the cell strainer, add 2mL of EC2 aiming at tissue remnants 1m
- 17 Incubate the undigested tissue at a 37 °C incubator with a 5% CO<sub>2</sub> supplementation for 00:15:00 15m
- 18 Remove a piston from a sterile, 10 mL syringes **Becton Dickinson (BD) Catalog #BD 309695** to smash the partly-digested tissues in the plate 1m





- 19 Place a sterile  Falcon® 70 µm cell strainer **Corning Catalog #352350** into the 50mL tube obtained at Step 15 to pool and sift through the cell/tissue mixture 1m 
- 20 Wash the trypsin-digested tissue on the cell strainer with WS1 thrice 1m 
- 21 Place a sterile  Falcon® 40 µm Cell Strainer **Corning Catalog #352340** into a clean 50mL tube and sift through the cell suspension 1m
- 22 Centrifuge the pooled cell suspensions at  300 x g, 20°C, 00:10:00 10m
- 23 Carefully aspirate the supernatant and resuspend the pellet in 1mL of WS2 2m 
- 24 Centrifuge the pooled cell suspensions again at  300 x g, 20°C, 00:05:00 5m
- 25 Carefully aspirate the supernatant and resuspend the pellet in 500µL of WS2 1m 
- 26 Determine the cell viability by mixing 10µl of the cell suspensions with 10µl of trypan blue 2m 

### Equipment

#### Countess 3 FL Automated Cell Counter

Automated Cell Counter

ThermoFisher scientific

AMQAF2000

<https://www.thermoFisher.com/th/en/home/life-science/cell-analysis/cell-analysis-instruments/automated-cell-counters/models/countess-3-fl.html>

NAME

TYPE

BRAND

SKU

LINK



26.1 Take at least four counts per sample

8m



26.2 Remember to tick "Trypan Blue correction" on the Countess cell counter

1m



27 Average cell numbers for each sample and dilute each cell suspension with WS2 to adjust the cell number to 700-1200 cells/ $\mu$ L as indicated in **10X Genomics Chromium Next GEM Single Cell 3' Reagent Kits v3.1 Protocol (CG000315)**

5m

27.1 Aim for cell viability above 70% for each sample



Expected result

Cell viability: 70% per count

27.2 If the cell viability is below 70%, you may choose to remove dead cells using a dead cell removal kit



28 Immediately proceed with **Chromium Next GEM Single Cell 3' Reagent Kits v3.1 Protocol (CG000315)**