

May 07, 2025

🌐 Clarification of wastewater samples

DOI

<https://dx.doi.org/10.17504/protocols.io.5jyl8qo99l2w/v1>

Victor Mabasa¹, Natasha Singh¹, Emmanuel Phalane¹, Mokgaetji Macheke¹, Sipho Gwala¹, Thabo Mangena¹, Lethabo Monamets¹, Lebohang Rabotapi¹, Nkosenhle Ndlovu¹, Fiona Els^{1,2}, Sibonginkosi Maposa¹, Said Rachida¹, Kerrigan McCarthy^{1,3}, Mukhlid Yousif^{1,3}

¹National Institute for Communicable Diseases; ²Gauteng City-Region Observatory; ³University of Witwatersrand

Wastewater Genomics S...



Victor Mabasa

National Institute for Communicable Diseases

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5jyl8qo99l2w/v1>

Protocol Citation: Victor Mabasa, Natasha Singh, Emmanuel Phalane, Mokgaetji Macheke, Sipho Gwala, Thabo Mangena, Lethabo Monamets, Lebohang Rabotapi, Nkosenhle Ndlovu, Fiona Els, Sibonginkosi Maposa, Said Rachida, Kerrigan McCarthy, Mukhlid Yousif 2025. Clarification of wastewater samples. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.5jyl8qo99l2w/v1>



License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 07, 2025

Last Modified: May 07, 2025

Protocol Integer ID: 217828

Keywords: Wastewater solids, wastewater, NICD, South Africa, WBE, WBS, wastewater genomics syndicate, wastewater sample processing protocol, wastewater sample, clarification of wastewater sample, supernatant components of wastewater, wastewater supernatant, nucleic acid extraction, different pathogen, such as ultrafiltration, ultrafiltration, bacteria, extraction, virus, national institute of communicable disease, sample

Abstract

This wastewater sample processing protocol by the Wastewater Genomics Syndicate at the National Institute of Communicable Diseases (NICD) is designed to enrich bacteria and viruses from the settled solids and supernatant components of wastewater, as different pathogens partition differently within this matrix. The end products are waterless solids and clarified wastewater supernatants, ready for downstream applications such as ultrafiltration, bead enrichment and nucleic acid extraction.

Materials

Equipment

–80°C freezer to store settled solids samples

4°C refrigerator storage of samples

Biological safety cabinet (Class II)

Refrigerated centrifuge (4°C) with adapters to accommodate 500 mL centrifuge bottles

Consumables

2 mL cryo-tubes

250 mL or 500 mL Conical centrifuge bottles

Waterproof Labels

Waste container

Plain wooden sticks

Solutions/Reagents

70% ethanol

1% Virkon

Sample Collection & Handling

- 1 Wearing appropriate personal protective equipment (PPE), collect a 1 000 mL grab wastewater sample from each location using a sterile bottle. Transport the samples to the laboratory under cold chain conditions and store them at 4 °C. To minimise microbial degradation, proceed to the next step within 23 hours of sample collection.

Preparing for centrifugation

- 2 Conduct all work in a Class II biological safety cabinet while wearing appropriate PPE. Clean the workspace using 1% Virkon followed by 70% ethanol. Double-line the waste container with biohazard waste bags, ensure all filter tip racks are refilled, and label a 250 mL conical centrifuge bottle for each wastewater sample.

Centrifugation

- 3 Transfer each wastewater sample to its corresponding conical centrifuge bottle, taking care not to overfill to prevent spillage during centrifugation.
- 4 Centrifuge the samples at $4\,700 \times g$ for 10 minutes at 4 °C.
- 5 Carefully transfer the supernatant to a clean, marked 1 000 mL bottle without disturbing the settled solids.
- 6 Repeat steps 3 to 5 until the original sample is finished.

Storage

- 7 Aliquot 6×5 mL cryotubes of the clarified supernatant and store them at -80 °C for long-term preservation. Keep the remaining clarified wastewater at 4 °C for other downstream applications.
- 8 Using clean, plain wooden sticks, collect settled solids and transfer them into labelled 2 mL cryotubes. Briefly centrifuge to remove excess liquid, then store the tubes at -80 °C until needed for downstream applications.