



May 16, 2023

Version 2

Complete Medium or Complete Medium Xylose (from Leach, Lang and Yoder 1982) V.2

DOI

dx.doi.org/10.17504/protocols.io.kqdg39mxqg25/v2

Megan Mcdonald¹

¹University of Birmingham, School of Biosciences



Megan Mcdonald

University of Birmingham

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kqdg39mxqg25/v2>

Protocol Citation: Megan Mcdonald 2023. Complete Medium or Complete Medium Xylose (from Leach, Lang and Yoder 1982). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.kqdg39mxqg25/v2> Version created by **Megan Mcdonald**

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: May 16, 2023

Last Modified: May 16, 2023

Protocol Integer ID: 81948

Keywords: maintenance of cochliobolus carbonum, cochliobolus carbonum, complete medium xylose

Abstract

For the growth and maintenance of Cochliobolus carbonum and Cochliobolus victoriae

Materials

Trace Minerals

Boric Acid (H_3BO_3) CAS: 10043-35-3

Cupric Sulfate (CuSO_4) CAS: 7758-98-7

Potassium Iodide (KI) CAS: 7681-11-0

Manganese (II) sulfate monohydrate $\text{MnSO}_4 \cdot \text{H}_2\text{O}$ CAS: 10034-96-5

Sodium permanganate monohydrate $\text{NaMoO}_4 \cdot \text{H}_2\text{O}$ CAS: 79048-36-5

Zinc sulfate heptahydrate $\text{ZnSO}_4 \cdot 7 \text{H}_2\text{O}$ CAS: 7446-20-0

Iron (III) chloride hexahydrate $\text{FeCl}_3 \cdot 6\text{H}_2\text{O}$ CAS: 10025-77-1

Salts

Calcium nitrate tetrahydrate $\text{Ca}(\text{NO}_3)_2 \cdot 4 \text{H}_2\text{O}$ CAS: 13477-34-4

Potassium Phosphate monobasic KH_2PO_4 CAS: 7778-77-0

Magnesium Sulfate heptahydrate $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$ CAS: 10034-99-8

Sodium Chloride NaCl CAS: 7647-14-5

Media

Yeast Extract CAS: 8013-01-2 (some people report differences between difference sources)

Casein Digests: Still a bit of work to do here to determine, which of these is best (including new plant-based sources). Big differences in prices between different digest types.

Acid hydrolysed Casein CAS: 65072-00-6

Peptone from casein, tryptic digest OR pancreatic digest (tryptone) CAS: 91079-40-2

Troubleshooting



Make Micronutrients Solution

- 1 9 mg H_3BO_3
58.5 mg $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$
1.95 mg KI (Potassium Iodine)
9 mg MnSO_4
7.6 mg NaMoO_4
822 mg $\text{ZnSO}_4 \cdot 6\text{H}_2\text{O}$
139.8 mg $\text{FeCl}_3 \cdot 6\text{H}_2\text{O}$

in 300 mL ddH₂O and filter sterilise

Citation:

Heterokaryosis and Parasexuality in the Fungus Ascochyta Imperfecta Author(s): K. E. Sanderson and A. M. Srb Source: American Journal of Botany , Jan., 1965, Vol. 52, No. 1 (Jan., 1965), pp. 72-81 Published by: Wiley Stable URL: <https://www.jstor.org/stable/2439977>

Make 100x Salt Solutions A and B

2 100X Salt Solution A

10g $\text{Ca}(\text{NaO}_3)_2 \cdot 4\text{H}_2\text{O}$
100 mL ddH₂O

Autoclave

Citation:

Leach, J., Lang, B. R. & Yoder, O. C. *Microbiology* **128**, 1719-1729, [doi:https://doi.org/10.1099/00221287-128-8-1719](https://doi.org/10.1099/00221287-128-8-1719) (1982).

3 100X Salt Solution B

2 g $\text{KH}_2\text{PO}_4 \cdot 7\text{H}_2\text{O}$
1.5g NaCl

100 mL H₂O
pH 5.3



Autoclave

Citation:

Leach, J., Lang, B. R. & Yoder, O. C. *Microbiology***128**, 1719-1729,
[doi:https://doi.org/10.1099/00221287-128-8-1719](https://doi.org/10.1099/00221287-128-8-1719) (1982).

Make Complete Medium (CM) or Complete Medium Xylose (CMX)

4 Complete Medium Base

10 g glucose OR xylose (substitute glucose for xylose for CMX medium)
1 g Yeast Extract
0.5 g Acid-hydrolysed Casein
0.5 g Enzyme-hydrolysed casein
20g Agar
10 mL Salt A
10 mL Salt B

Make up to 1000 mL with ddH₂O and Autoclave

After autoclaving add:

1 mL sterilised micronutrient solution

Citation:

Leach, J., Lang, B. R. & Yoder, O. C. *Microbiology***128**, 1719-1729,
[doi:https://doi.org/10.1099/00221287-128-8-1719](https://doi.org/10.1099/00221287-128-8-1719) (1982).

Complete Medium for Sporulation

5 Complete Medium Base

0.5 g glucose OR xylose (substitute glucose for xylose for CMX medium)
20 g Sorbose
1 g Yeast Extract
0.5 g acid hydrolysed casein
0.5 g enzyme hydrolysed casein
20g Agar

10 mL Salt A
10 mL Salt B

Make up to 1000 mL with ddH₂O and Autoclave



After autoclaving add:
1 mL sterilised micronutrient solution

Protocol references

Citation:

Heterokaryosis and Parasexuality in the Fungus Ascochyta Imperfecta Author(s): K. E. Sanderson and A. M. Srb
Source: American Journal of Botany , Jan., 1965, Vol. 52, No. 1 (Jan., 1965), pp. 72-81 Published by: Wiley Stable
URL: <https://www.jstor.org/stable/2439977>

Leach, J., Lang, B. R. & Yoder, O. C. *Microbiology***128**, 1719-1729, [doi:https://doi.org/10.1099/00221287-128-8-1719](https://doi.org/10.1099/00221287-128-8-1719) (1982).