

Apr 16, 2025

Version 1

Concentration and total RNA extraction of viruses from nasal swabs V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.261geepkdg47/v1>

Vanessa Smilansky¹, Simon Grimm¹, Olivia S Hershey¹

¹SecureBio

SecureBio Detection



Vanessa Smilansky

SecureBio

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261geepkdg47/v1>

External link: <https://naobservatory.org>

Protocol Citation: Vanessa Smilansky, Simon Grimm, Olivia S Hershey 2025. Concentration and total RNA extraction of viruses from nasal swabs. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.261geepkdg47/v1>

Manuscript citation:

Vanessa Smilansky, Simon Grimm, Olivia S Hershey 2025. Concentration and nucleic acid extraction of viruses from nasal swabs. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.261geepkdg47/v1>

Version created by [Vanessa Smilansky](#)

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 21, 2025

Last Modified: April 16, 2025

Protocol Integer ID: 124806

Keywords: nasal swab sampling, respiratory sampling, viral concentration, viromics, metagenomics, total rna extraction of virus, total rna per pooled swab, viruses present in nasal secretion, nucleic acid extraction of virus, total rna extraction, qiaamp viral rna mini kit, nasal swab sample, extracting total rna, rna extraction, nasal swabs this protocol, nasal swab, total rna, nucleic acid extraction, viral particle, rna, untargeted rna, nasal secretion, pooled swab, nucleic acid, virus, swab, nao protocol concentration, extraction, concentrated product

Funders Acknowledgements:

Open Philanthropy

Grant ID: NA

Abstract

This protocol describes the workflow for concentrating and extracting total RNA from a pool of nasal swab samples, aimed at untargeted RNA sequencing to identify viruses present in nasal secretions. It is based on the NAO protocol **Concentration and nucleic acid extraction of viruses from wastewater influent V.2**. To summarize, nasal secretions are dissociated from swabs ($n \geq 25$) into solution, bacteria and human cells are removed, and viral particles are concentrated to a final volume of $\sim 600 \mu\text{L}$. Total RNA is extracted from the concentrated product using the QIAamp Viral RNA Mini Kit according to the manufacturer's instructions. We estimate that this protocol yields an average of 84 ng (SE = 15 ng) of total RNA per pooled swab.

Attachments



[Estimate of total RN...](#)

47KB

Guidelines

RNA processing and handling: Please review **Protocol Note: Working with RNA Samples** before handling RNA samples.

Materials

Materials:

(for one pool of nasal swabs and one negative control)

- Nasal swabs in **Falcon 50 mL Conical Centrifuge Tubes** (VWR No. 21008-178): ≤ 20 swabs per tube
- Clean centrifuge tubes equal to number of tubes containing nasal swabs + one negative control
- 2 x **InnovaPrep Ultra CPT (Unirradiated)** tips for the InnovaPrep Concentrating Pipette Select
- 2 x **0.45 μ m PES 75 mm vacuum filtration tops** (VWR No. 10040-470)
- 2 x 100 mL pyrex bottles
- 4 x 5 mL centrifuge tubes
- 1 x 50 mL serological pipette (VWR No. 170358N)
- 1.5 mL microcentrifuge tubes (VWR No. 1420-2600)
- PREempt wipes (VWR No. 10822-456)
- Parafilm (Millipore Sigma No. HS234526C)
- RNaseZap
- Kimwipes (VWR No. 34120)
- Paper towels
- Filtered micropipette tips

Reagents:

- Ice
- **HBSS (Hank's Balanced Salt Solution)**
- **CP Select Elution buffer (Tris)**
- **CP Select Storage Fluid**
- All buffers in **QIAamp Viral RNA Mini Kit (50)**
- 100% ETOH
- 80% ETOH
- Nuclease-free water (NFW)

Equipment:

- **InnovaPrep Concentrating Pipette Select**
- **Attachment for vortexer** (Cole-Parmer No. UX-04724-89)
- Scientific Industries Vortex Genie 2
- Vacuum Line
- Micropipettes (1000 μ L, 200 μ L, 10 μ L, 2 μ L, 1 μ L) and holder
- Timer
- Floor Centrifuge (Ex: Beckman Coulter Avanti J series)
- Rotor compatible with 50 mL tubes (Ex: Beckman Coulter JA-14.50)

Safety warnings

! Biosafety precautions: Nasal swab samples will be received and stored using both primary and secondary containment. The primary container should remain inside the secondary container until processing begins. Samples must only be handled within a designated biosafety cabinet or fume hood. All laboratory personnel handling these samples are required to wear safety glasses, gloves, lab coats, and N95 face masks to minimize biological aerosol risk. When handling samples in the During transportation between processing stations, samples will be placed in a secondary container that has been disinfected with PreEmpt. In addition, all surfaces—including the outside of the fume hood/sash, centrifuge lid, and rotor—will be wiped down with PreEmpt. Any containers that come into contact with nasal swab samples will be disinfected by soaking with a 10% bleach solution for at least 20 minutes.

Before start

Read the 'Warnings' section for biosafety precautions necessary for handling nasal swab samples. Prepare the biosafety cabinet and fume hood for nasal swab handling. Gather materials and reagents. Ensure proper PPE.



Stage 1: Sample preparation


7m

- 1 In a biosafety cabinet, remove the nasal swabs ($n \geq 25$) from their secondary containers. Ensure that the nasal swabs are evenly distributed among Falcon 50 mL Conical Centrifuge Tubes, with no tube containing more than 20 swabs.

5m

Note



Falcon 50 mL Conical Centrifuge Tubes should be used as primary containers for collecting nasal swabs to streamline sample processing.

- 2 Using a 50 mL serological pipette, add the equivalent of  1 mL of HBSS per swab to each tube, taking care not to come into direct contact with the tube. For a negative control, add an equal volume of HBSS to a clean centrifuge tube. Cap and parafilm each tube.

2m

Stage 2: Dissociation of viral particles and removal of large solids

25m

- 3 Place all tubes on a vortexer using a 50 mL tube adapter. Shake for  1000 rpm ,  00:00:30 .




1m

- 4 Excluding the negative control, transfer the liquid contents into clean centrifuge tubes. Leave the swabs in their original tube, cap, and discard. Cap and parafilm each new tube and place in a secondary container.

2m

Note

For efficient centrifugation and pelleting in 50 mL centrifuge tubes, aim for a volume between 10-20 mL.

- 5 Transfer all tubes, including the negative control, from the biosafety cabinet to a centrifuge equipped with the appropriate rotor and/or adapters (e.g., Beckman Coulter Avanti J series with JA-14.50 rotor). Centrifuge the samples at  1200 x g, 4°C,  00:10:00 . After centrifugation is complete, wait  00:10:00 (as recommended by EHS) to allow aerosols to settle before opening the centrifuge. Remove tubes from the centrifuge, place them in a secondary container, and transfer them to the fume hood.

22m



Stage 3. Removal of suspended solids

22m

- 6 In the fume hood, prepare two separate 0.45 μm vacuum filtration apparatuses by attaching each filtration unit to a clean 100 mL pyrex bottle. 2m
- 7 Decant the supernatant from all tubes containing nasal swab samples directly into a 0.45 μm vacuum filtration top, taking care not to dislodge the pellets. 2m
- 8 Begin vacuum filtration by capping the vacuum filtration top and opening the vacuum line. When complete, cap the pyrex bottle and set aside. 15m

Note

The duration of this step will vary depending on the consistency of the supernatant.

- 9 For the negative control, decant directly into a 0.45 μm vacuum filtration top. There should not be a pellet. Perform vacuum filtration as was done for the nasal swab samples, cap the pyrex bottle, and set aside. 3m

Stage 4: Concentration of viral particles

25m

- 10 Perform the "Start-up" protocol for the InnovaPrep Concentrating Pipette Select.
 - Turn on the Concentrating Pipette, and navigate to "Maintenance" and then "Start-up". Follow the prompts.
 - Check that the maintenance tip is in place.
 - Place the waste line in the proper position.
 - Remove the storage fluid line and insert the foam elution canister.
 - Ensure that the screen reads "ULTRA".
- 11 Run the concentration protocol for the filtered nasal swab sample.
 - Remove the maintenance tip and place a fresh Ultra CPT into the tip port. Lower the tip into the sample. 10m

Note

Ensure that the CPT is as close to the bottom of the sample bottle as possible. If necessary, the bottle can be balanced on its edge while a weighted object holds down the top of the Concentrating Pipette. A bottle of HBSS can be used as the weighted object.



- Press "Start Run", allow the Concentrating Pipette to aspirate the waste, and wait until it beeps to signal the end of the run.
- While holding a 5 mL Eppendorf tube under the Tip, press "Elute" and catch the foam that is dispensed. The eluate should be ~ $600\ \mu\text{L}$; this can vary slightly depending on the volume of unconcentrated sample.
- About $1\ \text{mL}$ of unconcentrated filtrate should remain in the pyrex bottle. Pipette $600\ \mu\text{L}$ and transfer to a 5 mL Eppendorf tube.

12 Run the concentration protocol for the filtered negative control. Perform in the same manner as the nasal swab sample, including removal of a $600\ \mu\text{L}$ unconcentrated aliquot. Make sure to use a fresh tip.

5m

13 Perform the "Shut Down" protocol for the InnovaPrep Concentrating Pipette Select.

- Navigate to "Maintenance" and then "Shut Down".
- Place the maintenance tip into the tip port.
- Remove the elution canister.
- Check to ensure that there is adequate storage fluid and insert the storage fluid line.
- Turn off the device and remove the waste line.

5m

Stage 5: Total RNA extraction - QIAamp® Viral RNA Mini Kit

56m

14 Gather the materials and reagents for the QIAamp® Viral RNA Mini Kit in a biosafety cabinet. Equilibrate samples to room temperature.

- Add $25\ \text{mL}$ of 100% ETOH to Buffer AW1.
- Add $30\ \text{mL}$ of 100% ETOH to Buffer AW2.
- Check Buffer AVL for precipitate, and if necessary incubate at $80\ ^\circ\text{C}$ until the precipitate is dissolved. **It is not necessary to add carrier RNA.**




5m

Note

Refer to the "Guidelines" section for instructions on processing and handling RNA samples.



15 Add $2.4\ \text{mL}$ (4x volume) of Buffer AVL to each of the four 5 mL tubes containing $600\ \mu\text{L}$ of concentrated sample and control, and unconcentrated sample and

13m

control. Mix by pulse-vortexing for  00:00:15 . Incubate at  Room temperature for  00:10:00 ; spin down.

Note



The volume of Buffer AVL is scaled according to the manufacturer's recommendation.

16 Add  2.4 mL 100% ETOH to each tube. Mix by pulse-vortexing for  00:00:15 ; spin down.



3m

Note




The volume of 100% ETOH should be equal to the volume of Buffer AVL from step 15.

17 For each sample and control, transfer  700 μ L to a designated spin-column tube. Centrifuge at  6000 x g, 00:02:00 . Discard flow-through. Repeat 8x, or until entire volume has passed through spin-column. Take care to transfer to the correct tube each time. Once finished, transfer the spin-columns to clean collection tubes.






20m

18 Add  500 μ L Buffer AW1 to each tube. Centrifuge at  6000 x g, 00:01:00 . Discard flow-through.

3m

19 Add  500 μ L Buffer AW2. Centrifuge at  6000 x g, 00:03:00 . Transfer spin-columns to clean collection tubes. Centrifuge at  max rpm, 00:01:00 to remove residual buffer. Transfer spin-columns to designated 1.5 mL Eppendorf tube.

7m

20 Add  50 μ L NFW to spin-columns and incubate at  Room temperature for  00:02:00 . Centrifuge at  6000 x g, 00:01:00 to collect eluate. Discard spin-columns. Store eluates at  -80 $^{\circ}$ C .

5m

Note

Although the QIAamp® Viral RNA Mini Kit is designed for RNA isolation, some DNA co-precipitation is unavoidable. To eliminate residual DNA, DNase treatment of the eluate is necessary.

Expected result

We estimate that this protocol yields an average of 84 ng (SE = 15 ng) of total RNA per swab from the concentrated products (see Attachment for details). Total RNA per swab from the unconcentrated aliquots is not estimated because these yields are typically below the limit of detection for the Qubit HS RNA assay. Yields are expected to vary depending on the amount of the raw material present on each swab in the pool.

Protocol references

[Concentration and nucleic acid extraction of viruses from wastewater influent V.2](#)

[Innovaprep FluidPrep Concentrating Pipette Select](#)

[QIAamp® Viral RNA Mini Kit](#)