

May 09, 2025

Confocal microscopy to assess TelC and iGluSnFR expression in fixed mouse brain slices

 [Nature Communications](#)

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.eq2ly6q4wgx9/v1>

Safa Bouabid^{1,2}, Mark Howe^{1,2}

¹Department of Psychological & Brain Sciences, Boston University, Boston, MA, USA;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815, USA



Cláudia C. Mendes

University of Oxford



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.eq2ly6q4wgx9/v1>

Protocol Citation: Safa Bouabid, Mark Howe 2025. Confocal microscopy to assess TelC and iGluSnFR expression in fixed mouse brain slices. [protocols.io https://dx.doi.org/10.17504/protocols.io.eq2ly6q4wgx9/v1](https://dx.doi.org/10.17504/protocols.io.eq2ly6q4wgx9/v1)

**Manuscript citation:**

Bouabid *et al.*, (2025) Distinct spatially organized striatum-wide acetylcholine dynamics for the learning and extinction of Pavlovian associations.

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 05, 2025

Last Modified: May 09, 2025

Protocol Integer ID: 119623

Keywords: mouse, brain, confocal microscopy, glutamate recording, fixed mouse brain slice, mouse brain slices this protocol, encoded glutamate sensor, glutamate sensor, medial striatum, following tetanus toxin light chain, tetanus toxin light chain, viral injection in the dorsal, confocal microscopy, telc, dorsal, striatum, brain

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP)

Grant ID: ASAP-020370

Abstract

This protocol describes the tissue processing steps following (a) glutamate recordings using the genetically encoded glutamate sensor (iGluSnFR) and (b) behaviour experiments following tetanus toxin light chain (TelC) viral injection in the dorsal medial striatum.

Before start

This protocol was performed in **ChAT-IRES-Cre** mice (B6;129S6-Chattm2(cre)Lowl/J, JAX stock number 006410).

To block acetylcholine release from cholinergic interneurons, we injected **Cre-dependent Tetanus toxin light chain (TelC)** viral vector (RRID:Addgene_135391) in the dorsal medial striatum (aDMS) of ChAT-IRES-Cre mice.

To measure rapid changes in glutamate release, we expressed the **genetically encoded glutamate sensor iGluSnFR** selectively in anterior dorsal striatum (aDS) cholinergic interneurons of ChAT-IRES-Cre mice.



PFA Fixation by Intracardiac Perfusion & Brain Dissection

- 1 Perfuse intracardially with phosphate-buffered saline (PBS 1%, Fisher Scientific), followed by paraformaldehyde (PFA 4% in 1% PBS, Fisher Scientific).
- 2 Carefully extract brains and fix them overnight in 4% PFA dissolved in PBS.
- 3 Transfer brains to solution of 40% sucrose in 1% PBS (until the brains sank).

Cryosectioning & Mounting

- 4 Slice brains (50µm thickness) with a cryostat (Leica CM3050 S).
- 5 Mount coronal sections on super frost slides and cover slipped with Vectashield antifade mounting medium (Vector Laboratories, H-1900).

Image Acquisition using Confocal Microscopy

- 6 Acquire images on a Zeiss LSM 800 laser scanning confocal microscope.

Note

TelC and iGluSnFr fluorescence were not immuno-enhanced.