

Sep 03, 2020

Version 1

🌐 Coronavirus Lateral Flow Assay (LFA) sample preparation protocol V.1

📁 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.bksfkwbm>

peijun he¹

¹Highfield Diagnostics, UK


Coronavirus Method De...

XPRIZE Rapid Covid Tes...

1 more workspace



peijun he



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bksfkwbm>

Protocol Citation: peijun he 2020. Coronavirus Lateral Flow Assay (LFA) sample preparation protocol. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bksfkwbm>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol in our laboratory and it is working

Created: September 03, 2020

Last Modified: September 03, 2020

Protocol Integer ID: 41511

Keywords: sample preparation protocol coronavirus lateral flow assay, coronavirus lateral flow assay, lfa, assay

Abstract

Coronavirus Lateral Flow Assay (LFA) sample preparation protocol

Materials

MATERIALS

 Triton X-100 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #93426**

 DPBS 10X with Calcium and Magnesium **Gibco - Thermo Fisher Scientific Catalog #14-080-055**

 EDTA **Merck MilliporeSigma (Sigma-Aldrich) Catalog #ED2SS**

 Tributyl phosphate **Catalog #240494**



- 1 For each clinical sample to be processed, label 1 screwcap tube with the patient's ID.
Make  400 µL
aliquots of Lysis buffer in each tube.

Lysis buffer:

[M] 1 % volume Triton X - 100 [M] 0.3 % volume Tributyl phosphate 1x complete
EDTA-free [M] 5 millimolar (mM) EDTA in 1x PBS

- 2 Take out the swabs out of the  -80 °C freezer immediately before the analysis.
Leave to thaw at RT for  00:05:00 and insert into the respective screw-cap tube.
Note: handling times within a CL3 lab may be longer
than expected, so do not take the swabs out of the freezer until all tubes are properly
labelled and
containing  400 µL Lysis buffer in them.

- 3 'Plunge' each swab up and down in Lysis Buffer for  00:15:00 , constantly 'twisting'
the swab so its
sides are brushed against the side of the tube. Note: Ideally, do not handle more than 3-
4 swabs at a time
so as to spend enough time twisting each of them to 'encourage' maximum dissociation
of Nucleocapsid
protein from the swab and lysis of virions.

- 4 After  00:15:00 , discard the swabs appropriately (e.g. in a second container that will
then need to be
autoclaved before disposal), making sure to recover as much volume out of the swab as
possible ( 100 µL
will be lost as absorbed by the swab foam) and re-cap the tube.



- 5 Samples can be tested immediately by lateral flow device or stored at  -80 °C for further analysis.
- 6 Lateral flow assay could be performed directly inside a CL3 lab or within a Class I cabinet after decontamination of the outside of the tubes before taking them out of the CL3 lab. For the lateral flow assay, follow the lateral flow assay operation protocol.