

Oct 28, 2024

CRISPR/Cas9 genome editing in Candida albicans

DOI

dx.doi.org/10.17504/protocols.io.3byl4we1rvo5/v1

Chien-Der Lee¹, Alexander johnson¹

¹University of California San Francisco

protocols.io Ambassadors



Chien-Der Lee

UTSW

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl4we1rvo5/v1>

Protocol Citation: Chien-Der Lee, Alexander johnson 2024. CRISPR/Cas9 genome editing in Candida albicans. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.3byl4we1rvo5/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 26, 2024

Last Modified: October 28, 2024

Protocol Integer ID: 110965

Keywords: cas9 genome editing in candida, cas9 genome editing, homozygous knockout strain in diploid candida albican, generated cas9, digested cas9, diploid candida albican, transient cas9 expression method, crispr, candida, homozygous knockout strain, genome, plasmid, modification, knockout strain

Abstract

This protocol aims to create a homozygous knockout strain in diploid *Candida albicans*, based primarily on the transient Cas9 expression method by Min et al. in *mSphere* (2016), with some modifications. Notably, it uses plasmid-digested Cas9 (pp2280) instead of PCR-generated Cas9, optimizing both time and cost. The knockout strain should be successfully obtained within one week of the transformation reaction.

Guidelines

This protocol combines the classic SAT1-Flipper system with Cas9-induced double-strand breaks to streamline the creation of homozygous knockouts in diploid *Candida albicans*.


Materials

 pADH34_protocol_io.gbk 18KB

 pp2280_protocol_io.gbk 72KB

Troubleshooting

Safety warnings

 Due to the nature of the SAT1-Flipper system, the knockout strain will retain an FLP site scar. However, this scar can serve as a CRISPR/Cas9 target site for further modifications if needed.



Prepare for DNA required for transformation

1 There are only three DNA components required for transformation:

1. P_{ENO1}-Cas9-t_{CYC1} construct (Stul digested pp2280)
2. P_{SNR52}-gRNA-scaffold-terminator
3. SAT1-Flipper donor DNA with homologous end

2 **Component 1.** P_{ENO1}-Cas9-t_{CYC1} construct (Stul digested pp2280)

This restriction enzyme digested fragment provides the transient Cas9 expression in *Candida albicans* without integrating into the genome, as it does not have homology sequence to the genome.

I modified the plasmid **pADH140** (Aaron Hernday) to pp2280 (see materials) so that it can be digested by Stul for transient Cas9 expression.

Gel extraction after Stul digestion is optional.

3 **Component 2.** P_{SNR52}-gRNA-scaffold-terminator

Two independent PCRs to amplify P_{SNR52}-gRNA, and gRNA-scaffold-terminator from **pADH143** (Aaron Hernday). Then followed by fusion PCR (alternatively, Gibson assembly or isothermal assembly also work) to fuse these two fragment together.

4 **Component 3.** SAT1-Flipper donor DNA with homologous end

PCR the classic SAT1-Flipper cassette with 80 bp flanking region to the targeted gene.

Note: Generally speaking, the longer the flanking region, the higher the knockout efficiency. 60~100 bp flanking region should all work.

Transformation

1w 2d

5 Streak out strain from -80°C to YPD plate, 30°C sit for 2 days

2d

6 Pick single colony and inoculate to 3 mL YPD, roller drum 30°C O/N

1d

7 In the morning, back dilute the O/N culture to 10 mL YPD per transformation reaction (~0.2 OD/mL to start). 30°C 200rpm shaking in flask



8 About 3~4 hours, the cell density reaches to ~0.8 OD/mL

4h

9 Spin down the cell by 3000xg 2min, wash with 10 mL water, and then spin down again

10m

10 Resuspend all cell with: (all need to be sterilized)

	A	B
	50 % PEG 3350	800 µL
	1 M Lithium Acetate	100 µL
	10 mg/mL heat denatured salmon sperm DNA	50 µL
	10X TE pH 8.0	100 µL
	Component 1	50 µL
	Component 2	50 µL
	Component 3	50 µL

11 44°C heat shock for 15min (we use water bath here)

15m

12 Spin down cell at 3000xg 1min, wash twice with 1 mL YPD

10m

13 Resuspend all cell in 2 mL YPD, 30°C roller drum for at least 4 hours (O/N is fine).

4h

14 Spin down all cell at 3000xg 1min, wash twice with water

10m

15 Resuspend all cell in 1 mL water and take 100 µL to YPD+NAT plate. 30°C sit for >2days

2d

16 Pick several colonies to new YPD+NAT plate (Due to the high efficiency, 4 colonies should be enough)



- 17 Inoculate cells to 2 mL YP+2% maltose 30°C O/N to Flip out SAT1 cassette 1d
- 18 Streak small amount of O/N culture to YPD. 30°C sit for >2days 2d
- 19 Patch single colonies to YPD. Then test if it lack the NAT-resistance. (indicating the SAT1-Flipper has been removed) 1d
- 20 Store strains lacking NAT-resistance.

Appendix 1. PCR primer design information: Ade2 Knockout as an example

21 **Component 2.** P_{SNR52}-gRNA-scaffold-terminator (pADH143 as template)

Component 2 PCR-A (P_{SNR52})

	A	B
	templat e	pADH143
	F1	catctaactcaactcccagat
	R1	caaattaaaaatagtttacgcaag

Component 2 PCR-B (gRNA-scaffold-terminator)

	A	B
	templat e	pADH143
	F2	cttgcgtaaactatttttaatttgCAACAATCATACGACCTAATgttttagagctagaaatagc a
	R2	taaaaaaaCTCGAGAAAAAAGCAC

F1, R1 and R2 are universal for different KO. Only F2 (contain gRNA) is gene specific

Component 2 PCR-C (fusion PCR of A and B)

This PCR product will be used for transformation)

A	B
template	PCR-A and PCR-B product, equal molar ratio (50ng each is good for 50 µL PCR reaction)
F1	catctaataactcccagat
R2	taaaaaaaCTCGAGAAAAAAGCAC

The PCR-A and PCR-B product need to be DpnI digested and PCR clean up to remove both the template plasmid and primers.

Make sure to use the correct primer pair here: F1 and R2

Component 3. Primer for SAT1-Flipper donor DNA with homologous end (pADH34 as template)

A	B
template	pADH34 (or any plasmid with SAT1-Flipper should work)
F	AAAAACAACCAACCAACCCTTAACCCATTAACGAATTAAGATTTGTT CTATTTGACTACCAAGAATATAACCCATATTactagtgaattcgctcgag
R	TTCTTTTATATATTAATTATAATTTTATTGCACCATAACGTTTACTTGT TTAATATGCTATTGATATCTATATTTTTTcgctctagaactagtggatc

Appendix 2. gRNA design

- 22 We design the gRNA using Benchling tools following the protocol in Nguyen *et. al.*, *mSphere* 2017
https://journals.asm.org/doi/suppl/10.1128/mspheredirect.00149-17/suppl_file/sph002172275s1.pdf
- 23 Alternatively, *Candida albicans* gRNA sequence for each gene can be found in the supplementary information in Vyas *et. al.*, *Sci Adv.* 2015



Protocol references

Vyas VK, Barrasa MI, Fink GR. A *Candida albicans* CRISPR system permits genetic engineering of essential genes and gene families. *Sci Adv*. 2015;1(3):e1500248. doi: 10.1126/sciadv.1500248. PMID: 25977940; PMCID: PMC4428347.

Min K, Ichikawa Y, Woolford CA, Mitchell AP. *Candida albicans* Gene Deletion with a Transient CRISPR-Cas9 System. *mSphere*. 2016 Jun 15;1(3):e00130-16. doi: 10.1128/mSphere.00130-16. PMID: 27340698; PMCID: PMC4911798.

Nguyen N, Quail MMF, Hernday AD. An Efficient, Rapid, and Recyclable System for CRISPR-Mediated Genome Editing in *Candida albicans*. *mSphere*. 2017 Apr 26;2(2):e00149-17. doi: 10.1128/mSphereDirect.00149-17. PMID: 28497115; PMCID: PMC5422035.

Acknowledgements

We thank Aaron Hernday for providing the plasmids, Manning Huang and Carrie Graham for helpful discussion.