



Mar 15, 2024

## CRISPR inhibition (CRISPRi) of LINC01876 in hiPSCs and fbNPCs

DOI

[dx.doi.org/10.17504/protocols.io.j8nlkwnq6l5r/v1](https://dx.doi.org/10.17504/protocols.io.j8nlkwnq6l5r/v1)

anita.adami<sup>1</sup>

<sup>1</sup>Laboratory of Molecular Neurogenetics, Department of Experimental Medical Science, Wallenberg Neuroscience Center and Lund Stem Cell Center, BMC A11, Lund University, 221 84 Lund, Sweden.



**Raquel Garza**

Lund University

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlkwnq6l5r/v1>

**Protocol Citation:** anita.adami 2024. CRISPR inhibition (CRISPRi) of LINC01876 in hiPSCs and fbNPCs. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.j8nlkwnq6l5r/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** February 07, 2023

**Last Modified:** May 31, 2024

**Protocol Integer ID:** 76525

**Keywords:** ASAPCRN, crispr inhibition, crispr, crispri in hipsc, linc01876 in hipsc, crispri, hipsc, linc01876

**Funders Acknowledgements:**

**Aligning Science Across Parkinson's through the Michael J. Fox Foundation for Parkinson's Research**

Grant ID: ASAP-000520

**Swedish Research Council**

Grant ID: 2018-02694

**Swedish Brain Foundation**

Grant ID: FO2019-0098

**Cancerfonden**

Grant ID: 190326

**Barncancerfonden**

Grant ID: PR2017-0053

**NIHR Cambridge Biomedical Research Centre**

Grant ID: NIHR203312

**Swedish Society for Medical Research**

Grant ID: S19-0100

**National Institutes of Health**

Grant ID: HG002385

**Swedish Research Council**

Grant ID: 2021-03494

**Swedish Research Council**

Grant ID: 2020-01660

## Abstract

This protocol describes how to perform CRISPRi in hiPSCs and fbNPCs

## Troubleshooting



## gRNA design

- 1 To **design** the silencing of the expression of LINC01876 in hiPSCs, the protocol detailed in [Johansson et al., 2022](#) was adapted.

Single guide sequences were designed to recognise DNA regions near the transcription start site (TSS) of the chosen locus according to the [GPP Portal \(Broad Institute\)](#).

## CRISPR gRNA cloning and virus production

1h 30m

- 2 The guide sequences were **cloned** into a deadCas9-KRAB-T2A-GFP lentiviral backbone, pLV hU6-sgRNA hUbC-dCas9-KRAB-T2a-GFP, a gift from Charles Gersbach ([Addgene plasmid #71237 RRID:Addgene\\_71237](#)), using annealed oligos and the BsmBI cloning site.

To see the detailed protocol for the gRNA cloning, see [CRISPR gRNAs cloning protocol](#).

- 3 **Lentiviral vectors** were produced from the plasmids according to the protocol published by [Zufferey et al., 1997](#).

1h 30m

Briefly, HEK293T cells were grown to a confluency of 70 – 90% for lentiviral production. Third-generation packaging and envelope vectors (pMDL, psRev, and pMD2G) together with Polyethyleneimine (PEI Polysciences PN 23966, in DPBS (GIBCO) were used in conjunction with the lentiviral plasmids previously generated.

The lentivirus was harvested 2 days after transfection. The media was collected, filtered and centrifuged at 25000 x g, 4 °C for 01:30:00 . The supernatant was removed from the tubes and the virus was resuspended in PBS and left at 4 °C . The resulting lentivirus was aliquoted and stored at 80 °C .

## CRISPRi in iPSCs

- 4 The iPSCs were then **transduced** with MOI 10 of LacZ control and two LINC1876-targeting gRNAs, which led to three separate conditions:

1. LacZ control-,
2. LINC1876 g1-,
3. and LINC1876 g2- transduced cells.

- 5 The **gRNAs efficiency** was validated using standard quantitative real-time RT-PCR.


Briefly, total RNA was first extracted using the miniRNeasy kit (QIAGEN), and cDNA was generated using the Maxima First Strand cDNA Synthesis Kit (Thermo Scientific).

Quantitative qPCR was performed using SYBR Green I master (Roche) on a LightCycler 480 (Roche). The  $2^{-\Delta\Delta C_t}$  method was used to normalise expression to control, relative to housekeeping genes GAPDH and B-ACTIN expression.

#### Note

Gene Primer sequences (5' to 3'):


- LINC01876 Fw AATCCGTGCCAGCAGTAAGT Rev GGACCTCTTCAAGTCCCAGG
- ACTB Fw CCTTGCACA TGCCGGAG Rev GCACAGAGCCTCGCCTT
- GAPDH Fw TTGAGGTCAARGAAGGGGTC Rev GAAGGTGAAGGTCGGAGTCA

- 6 To select the successfully transduced cells expressing the vector we leveraged the GFP gene expressed by the employed vector as follows. 

7 days post transduction, cells were detached with Accutase (75 ml/cm<sup>2</sup>; GIBCO), resuspended in **iPS media** (StemMACS iPS-Brew XF and 0.5% penicillin/streptomycin (GIBCO)) containing RY27632 (10 µM) and Draq7 (1:1000) and strained with a 70µm filter. Gating parameters were determined by side and forward scatter to eliminate debris and aggregated cells. The GFP-positive gates were set using untransduced iPSCs. The sorting gates and strategies were validated via reanalysis of sorted cells (> 95% purity cut-off). 200.000 GFP-positive/Draq7-negative cells were collected per sample, spun down at 400g for 5 min and resuspended in iPS media containing RY27632 (10 µM) and either expanded or frozen down for further use.

## CRISPRi in neural progenitor cells (NPCs)

- 7 The selected iPSCs were differentiated into NPCs according to the differentiation protocol described in **Forebrain neural progenitor cells differentiation from iPSCs.**

- 8 NPCs that were efficiently expressing the vector of interest were selected as follows. 

At day 14 cells were detached with Accutase (75 ml/cm<sup>2</sup>; GIBCO), resuspended in **B27 media** (Neurobasal supplemented with 1% B27 without vitamin A (GIBCO), 2 mM L-glutamine and 0.2% penicillin/streptomycin Y27632 (10 µM), BDNF (20 ng/ml; R&D), and L-ascorbic acid (0.2 mM; Sigma) containing RY27632 (10 µM) and Draq7 (1:1000, BD Bioscience) and strained with a 70µm (BD Bioscience) filter. Gating parameters were determined by side and forward scatter to eliminate debris and aggregated cells. The GFP-positive gates were set using untransduced fbNPCs. The sorting gates and strategies were validated via reanalysis of sorted cells (> 95% purity cut-off).



200.000 GFP-positive/Draq7-negative cells were collected per sample, spun down at 400g for 5 min and snap frozen on dry ice. Cell pellets were kept at  $-80^{\circ}\text{C}$  until RNA was isolated.