



Mar 13, 2025

CRISPRi Knockdown in NGN2 iPSCs

DOI

<https://dx.doi.org/10.17504/protocols.io.rm7vzk268vx1/v1>

Alexis Penverne¹

¹Gandhi Lab



Alexis Penverne

Gandhi Lab, University College London, University of London

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vzk268vx1/v1>

Protocol Citation: Alexis Penverne 2025. CRISPRi Knockdown in NGN2 iPSCs. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.rm7vzk268vx1/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 22, 2025



Last Modified: March 13, 2025

Protocol Integer ID: 118874

Keywords: crispri knockdown in ngn2 ipscs crispri knockdown, ngn2 ipscs crispri knockdown, ngn2 ipscs with dcas9, using sgrna, ngn2 ipsc, crispri knockdown, dcas9, ipsc

Abstract

CRISPRi Knockdown in NGN2 iPSCs with dCas9 using sgRNAs

- 1 Doxycycline-inducible NGN2 iPSCs are maintained on Geltrex in mTeSR (STEMCELL) , and passaged using  0.5 micromolar (μM) EDTA. iPSCs should be between 50-70% confluent to increase effectiveness of transduction.
- 2 Prepare lentiviral stocks with sgRNAs targeting the genes for Knockdown.
- 3 Dilute aliquot of lentivirus preparation in enough media to cover cells in 6-well plate. We diluted our stocks 1:100 in 1mL of mTeSR.

Note

Using 1mL of media to save lentiviral stock, we transduced the cells overnight and changed the media early the next day to ensure cell survival.

- 4 Change media with  2 mL of mTeSR per well of a 6-well plate. Leave the cells to recover for 2 days, changing the media daily.
- 5 Once cells look healthy after transduction step, add 0.5ug/mL of Puromycin to kill untransduced cells. Keep puromycin as long as control untransduced cells are not dead.

Note

Include antibiotic to which lentiviral construct presents the relevant resistance. Keep one well of untransduced cells to control whether the antibiotic accurately gets rid of cells that haven't taken up the plasmid. It is possible to increase the antibiotic concentration to 1ug/mL if cells are staying alive.