



Apr 25, 2025

Version 2

🌐 Crystallisation of Enterovirus coxsackievirus A16 2A protease in space group C2 used for apo, fragment screen and follow compounds (PDB 8POA) V.2



DOI

<https://dx.doi.org/10.17504/protocols.io.3byl49kdzgo5/v2>Ryan Lithgo^{1,2,3}, Peter Marples^{1,2,3}, Lizbé Koekemoer^{4,5,3}, Daren Fearon^{1,2,3}¹Diamond Light Source; ²Research Complex at Harwell; ³ASAP Discovery Consortium;⁴Centre of Medicines Discovery; ⁵University of Oxford

Ryan Lithgo: The principle crystallographer on the Coxsackievirus A16 project.;

ASAP Discovery

OpenBind Consortium

**Lizbé Koekemoer**

University of Oxford

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

[Create free account](#)OPEN  ACCESSDOI: <https://dx.doi.org/10.17504/protocols.io.3byl49kdzgo5/v2>

Protocol Citation: Ryan Lithgo, Peter Marples, Lizbé Koekemoer, Daren Fearon 2025. Crystallisation of Enterovirus coxsackievirus A16 2A protease in space group C2 used for apo, fragment screen and follow compounds (PDB 8POA).

protocols.io <https://dx.doi.org/10.17504/protocols.io.3byl49kdzgo5/v2> Version created by **Mary-Ann Xavier**

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 25, 2025

Last Modified: April 25, 2025

Protocol Integer ID: 162007

Keywords: crystallisation, XChem, ASAP, AViDD, CMD, Diamond Light Source, i04-1, Coxsackievirus, A16, crystallisation of enterovirus coxsackievirus a16, coxsackievirus a16 crystal, 2a protease of the virus, enterovirus coxsackievirus a16, picornaviridae coxsackievirus a16, assembly of capsid protein, final stages of viral replication, cleavage from the poly protein, capsid protein, antiviral activity, virus, crystallographic fragment screening, poly protein, throughput crystallographic fragment screening, viral replication, protein production protocol, pdb 8poa, 2a protease in space group c2, 2a protease, formation of mature virion, protein production, protease, crystallisation, target for pandemic preparedness, solvent, protein

Funders Acknowledgements:

National Institutes of Health/National Institute Of Allergy and Infectious Diseases (NIH/NIAID)

Grant ID: Grant ID: U19AI171399

Disclaimer

The content is solely the responsibility of the authors and does not necessarily represent the official views of the National Institutes of Health.

Acknowledgements:

Diamond Light Source Ltd, Harwell Science and Innovation Campus, Didcot OX11 0QX, UK
Research Complex at Harwell, Harwell Science and Innovation Campus, Didcot OX11 0FA, UK
Oxford Lab Technologies crystal shifter <https://doi.org/10.1107/S2059798320014114>

Abstract

Picornaviridae coxsackievirus A16 is the causative agent of paediatric hand-foot-and-mouth disease, and a target for pandemic preparedness due to the risk of higher order complications in a large-scale outbreak. The 2A protease of the virus is responsible for self-cleavage from the poly protein, allowing for correct folding and assembly of capsid proteins in the final stages of viral replication. Inhibition deranges capsid folding and assembly, preventing formation of mature virions in host cells and making the protease a valuable target for antiviral activity. This protocol was used to grow coxsackievirus A16 crystals (PDB 8POA) that were used in high-throughput crystallographic fragment screening, and follow up compounds on the target. In this new version we added: the group deposition code; details about the fragment screen and solvent tolerance; also the protein production protocol.

Materials

<https://swissci.com/product/3-lens-crystallisation-plate/> Codes:

Midi: UVXPO-3LENS 3W96T-PS 3W96T-UVP

[M] 1 Molarity (M) MES  6.7 , Molecular Dimensions, Catalog # MD2-013-PH 6.7
50% w/v PEG 20000, Molecular Dimensions, Catalog # MD2-250-16

Purified SARS CoV-2 Coxsackievirus A16 protein ([M] 20 mg/mL) in [M] 10 millimolar (mM) HEPES,  7.5 ,
[M] 0.5 Molarity (M) NaCl, 5% glycerol, [M] 0.5 millimolar (mM) TCEP

Protein construct <https://www.addgene.org/204809/>

Protocol materials

 Human Coxsackievirus A16 strain G10 2A protease addgene Catalog #228632

Safety warnings

 Follow all handling warning for the chemicals used in the crystallisation screen composition.



Protein expression and purification

- 1 We used the following protein expression and purification protocol for the plasmid

 Human Cocksackievirus A16 strain G10 2A protease **addgene Catalog #228632**

Protocol

NAME

Enterovirus coxsackievirus A16 2A protease small scale expression and purification protocol

CREATED BY

korvus.wang

Preview

Equipment needed

- 2 Formulatrix Rock Imager (or incubator of choice)
SPT mosquito

P100 8 multi-channel pipette

SwissCI 3 lens plate

Crystallization experiment

1d

- 3 Prepare seed stock:



Protocol



NAME

Diamond XChem Seeding Protocol

CREATED BY

Peter Marples

[Preview](#)

1: 1000 dilution Sample seeds

4 Protein and buffer requirements: 43.2 μ L [M] 20 mg/mL Sample

3.36 mL Crystallization screen

 14.4 μ L Sample seeds, dilution 1:1000**5 Crystallisation screen composition:**

13.5 % PEG 20000

[M] 0.1 Molarity (M) MES 6.7

Stock solutions used:

[M] 1 Molarity (M) MES 6.7

50% w/v PEG 20000

Note

The crystallisation screen can be stored in a duran bottle or aliquoted into 96 deep well block for easy dispensing into SwissCI 3 lens plates.

For long term storage keep the Crystallisation screen in the fridge at 4°C.

6 Dispense 35 μ L Crystallisation screen into SwissCI 3 lens plate reservoir wells using a 100 μ l multi-channel pipette.

Dispense 150 μ L [M] 20 mg/mL Sample to each lens using the SPT mosquito.Dispense 150 μ L Crystallisation screen to each lens using the SPT mosquito.Dispense 50 μ L Seeds to each lens using the SPT mosquito.



Drop ratio: 3:3:1 ratio (150 nl  Sample : 150 nl reservoir solution: 50 nl seeds)

Final drop volume: 350 nl

- 7 Incubate at  20 °C for  24:00:00 h in Formulatrix Rock Imager.

1d

Imaging Schedule: The first images are taken after 12hrs and the imaging schedule follows a Fibonacci sequence of days for further collections.

- 8 Crystal typically form after ~24hrs

Expected result

Crystals typically reach their maximum size after ~24 h.


Morphology: typically rectangles.

Size: ~75 μm in length and ~10 μm in width, depth of the crystals is ~10 μm , giving a rectangular appearance

Average resolution: 1.6 $\text{\AA}$

Space group: C2

Unit cell: 86 $\text{\AA}$, 57 $\text{\AA}$, 32 $\text{\AA}$
90°, 95°, 90°



An example of a drop containing Cocksackievirus A16 crystals.



Data collection at Synchrotron

9 Diamond Light Source
Unattended Data Collection (UDC)
Data Collection Temperature: 100K
Detector: DECTRIS EIGER2 X 9M
Beamline: I04-1
Wavelength: 0.9212 Å
Resolution (Å): 1.21
Beam Size (μm): 60 X 50
Number of images: 3600
Oscillation: 0.10°
Exposure (s): 0.0020
Transmission (%): 100
Flux (ph/s): 9.50e+11

Solvent tolerance

10 Best conditions were: 10%DMSO, for 1h.

Protocol



NAME

XChem solvent test

CREATED BY

Peter Marples

Preview

Fragment and compound soaking

11 Conditions used 5-10% DMSO, for 1h



Protocol



NAME

XChem crystallographic fragment screening

CREATED BY

Peter Marples

Preview

Results

- 12 The results of this fragment screen was deposited under the PDB code G_1002288
https://www.rcsb.org/groups/summary/entry/G_1002288

Protocol references

https://www.rcsb.org/groups/summary/entry/G_1002288

<https://www.rcsb.org/structure/8POA>

<https://doi.org/10.1101/2024.04.29.591684>