

Oct 01, 2020

CTAB extraction of DNA and RNA of respiratory samples for microbial work

 In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.bf28jqhw

Leah Cuthbertson¹, Michael Cox¹, Phillip James¹, Colin P Churchward¹, Bill Cookson¹, Miriam Moffatt¹

¹National Heart and Lung Institute, Imperial College, London, UK, SW3 6LY



Leah Cuthbertson

University of Leicester



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bf28jqhw>

Protocol Citation: Leah Cuthbertson, Michael Cox, Phillip James, Colin P Churchward, Bill Cookson, Miriam Moffatt 2020. CTAB extraction of DNA and RNA of respiratory samples for microbial work. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bf28jqhw>



Manuscript citation:

Cowman SA, James P, Wilson R, Cookson WOC, Moffatt MF, Loebinger MR (2018) Profiling mycobacterial communities in pulmonary nontuberculous mycobacterial disease. PloS one. <https://doi.org/10.1371/journal.pone.0208018>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 06, 2020

Last Modified: October 01, 2020

Protocol Integer ID: 36672

Keywords: DNA, Extraction, respiratory, microbial, microbiome, metagenomics, Respiratory, upper respiratory tract swab, ctab extraction of dna, rna of respiratory sample, ctab extraction, throat swab, microbial work this extraction method, lower respiratory tract sample, comparison with dna extraction kit, respiratory sample, dna extraction kit, extracting dna, upper respiratory tract lavage, dna from cultured isolate, bronchoalveolar lavage, cultured isolate, lower respiratory tract, extraction, resulting dna pellet, lower respiratory tract synthetic absorptive matrix, extraction method, nasal lavage, nucleic acid yield, saliva sputum, bronchoscopy, precipitation of nucleic acid, dna pellet, nasopharyngeal aspirate, microbial work, mouth wash, nucleic acid, strain, ctab, isoamyl alcohol to each sample, double extraction protocol,



Abstract

This extraction method is modified from Griffiths et al 2000 Applied Environmental Microbiology 66(12): 5488–5491 and DeAngelis et al 2009 ISME Journal 3 pp 168–178.

Citation

Griffiths RI, Whiteley AS, O'Donnell AG, Bailey MJ (2000)

. Rapid method for coextraction of DNA and RNA from natural environments for analysis of ribosomal DNA- and rRNA-based microbial community composition..

Applied and environmental microbiology.

Citation

DeAngelis KM, Brodie EL, DeSantis TZ, Andersen GL, Lindow SE, Firestone MK (2009)

. Selective progressive response of soil microbial community to wild oat roots.. The ISME journal.

<https://doi.org/10.1038/ismej.2008.103>

LINK

In comparison with DNA extraction kits, nucleic acid yields tend to be higher, though it is more laborious for larger numbers of samples. The extraction uses CTAB and Phenol:Chloroform:Isoamyl alcohol for lysis and the precipitation of nucleic acid is PEG based. An alternative to bead-beating is described, which should only be used for extracting DNA from cultured isolates. Bead-beating is avoided in this case in order to minimise shearing of DNA prior to genomic DNA sequencing. For all other samples use bead-beating is recommended.

This Protocol will cover the extraction of the following sample types;

Bacterial isolates

Upper respiratory tract swabs (nasal and throat swabs)

Upper respiratory tract lavage and aspirates (nasopharyngeal aspirate, endotracheal aspirate, nasal lavage, Mouth wash)

Saliva

Sputum (expectorated spontaneously produced, and induced sputum)

Upper and lower respiratory tract synthetic absorptive matrix (SAM) strips



Lower respiratory tract samples obtained via bronchoscopy - lung brushes, bronchoalveolar lavage and pleural fluid

Lung tissue – human or mouse
Lung biopsies
Cerebral spinal fluid (CSF)



This is a double extraction protocol, you are adding a second volume of Phenol:Chloroform:Isoamyl alcohol to each sample, precipitating two aqueous phases per sample and then recombining the resulting DNA pellets. This approach significantly increases the yield of DNA from samples and strains.

Materials

MATERIALS

-  Phenol:Chloroform:Isoamyl alcohol (25:24:1) Saturated with 10 mM Tris, pH 8.0
-  Hexadecyltrimethylammonium bromide **Merck MilliporeSigma (Sigma-Aldrich) Catalog #H6269**
-  Ethanol (molecular biology grade, ≥99.8%) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #51976-500ML-F**
-  Sodium Chloride **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S9888**
-  Chloroform:Isoamyl alcohol 24:1 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #C0549**
-  lysing matrix E **Catalog #116914050**
-  Phase Lock Gel Separation tube Heavy QuantaBio **VWR International (Avantor) Catalog #733-2478**
-  FastPrep-24 Homogenizer **MP Biomedicals Catalog #116004500**
-  Linear Polyacrylamide **Merck MilliporeSigma (Sigma-Aldrich) Catalog #56575-1ml**
-  Aluminum ammonium sulfate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #402816**
-  Polyethylene glycol 6000 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #8074911000**
-  0.5ml Plain Skirted Tube **StarLab Catalog #E1405-2142**

STEP MATERIALS

-  0.5ml Plain Skirted Tube **StarLab Catalog #E1405-2142**
-  Phase Lock Gel Separation tube Heavy QuantaBio **VWR International (Avantor) Catalog #733-2478**

Protocol materials

- ✕ Polyethylene glycol 6000 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #8074911000**
- ✕ Sodium Chloride **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S9888**
- ✕ lysing matrix E **Catalog #116914050**
- ✕ Linear Polyacrylamide **Merck MilliporeSigma (Sigma-Aldrich) Catalog #56575-1ml**
- ✕ 0.5ml Plain Skirted Tube **StarLab Catalog #E1405-2142**
- ✕ 0.5ml Plain Skirted Tube **StarLab Catalog #E1405-2142**
- ✕ Phase Lock Gel Separation tube Heavy QuantaBio **VWR International (Avantor) Catalog #733-2478**
- ✕ Ethanol (molecular biology grade, ≥99.8%) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #51976-500ML-F**
- ✕ Chloroform:Isoamyl alcohol 24:1 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #C0549**
- ✕ Phase Lock Gel Separation tube Heavy QuantaBio **VWR International (Avantor) Catalog #733-2478**
- ✕ FastPrep-24 Homogenizer **MP Biomedicals Catalog #116004500**
- ✕ Hexadecyltrimethylammonium bromide **Merck MilliporeSigma (Sigma-Aldrich) Catalog #H6269**
- ✕ Aluminum ammonium sulfate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #402816**
- ✕ Phenol:Chloroform:Isoamyl alcohol (25:24:1) Saturated with 10 mM Tris, pH 8.0
- ✕ 0.5ml Plain Skirted Tube **StarLab Catalog #E1405-2142**
- ✕ Phase Lock Gel Separation tube Heavy QuantaBio **VWR International (Avantor) Catalog #733-2478**

Troubleshooting

Safety warnings



- Informal training by a competent person before task can be undertaken unsupervised
- Appropriate PPE must be worn at all time
- Recommend doubling up nitrile gloves for handling phenol
- Work conducted at Containment Level 2
- Primary samples may be handled, but phenol and chloroform and to a lesser extent CTAB are the greater risks, so perform **work in the fume cupboard**, rather than the biological safety cabinet. Do not handle samples that are suspected to contain pathogens above ACDP category 2.
- Ensure correct waste disposal procedures are used. Waste solvents should be collected in labelled glass waste.
- The CTAB buffer cannot be disposed of down the sink as it is an environmental hazard. Collect waste solution as above.
- Note, when pipetting solvents, pre-wet the pipette tip by very gently pipetting up and down. This saturates the headspace within the tip with the volatile solvent and prevents vapour pressure from causing the tip to drip excessively when pipetting.
- Sealed buckets/rotors should be used for centrifugation steps
- Sample spillages should be disinfected with 70% ethanol followed by cleaning with Surfanios (Biolab, 059840)

Ensure clear and appropriate labelling of all stored samples



Before start

Prepare reagents:

CTAB extraction buffer:

1Part: 10% (w/v) CTAB (hexadecyltrimethylammonium bromide) in 1M NaCl

1 part: 0.5 M phosphate buffer (pH 8.0) in 1 M NaCl

Prepare 1M 1 Molarity (M) NaCl in milliQ water.

Use the 1M 1 Molarity (M) NaCl solution you have prepared in place of water for preparing the 10 % CTAB stock solution.

Prepare the phosphate buffer at pH 8.0.

- Prepare 1M 1 Molarity (M) NaH_2PO_4 (Monobasic phosphate) and 1M 1 Molarity (M) Na_2HPO_4 (Dibasic phosphate) stocks using 1M 1 Molarity (M) NaCl in place of water. They may require gentle heating and stirring in order to dissolve.

- Combine 15.9 mL of monobasic phosphate with 284.1 mL of dibasic phosphate and make up to 600 mL with 1M 1 Molarity (M) NaCl to achieve pH 8.0

Combine the phosphate buffer and CTAB solution 1:1 to complete the CTAB extraction buffer

Sterilise by autoclaving

PEG/NaCl precipitation solution:

30 % (w/v) polyethylene glycol 6000 in 1.6 Molarity (M) NaCl

Sterilise by autoclaving


0.1M Aluminum ammonium sulphate:

0.1 Molarity (M) Aluminum ammonium sulfate ($\text{AlNH}_4(\text{SO}_4)2 \cdot 12\text{H}_2\text{O}$, Sigma 402816)

Filter-sterilize through 0.2mm filter






Before you start

- 1 Add working stock of Phenol:Chloroform:Isoamyl (P:C:I (25:24:1)) alcohol to a  50 mL Falcon tubes in the fume cupboard.



Note

Use a 10ml pipette to remove required volume of P:C:I from below the the top buffered layer in the stock bottle.
*For a new bottle ensure the buffer has been added and allowed to settle prior to starting the extraction.


- 2 Add working stock of Chloroform:Isoamyl (C:I (24:1)) alcohol to a  50 mL Falcon tubes in the fume cupboard.
- 3 Pre-spin the 2 Phase-lock gel tubes ( 2 mL Hard gel) per samples .Phase lock gel should be pelleted at the bottom of each tube (often on the sides when they arrive from manufacturer). Centrifuge at 16,000 x g for  00:05:00 .



Phase Lock Gel Separation tube Heavy QuantaBio **VWR International**
(Avantor) Catalog #733-2478

- 4 Add  1 μ L of linear polyacrylamide (LPA) to  2 mL eppendorfs, 2 per sample.This is a carrier and precipitates along with DNA increasing the yield. Unlike glycogen which can also be used it does not affect later sample use.

Extraction


- 5 Add  50 μ L aluminium ammonium sulfate to each Lysing Matrix E (LME) tube

Note





It is helpful to write the sample numbers on the lid and side of the LME tubes as the beadbeater can rub the numbers off.






6

Aseptically transfer sample to the LME tube and add 500ul of CTAB to the LME tube and incubate for  00:15:00 .




**Note****Bacterial Isolates**

- Grow up bacterial isolates overnight in appropriate culture broth
- Spin down  2 mL of broth for  00:10:00 at 16,000 x g
- In a safety cabinet add 500ul of CTAB and resuspend the pellet. Incubate at room temperature for  00:15:00 .
- Bacteria and CTAB can be stored frozen prior to extraction or extracted after  00:15:00 incubation. **No additional CTAB is required.**

Upper respiratory tract swabs (nasal and throat swabs)

- Prepare 1 sterile spin baskets in a  1.5 mL trefflab tube per sample.
- 1 pair of autoclaved scissors will be required per swab to transfer swab tips into bead beating tubes
- Keep swabs on ice until transferred into LME tubes
- Aseptically transfer swab tips into the LME tubes using sterile scissors.
- Add  500 µL of CTAB to the LME tube and incubate for  00:15:00 .

Upper respiratory tract lavage and aspirates (nasopharyngeal aspirate, endotracheal aspirate, nasal lavage)


- Spin a minimum volume of  2 mL , optimal volume 5mls of lavage fluid for  00:20:00 at full speed
- Carefully resuspend pellet in  500 µL CTAB extraction buffer and transfer to LME tube.

NOTE:For saliva or mouth wash samples and sputum samples pre-alequoting using wide bore tips prior to storage is recommended.

Saliva or Mouth wash samples

- Transfer a known volume of no more than  300 µL into a LME tube.

Sputum (expectorated spontaneously produced, and induced sputum)

- Transfer a known volume of no more than  300 µg or individual sputum plug, into a LME tube.





Upper and lower respiratory tract synthetic absorptive matrix (SAM) strips

- 1 pair of autoclaved scissors will be required per SAM to transfer SAM into LME tubes
- Keep SAMs on ice until transferred into LME tubes

Lung brushes

- Keep lung brushes on ice until transferred to LME tubes

Bronchoalveolar lavage or Pleural fluid

- Spin a minimum volume of  2 mL , optimal volume  5 mL of lavage fluid for  00:20:00 at full speed
- Carefully resuspend pellet in  500 µL CTAB extraction buffer before transferring to LME tube. **NO further CTAB is required**



**Lung tissue- Human or Mouse**

- a. Keep tissue on ice until transferred to LME tubes

Lung biopsies


- a. Keep tissue on ice until transferred to LME tubes.

CSF

- a. Maximise and standardise sample volume across the study.
- b. Spin fluid for  00:20:00 at 16,000 x g
- c. Carefully suspend pellet in  500 µL CTAB extraction buffer before transferring to LME tube. **NO further CTAB is required**

Safety information

Transfer of patient samples should be performed in a class 2 safety hood.

- 7 Moving to a fume cupboard, carefully and immediately add  500 µL of P:C:I (25:24:1).

Note

Pre-wet the pipette to avoid drips due to vapour pressure in tip.

Safety information

Perform this step in a fume hood


Safety information

Wear double gloves when handling P:C:I (25:24:1). Should a small spillage occur, outer gloves can then be disposed of without risk.



Safety information

DO NOT leave tubes with P:C:I (25:24:1) for extended periods, it can degrade the plastic

- 8 Transfer to the bead-beater and beat using the pre-programmed CTAB setting for  00:01:00 . Return tubes to ice immediately after beating.


Note

CTAB Bead beating Program:

Speed: 5.5m/sec
Adapter: Quickprep
Time: 60 sec
Lysing Matrix: E
Quantity: 1ml
Cycles: 1
Pause time: 300sec


Note



Ensure the lids of the LME tubes are securely fastened with no beads in the seal, and that the tubes are labelled on the top and on the side as the

- 9 Centrifuge LME tubes at 16,000 x g for  00:05:00 .

- 10 Transfer all liquid to phase lock gel tube and keep the LME tube on ice.



Note

Phase lock tubes should have been pre-spun. Phase lock gel should be pelleted at the bottom of each tube (often on the sides when they arrive from manufacturer). Centrifuge at 16,000 x g for  00:05:00

- 11 Centrifuge the phase lock tube at 16,000 x g for  00:05:00 at  4 °C .




Note

Gel will form a barrier between the aqueous and P:C:I (25:24:1)/C:I (24:1) phases.



- 12 Add 1 volume of C:I (24:1) to each phase lock gel tube, shake briefly to mix. Centrifuge at 16,000 x g for  00:05:00 at  4 °C .

Note

If barrier does not form, extend centrifugation time.



- 13 **Second extraction;** In the fume hood, add  50 µL of Aluminium ammonium sulphate,  500 µL of CTAB extraction buffer and  500 µL of P:C:I (25:24:1) to each bead beating tube.
- 14 Repeat steps 8 to 13 then continue with the precipitation step below.

Precipitation

- 15 Transfer the aqueous phase from each tube to the pre-prepared eppendorfs with  1 µL of LPA in them. As the solvent is locked beneath the phase lock gel, you should be able to pour this. You may need 2 eppendorfs per sample to fit allow room for PEG.
- 16 Add 2 volumes of PEG/NaCl solution and mix well (solution is viscous). Leave overnight in the fridge at  4 °C to precipitate.

**Note**



If you are in a hurry it can be left for 2 hours but this will reduce the yield.

17 Centrifuge all eppendorfs at 16,000 x g for  00:20:00 at  4 °C .





18 Carefully aspirate the PEG/NaCl solution from the pellets.

Note

Pellets are usually large and translucent. They tend to be fairly uniform in size as the LPA has also precipitated.

19 Wash pellets with  500 µL ice-cold 70 % ethanol to remove any precipitated salts and centrifuge at 16,000 x g for  00:05:00 .

20 Repeat the wash twice with  200 µL ice-cold 70 % ethanol.

21 After briefly air drying the pellet for about 5 minutes, resuspend in  30 µL low EDTA TE and combine the resuspended total nucleic acid from each sample into a single tube ( 60 µL total per sample). For sputum samples (or other high biomass samples) you may need to increase this volume to  50 µL per tube ( 100 µL per sample).

Note

We do not recommend storing DNA in cryo tubes as you are unable to spin them down.
Recommended tubes: Tethered O ring, sterile tubes E1405-2142, star labs

 0.5ml Plain Skirted Tube **StarLab Catalog #E1405-2142**

22 This extract can be stored at  -20 °C or  -80 °C until required.



Citations

Griffiths RI, Whiteley AS, O'Donnell AG, Bailey MJ. Rapid method for coextraction of DNA and RNA from natural environments for analysis of ribosomal DNA- and rRNA-based microbial community composition.

DeAngelis KM, Brodie EL, DeSantis TZ, Andersen GL, Lindow SE, Firestone MK. Selective progressive response of soil microbial community to wild oat roots.

<https://doi.org/10.1038/ismej.2008.103>