



Apr 15, 2025

Dengue Virus Serotype 2 NS2B-NS3 protease fusion protein Avi-tagged (Biotinylated) protein expression and purification: large scale 1 L cultures

Forked from a private protocol

DOI


<https://dx.doi.org/10.17504/protocols.io.rm7vzkzzxvx1/v1>Michael Fairhead¹¹University of Oxford, ASAP Discovery Consortium

ASAP Discovery



Mary-Ann Xavier

Diamond Light Source, Centre for Medicines Discovery, ASAP D...



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

[Create free account](#)OPEN  ACCESSDOI: <https://dx.doi.org/10.17504/protocols.io.rm7vzkzzxvx1/v1>External link: <https://orcid.org/0000-0001-5361-3933>

Protocol Citation: Michael Fairhead 2025. Dengue Virus Serotype 2 NS2B-NS3 protease fusion protein Avi-tagged (Biotinylated) protein expression and purification: large scale 1 L cultures. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.rm7vzkzzxvx1/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 21, 2025

Last Modified: April 15, 2025

Protocol Integer ID: 123198

Keywords: parallel protein purification, Recombinant protein, Escherichia coli, ZIKV, Zika, NS2B-NS3 protease, purification of dengue virus serotype, ns3 protease fusion protein avi, dengue virus serotype, ns3 protease fusion protein, recombinant protein, purification, protein expression, protein, biotinylation tag, avitag, virus

Funders Acknowledgements:

National Institutes of Health/National Institute Of Allergy and Infectious Diseases (NIH/NIAID)

Grant ID: U19AI171399

Disclaimer

The content is solely the responsibility of the authors and does not necessarily represent the official views of the National Institutes of Health.

Abstract

This protocol details the expression and purification of Dengue Virus Serotype 2 NS2B-NS3 protease fusion protein with a AviTag (biotinylation tag) for biophysical assay at a 1 L culture scale. Recombinant proteins are expressed in *Escherichia coli* using the autoinduction method and then purified in using a IMAC, desalt, tag cleavage, reverse IMAC and gel filtration work flow.

Attachments



PAGE24-00124 -

AVIDD...

935KB

Guidelines

Method overview

Standard workflow is expression via autoinduction followed by purification using IMAC/PD-10/revIMAC and serial gel filtration

Materials

Dengue Virus Serotype 2 NS2B-NS3 protease fusion protein

Vector: pNIC (Kanamycin resistance)

Cell line: E. coli strain BL21(DE3)-RR + pCDF-BirA plasmid GenBank: JF914075.1 (Streptomycin resistance)

Tags and additions: N-terminal His6 Twin Strep TEV Avitag

Construct protein sequence before tag cleavage

MHHHHHHSSGASWSHPQFEKGGGSGGGSGGSAWSHPQFEKSGVDLGTENLYFQSGLNDIFEAQKIEWHEMADLELERA
ADVKWEDQAEISGSSPILSITISEDGSMSEKNEEEETLGGGGSGGGGAGVLWDVPSPPPMGKAELEDGAYRIKQKGILGYSQ
IGAGVYKEGTFHTMWHVTRGAVLMHKGKRIEPSWADVKKDLISYGGGWKLEGWEKEEVEQVLALPGKNPRAVQTKPGL
FKTNAGTIGAVSLDFSPGTSGSPIIDKKGKVVGLYGNGVVTRSGAYVSAIAQTEKSIEDNPEIEDDIFRK

- MW = 27.116 kDa
- Extinction co-efficient = $43430 \text{ M}^{-1} \text{ cm}^{-1}$
- pI = 5.67

Construct protein sequence after tag cleavage

SGLNDIFEAQKIEWHEMADLELERAADVKWEDQAEISGSSPILSITISEDGSMSEKNEEEETLGGGGSGGGGAGVLWDVPSP
PPMGKAELEDGAYRIKQKGILGYSQIGAGVYKEGTFHTMWHVTRGAVLMHKGKRIEPSWADVKKDLISYGGGWKLEGWEK
EGEEVQVLALPGKNPRAVQTKPGLFKTNAGTIGAVSLDFSPGTSGSPIIDKKGKVVGLYGNGVVTRSGAYVSAIAQTEKSIED
NPEIEDDIFRK

- MW = 27.116 kDa, +Biotin MW = 27.886 kDa
- Extinction co-efficient = $47440 \text{ M}^{-1} \text{ cm}^{-1}$
- pI = 4.92

 PD-10 desalting columns packed with Sephadex G-25 resin **Cytiva Catalog #17085101**

 PD-10 Spin Adapter **Cytiva Catalog #28923245**

 LabMate PD-10 Buffer Reservoir **Cytiva Catalog #18321603**

 His GraviTrap **Cytiva Catalog #11003399**

 Super Broth **Formedium Catalog #SPB0102**

 AirOtop Enhanced Flask seals **Generon Catalog #899425**

 Nalgene®; Unwire®; Test Tube Racks: Resmer®; Manufacturing Technology, for 30mm tubes, white **Thermo Fisher Catalog #5970-0030**

 AIM – Terrific Broth Base including Trace elements **Formedium Catalog #AIMTB0210**

 Ultra Yield 2.5L Flask, Sterile **Generon Catalog #931136-B**

Optional but useful

 BENCHMIXER™ XL MULTI-TUBE VORTEXER Benchmark Scientific Catalog #BV1010

Materials (1 L cultures) for Expression:

- Plates with LB-agar+antibiotics
-  1 L of autoclaved autoinduction TB + 20 g/L glycerol + antibiotics
-  1 mL of 10 % Antifoam 204 (Sigma) in ethanol
-  2.5 L Ultra Yield flasks (fitted with loose foil cover**)

Materials (1 L cultures) for Purification:

- **1L of Base Buffer**

	A	B
	HEPES	10 mM
	Glycerol	5%
	NaCl	500 mM
	TCEP, pH 7.5	0.5 mM

-  100 mL of **IM1 3 Mass Percent** imidazole pH 7.5.
-  100 mL of 10 % Triton X-100 in water.
-  10 mL Lysozyme solution (100 x).
-  1 mL homemade benzonase (1000x). Maybe substituted for 10 mg/mL of commercial DNase I
- 2 x His GraviTrap column per litre of culture to be purified (Cytiva) fitted with LabMate extender (Cytiva) and PD-10 spin adapter (Cytiva) in 24 place Nalgene rack.
- 2 x PD-10 desalting column per litre of culture to be purified fitted with LabMate extender and PD-10 spin adapter in 24 place Nalgene rack.
- 2 × 50 mL centrifuge tubes per litre of culture to be purified in a 24 place Nalgene rack.

Protocol materials

 Dengue virus serotype 2 NS2B-NS3 protease fusion with non cleavable linker and non cleavable avi-tag **addgene Catalog #228649**



Safety warnings

- ❗ Triton x-100 is currently restricted for use in the EU and cannot be used without an exemption certificate REACH Annex XIV (Jan 2021). It can be readily substituted with IGEPAL CA-630 (which is likely to be subject to the same restrictions in the near future). Alternatives that also maybe used are Tergitol 15-S-9 or Tween-20 or octyl glucoside.



Expression

1d 19h

- 1 Transform *BL21 (DE3)* with the appropriate plasmid



Dengue virus serotype 2 NS2B-NS3 protease fusion with non cleavable linker and non cleavable avi-tag **addgene Catalog #228649**

and spread onto LB-agar plate + 50 µg/mL kanamycin + 50 µg/mL streptomycin and incubate Overnight 37 °C *.

- 2 Prepare a starter culture by using several colonies to inoculate 10 mL of superbroth + 1 % glucose+ 50 µg/mL kanamycin + 50 µg/mL streptomycin in a 50 mL tube and incubate at 250 rpm, 37°C, 16:00:00

16h

- 3 Use 10 mL started culture to inoculate 1 L AIM-TB + 50 µg/mL kanamycin + 50 µg/mL streptomycin + 0.01% Antifoam 204 in a 2.5 L Ultra Yield baffled flask (Thomson) fitted with an AirOtop enhanced flask seal (Thomson)

- 4 Grow 250 rpm, 37°C, 04:00:00 shaking.

4h

- 5 Add 1 mL 500 millimolar (mM) D-Biotin dissolved in DMSO

- 6 Gow 250 rpm, 18°C, 01:00:00 shaking.

1h

- 7 Add 1 mL 500 millimolar (mM) D-Biotin dissolved in DMSO

Add 1 mL 100 millimolar (mM) IPTG

- 8 Grow 250 rpm, 18°C, 22:00:00 shaking.

22h

- 9 Harvest at 4000 x g, 4°C, 00:20:00 .

- 10 Scrape out pellet and place in plastic polygrip bag and place in -80 °C freezer. Final wet cell weight is 45 g/L of culture



Cell lysis

3h 30m

- 11 Place the polygrip bag on a flat surface and smash cell pellet into small pieces and transfer the cell pellet into 500 mL beaker.
- 12 Add  3 mL Base Buffer/g cell pellet ( 10 millimolar (mM) HEPES,  500 millimolar (mM) NaCl, 5 % Glycerol,  0.5 millimolar (mM) TCEP,  7.5) +  0.5 mg /mL Lysozyme,  1 µg /mL Benzonase or  10 µg /mL DNase I, 1 % Triton X-100***,  30 millimolar (mM) Imidazole  0.5 millimolar (mM) D-Biotin



Note

***Triton X-100 is currently restricted for use in the EU and cannot be used without an exemption certificate REACH Annex XIV (Jan 2021). It can be readily substituted with IGEPAL CA-630 (which is likely to be subject to the same restrictions in the near future). Alternatives that could be used are Tergitol 15-S-9 or Tween-20 or octyl glucoside.

- 13 Use stripette to dissolve pellet and transfer  45 mL (maximum) to a 50 mL tube (4 tubes in total).
- 14 Leave  00:30:00  Room temperature .
- 15 Place tubes in  -80 °C freezer for 2h or overnight if preferred.
- 16 Thaw in  Room temperature water bath for  01:00:00 and mix.
- 17 Centrifuge at  4000 x g, 4°C, 01:00:00 .

30m



1h



1h



Purification

3h 30m




- 18 Apply supernatant from 2 × 50 mL tubes to 1 mL His GraviTrap column (Cytiva) fitted with a LabMate PD-10 Buffer Reservoir (Cytiva). His GraviTrap columns are pre-equilibrated in Base Buffer + [M] 30 millimolar (mM) Imidazole.
- 19 Wash with  20 mL Base Buffer + [M] 30 millimolar (mM) Imidazole
Repeat wash a total of 3 times
- 20 Slot His GraviTrap column into PD-10 desalting column (Cytiva) fitted with a LabMate PD-10 Buffer Reservoir. PD-10 desalting columns are pre-equilibrated in Base Buffer w/o Imidazole.
- 21 Elute protein, with  2.5 mL of Base Buffer + [M] 500 millimolar (mM) Imidazole, directly onto PD-10 column.
- 22 Remove His GraviTrap column.
- 23 Place PD-10 desalting column into 50 mL falcon tube and add  3.5 mL Base Buffer w/o Imidazole and collect the flow through from the column.
- 24 Run TEV digested sample over His GraviTrap column to remove protease, tag and un-cleaved protein
Wash column with 2.5 mL Base Buffer + [M] 30 millimolar (mM) Imidazole
- 25 Concentrate to 10 mg/mL using a 10,000 MWCO centrifugal filter (Amicon, Merck-Millipore) and inject the sample (5 mL) onto 125 mL supdex 75 pg column pre-equilibrated in Base Buffer. Collect 2 mL fractions using a flow rate at 1 mL/minute.
- 26 Pool the peak fraction(s) (e.g. 10-16 in the chromatogram below) and concentrate to 10 mg/mL using a 10,000 MWCO centrifugal filter (Amicon, Merck-Millipore)
Snap freeze protein using liquid nitrogen as 0.1 mL aliquots and store the sample at -80°C until use
Approximate final yield is 55 mg

Purification Results


3h 30m

26.1 Chromatogram using Base Buffer as the mobile phase




QQ01D2VNS2BA-c003 SEC profile on 125 mL Superdex 75 pg

26.2 SDS PAGE on 4-12% Bis-TRIS gel



QQ01-c003 SEC profile on 125 mL Superdex 75 pg SDS PAGE of samples during purification

27 Mass spectrometry of final sample to confirm biotinylation








QQ01-c003 Mass spectrometry showing the final sample has the expected MW for biotinylated protein

Column regeneration: PD-10

- 28 Wash PD-10 columns with 50 mL - 100 mL of Milli-Q water.
- 29 Store all columns in water at 4 °C . For long term storage, use 20 % Ethanol



Column regeneration: His GraviTrap

- 30 Wash the IMAC columns with  40 mL Milli-Q. 
- 31 Wash IMAC columns  10 mL 20 % Ethanol +  0.1 Mass Percent EDTA*. 
- *PUT NICKEL WASTE IN APPROPRIATE CONTAINER FOR DISPOSAL!
- 32 Wash IMAC columns  40 mL Milli-Q. 
- 33 Wash IMAC columns  10 mL  1 Mass Percent NaOH. 
- 34 Wash IMAC columns  40 mL Milli-Q. 
- 35 Wash IMAC columns  10 mL  1 Mass Percent Acetic Acid + 1 % Triton X-100. 
- 36 Wash IMAC columns  40 mL Milli-Q.
- 37 Wash IMAC columns with  0.5 mL  100 millimolar (mM) Nickel Sulfate +  20 millimolar (mM) Tris.HCl pH 8*. 
- Note

*PUT NICKEL WASTE IN APPROPRIATE CONTAINER FOR DISPOSAL!
- 38 Wash IMAC columns with  40 mL Milli-Q. 
- 39 Store all columns in water at  4 °C . For long term storage use 20 % Ethanol

Protocol references

dx.doi.org/10.17504/protocols.io.4r3l22jkxl1y/v1