



Aug 08, 2023

Version 1

Detection of arboviruses in salivary glands and midgut of mosquitoes vector V.1

DOI

dx.doi.org/10.17504/protocols.io.8epv5xjndg1b/v1

ISABELLA FERREIRA DA COSTA¹, CAROLINE PIRES SANTANA¹, Isabela Cinquini Junqueira¹, Juscelino Rodrigues², Christian Luz², Suleimy Marinho Fernandes³, Izabela Batista Melo³, Fabíola Fiaccadori³, Valéria Christina de Rezende Feres¹

¹Laboratory of Molecular Biology and Applied Technologies for Laboratory Diagnosis, Faculty of Pharmacy, Federal University of Goiás (UFG), Goiânia, GO, Brazil.;

²Laboratory of Invertebrate Pathology, Institute of Tropical Pathology and Public Health, Federal University of Goiás (UFG), Goiânia, GO, Brazil.;

³Laboratory of Virology and Cell Culture, Institute of Tropical Pathology and Public Health, Federal University of Goiás (UFG), Goiânia, GO, Brazil.

**ISABELLA FERREIRA DA COSTA**

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

[Create free account](#)OPEN  ACCESS

DOI: <https://dx.doi.org/10.17504/protocols.io.8epv5xjndg1b/v1>

Protocol Citation: ISABELLA FERREIRA DA COSTA, CAROLINE PIRES SANTANA, Isabela Cinquini Junqueira, Juscelino Rodrigues, Christian Luz, Suleimy Marinho Fernandes, Izabela Batista Melo, Fabíola Fiaccadori, Valéria Christina de Rezende Feres 2023. Detection of arboviruses in salivary glands and midgut of mosquitoes vector. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.8epv5xjndg1b/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 07, 2023

Last Modified: August 08, 2023

Protocol Integer ID: 86053


Keywords: Arbovirus, Mosquitoes, Salivary Glands, Midgut, Cell Culture, Detection, characterization of these arbovirus, mosquitoes vector dengue, arbovirus, circulation of these arbovirus, midgut of mosquitoes vector dengue, mosquitoes of the genus aede, experimental technical protocol in adult mosquito, adult mosquito, mosquito, chikungunya virus, aedes albopictus, viral isolation, aedes albopictus for amplification, zika, salivary gland, genus aede, viral replication, infections in human

Abstract

Dengue, Zika and chikungunya viruses are widely disseminated in Brazil and cause infections in humans, with similar clinical manifestations of varying intensity, with relevance in the context of public health. These diseases are mainly transmitted by mosquitoes of the genus *Aedes*. Therefore, monitoring the circulation of these arboviruses in vectors is a crucial step in epidemiological surveillance. In this protocol, adult mosquitoes were dissected to extract the salivary glands and midgut, which are structures directly associated with viral replication and transmission. These structures were inoculated into cells of the C6/36 lineage of *Aedes albopictus* for amplification and viral isolation. Our main objective was to develop and evaluate an experimental technical protocol in adult mosquitoes with the aim of enhancing the detection and characterization of these arboviruses.

Troubleshooting

Safety warnings

 Please use personal protective types of equipment (such as gloves, safety glasses, a lab coat and masks) during all experiment.

Before starting

- 1 Prepare pools of 2 to 10 live and engorged and non-engorged adult mosquito specimens, sorted by sex (male/female). Place the pools in 1.5 µl microtubes and store in a freezer at -80°C.

Mosquitoes dissection

- 2 Pinch one mosquito at a time and remove the legs.
- 3 Place the mosquito on a micro glass slide under the stereoscopic magnifying glass
- 4 Add a drop of sterile solution of phosphate buffer saline (PBS) 1% to the micro glass slide to create a conductive medium for the separation of the desired structures.
- 5 With a needle, immobilize the mosquito's thorax by pinning it against the micro glass slide.
- 6 Using tweezers, carefully remove the entire abdomen of the mosquito, ensuring that the midgut remains attached to the immobilized thorax.
- 7 Use a needle to separate and section the midgut from the thorax.
- 8 Transfer the midgut to microtubes already prepared with 200 µL of PBS solution and store them appropriately.
- 9 While the mosquito is still fixed on the slide, use another needle to separate the mosquito's head from the rest of the body. This is necessary because the salivary glands, located in the anterior portion of the thorax, may partially adhere to the head.
- 10 Gently press the head and thorax to extract the salivary glands.
- 11 Using a pipette, extract the salivary glands from the medium and transfer them to the respective microtube where the intestine has already been stored.

- 12 Each tube will be homogenized by repeated pipetting (60 seconds). Subsequently, vortexed for 60 seconds to disrupt the cells.

Cell culture

- 13 Maintain C6/36 cells (*Aedes albopictus*) in culture using Leibovitz medium (L-15 Sigma-Aldrich, St. Louis, MO, USA), supplemented with 10% fetal bovine serum (FBS) (Sigma-Aldrich, St. Louis, MO, USA) and 1% Penicillin/Streptomycin (Gibco, Waltham, MA, USA) at 28°C in an incubator.
- 14 Monitor cell growth and wait until the cells reach approximately 80% of confluence.
- 15 Completely discard the medium from the cell culture bottle and resuspend the cells in approximately 2 ml of L-15 medium.
- 16 Transfer the cells to a falcon tube (50 ml) and proceed to count them using a Neubauer chamber.
- 17 Prepare a solution combining the cell suspension and Trypan blue dye at a 1:4 dilution ratio.
- 18 Under an optical microscope, count all living cells present in the four lateral quadrants of the Neubauer chamber.
- 19 Calculate the average of the cells counted in the four quadrants. Then, multiply the result by the dilution factor and the Neubauer chamber factor (10^4).
- 20 Aliquot 500 μ L of the cell suspension into each well of a 24-well cell culture plate. Incubate the cell culture plate at 28°C for 48 hours.
- 21 After confirming the formation of a cell monolayer, inoculate 50 μ L of the samples processed as described in item 2. Then, add 500 μ L of L-15 medium to each well and incubate the plate for 7 days.
- 22 At the 7th day, freeze the plates at -80°C for at least overnight. Then, thaw the plates at room temperature to promote cell disruption and virus release.
- 23 Aspirate the contents of each well individually and transfer them to centrifuge tubes.



- 24 Centrifuge the tube at 14.000 rpm for 60 seconds at room temperature. After centrifugation, separate the supernatant and refrigerate it at -80°C for subsequent viral detection and identification through molecular biology techniques.
- 25 All stages of cell cultivation must be carried out in a sterile biological safety cabinet and in accordance with biosafety standards.

Arbovirus Molecular identification

- 26 Extract the viral RNA using the silica/gel membrane column methods, utilizing a commercial kit (QIAamp Viral RNA Mini Kit, Qiagen, Hilden, Germany).
- 27 Perform reverse transcription and amplify the genetic material using real-time polymerase chain reaction (PCR).
- 28 This protocol was standardized using two methodologies for real-time RT-PCR. The protocol described by Huhtamo et al. 2010 for Dengue virus (DENV) detection, which employed specific primers and a probe for the conserved 3'UTR region (Path-IDtm multiplex One-Step RT-qPCR - Applied Biosystems – Carlsbad, CA, USA). The final reaction volume was 10 µL, and the components included 2.5 µL of extracted material, 5.0 µL of Multiplex RT-PCR buffer, 1.0 µL of Multiplex Enzyme Mix, 100 nM of each primer, 200 nM of TaqMan probe and Nuclease-Free water provided in the kit to complete the volume.
- 29 The protocol described by Lanciotti et al. 2008 for Zika virus (ZIKV) detection was performed using specific primers and a probe targeting the conserved region of the protein E gene (GoTaq® Probe 1-Step RT-qPCR System, Promega – USA). The final reaction volume was 10 µL, and the components included 2.5 µL of extracted material, 5.0 µL of GoTaq® Probe Master Mix with dUTP, 0.4 µL of GoScript™ RT Mix, 100 nM of each primer, 200 nM of TaqMan probe and, Nuclease-Free water provided in the kit to complete the volume.
- 30 For the detection of nucleic acids (RNA) from ZIKV, DENV (4 serotypes) and Chikungunya virus, the commercial Kit IBMP Biomol ZDC-Zika, Dengue and Chikungunya (IBMP, Curitiba - PR - Brazil) was utilized, following the manufacturer's instructions.

Protocol references

Coleman J, Juhn J, James AA. Dissection of midgut and salivary glands from *Ae. aegypti* mosquitoes. *J Vis Exp.* 2007;(5):228. doi: 10.3791/228. Epub 2007 Jul 4. PMID: 18979026; PMCID: PMC2557086.