



May 07, 2025

Dinoflagellate Single Cell Transcriptome

DOI

<https://dx.doi.org/10.17504/protocols.io.261ge8dpjg47/v1>

Brittany Ott¹

¹US FDA



Brittany Ott

US FDA

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261ge8dpjg47/v1>

Protocol Citation: Brittany Ott 2025. Dinoflagellate Single Cell Transcriptome. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.261ge8dpjg47/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 02, 2025



Last Modified: May 07, 2025

Protocol Integer ID: 210830

Keywords: dinoflagellate single cell transcriptome this protocol, dinoflagellate single cell transcriptome, transcriptomes from dinoflagellate, dinoflagellate, transcriptome, small batch culture, cell

Abstract

This protocol provides instruction on how to generate transcriptomes from Dinoflagellates, either from single-cell or small batch cultures.

Before start

NOTE: PLEASE READ THROUGH THE ENTIRE PROTOCOL TO DETERMINE ALL OF THE NECESSARY ITEMS REQUIRED. THERE ARE MANY CONSUMABLE AND NON-CONSUMABLE ITEMS NEEDED, WITH SOME OPTIONS FOR REPLACEMENT.



Step 1. Isolate Cells in RNAlater

4h

- 1 Utilizing your own way to isolate single dinoflagellate cells, add a cell to  20 μL of RNAlater into a capped PCR tube (either strips or individual tubes).

Note: a PCR plate doesn't really work with this, as the films used allow the RNAlater to evaporate, even while frozen. Capped tubes or strips are best.

- 1.1 Allow the cells to sit in the RNAlater at room temperature for at least  04:00:00 , then freeze (either at  $-80\text{ }^{\circ}\text{C}$ or  $-20\text{ }^{\circ}\text{C}$. The RNAlater is also capable of preserving cells and RNA at  Room temperature for up to 2 weeks.

4h

Step 2. RNA Isolation

26m 30s

- 2 *Note: please also reference the protocol that comes with the extraction kit used in this step*

For downstream methods, RNAlater must be removed from the cells, otherwise it will inactivate sensitive enzymes. So, utilizing a **Macherey Nagel XS RNA extraction kit**, you will remove the RNAlater from the cell.

This technique works with both single cells and batch cells.

- 2.1 **Gather and prepare the necessary items/reagents:**

5m

a. One,  1.5 mL tube containing  1 mL of 70% ethanol (prepared fresh beforehand with 100% molecular grade ethanol and RNase-free H₂O;  700 μL of ethanol +  300 μL of water)

b. Remove the TCEP and thaw  On ice . Spin down:  7500 rpm, 4°C , 00:05:00
Remove the supernatant (~  16 μL) while leaving ~  4 μL in the tube

c. Pull out samples and place  On ice

d. Prepare a working solution of Carrier RNA (thaw on the bench): add  1 μL of stock solution to  99 μL of Buffer RA1. Do **NOT** place this working solution



On ice , as it will freeze due to the Buffer EA1. When done with the concentrated Carrier RNA, place this on ice until it can be returned to the freezer.

2.2 **Rupture the Cells:** Add 100 μL of Lysis Buffer RA1 and 2 μL of TCEP, vortex to mix (2 x, 00:00:05 each), then do a quick spin down.

2.3 **Add the Carrier RNA:** Add 5 μL of the Carrier RNA working solution (prepared beforehand) to each tube, vortex to mix (2x, 00:00:05 each), then do a quick spin down.

2.4 **Adjust the RNA binding conditions:** Add 100 μL of the 70% ethanol (prepared fresh beforehand) to the solution from step 2 (adding it to the 1.5 mL tube), vortex to mix (2x, 00:00:05 each), then do a quick spin down.

2.5 **Bind the RNA:** Obtain a nucleospin RNA XS column (light blue ring in individual packets - and should already be in a collection tube), mix the lysate (from Step 2.2) with a pipette and load the liquid into the top of the column (max: 650 μL), centrifuge 11000 x g, Room temperature, 00:00:30 , discard the used collection tube and place into a fresh tube.

30s

2.6 **Desalt the silica membrane:** Add 100 μL of the Membrane Desalting Buffer (MDB) to the top of the column, centrifuge 11000 x g, Room temperature, 00:00:30 , discard the flow-through and place the column back into the collection tube. *Note: if any additional liquid touches the bottom of the column, centrifuge again at* 11000 x g, 00:00:30 .

1m

2.7 **Digest the DNA:** Prepare the DNase mixture in a sterile, RNase-free (DEPC-treated as an option) 1.5 mL microcentrifuge tube. For each isolation, mix 10 μL of the reconstituted DNase to 90 μL of the reaction buffer (example for 2 RNA samples: 20 μL of DNase + 180 μL of buffer), mix by flicking the tube, quickly spin down, apply 95 μL of the solution directly to the center of the membrane (DO NOT PUNCTURE THE MEMBRANE), incubate at Room temperature for 00:15:00

15m



- 2.8 **Wash and dry the silica membrane:** add 100 μL of the Buffer RA2 to the top of the column, Room temperature incubate for 00:02:00 then centrifuge 11000 x g, Room temperature, 00:00:30 then discard the whole collection tube and flow-through and place the column in a new collection tube. Add 400 μL of the Buffer RA3 to the top of the spin column, pipetting around the edges to remove salts from the inner walls, then centrifuge 11000 x g, Room temperature, 00:00:30 then discard the flow-through and place back into the collection tube. Add 200 μL of the Buffer RA3 to the top of the spin column and centrifuge 11000 x g, Room temperature, 00:02:00 to dry the membrane. *Note: if any liquid touches the bottom of the column, then centrifuge again.* Then place the column into a 1.5 mL tube (provided by the kit). 3m
- 2.9 **Elute the RNA:** Add 5 μL of RNase-free water (provided by the kit). *Note: if you are extracting multiple cells for a total RNA isolation, utilize the maximum elution volume:* 10 μL . Let the tube incubate at room temperature for 00:01:00. Then centrifuge 11000 x g, Room temperature, 00:01:00 2m

Step 3. Check RNA quality and quantity (optional)

- 3 To check the RNA quantity and quality, we suggest using Agilent's High Sensitivity RNA ScreenTape Analysis for TapeStation Instruments (published 2025; this may become defunct later). Follow the instructions provided by the manufacturer.

Note: this will utilize 2 μL of your total 5-10 μL sample; if you believe that your procedure worked, we suggest that you move straight to Step 4, the SmartSeq v4 RNA kit procedure. However, for a first extraction attempt, it would be good to check to make sure your technique is reliable before moving on.

If you extracted from multiple cells for a total RNA isolation, it is suggested that you check the quantity and quality for downstream applications.

Store the rest of the sample at -80 $^{\circ}\text{C}$ until Step 4 (if Step 4 will not occur directly after the completion of Step 3).

Step 4. SmartSeq v4 RNA Kit

56m

4 *Note: please also reference the original protocol found on the manufacturer's website*

4.1 Step A: First Strand cDNA synthesis: Set up

1. Clean and Prep the PCR workstation: wipe down with 10% bleach, followed by 70% ethanol, and RNA-away (or Eliminate, whatever is the preferred Rnase and Dnase remover). Also, clean the vortexer and mini-centrifuge for inside the workstation.

2. Obtain 2 buckets of ice: one for reagents, one for the RNA samples (you want to keep these separate).

3. Preheat your thermal cyclers:

*One at 72 °C

*One with the following program:

i. 42 °C for 01:30:00

ii. 70 °C for 00:10:00

iii. 4 °C for ∞

4. Depending on the number of cells that you are using:

a. More than 1 cell (such as in a total RNA extraction), calculate the amount of RNA needed to make 10 ng total with a volume of 9.5 µL (RNase-free water can be used to make up this volume). If you have the results from the TapeStation (*Step 3: Check RNA quality and quantity*) post RNA isolation, you can use them here to do the calculations. Take the time to analyze these results and consolidate the information before beginning the procedure.

b. If using a single-cell isolation, it is suggested to proceed directly through the protocol so as not to lose sample volume.

5. Mix the 5x Ultra Low First Strand Buffer and SmartSeq V4 Oligonucleotide (hereafter termed Buffer+Oligo) as described in Table 1, and keep On ice .

	A	B	C	D
	Master Mix	1x (uL)	6x (uL) w/o ctrl	8x (uL) w/ ctrl
	5x Ultra Low First Strand	4	26.4	35.2

	A	B	C	D
	Buffer			
	SMARTseq v4 Oligonucleotide (48 uM)	1	6.6	8.8
	Total Volume	5	33	44

Table 1. SmartSeq V4 primer binding PCR.

4.2 Step A: First Strand cDNA synthesis: Begin Procedure

5m

1. Thaw the 5x Ultra Low First Strand Buffer at  Room temperature . Store at  Room temperature after the first thaw, and vortex well before using.

2. Thaw the 10x Lysis buffer, RNase Inhibitor, Diluted Control RNA (optional, only if this is being used as a control), nuclease-free water, sample RNA, 3' SmartSeq CDS Primer IIA, and SmartSeq V4 Oligonucleotide ON ICE. *Note: remember that the sample RNA goes into its own bucket.* Vortex to mix each reagent (briefly and gently).

3. Prepare a stock solution of the 10x Reaction Buffer (keep  On ice until 3' CDS Primer IIA is thawed). Mix briefly and spin down (to avoid bubbles).

 19 μL 10x Lysis Buffer

 1 μL RNase inhibitor

 20 μL total volume 10x reaction buffer

4. Prepare each sample ( 10.5 μL total volume) in individual  0.2 mL RNase-free PCR tubes. If using control RNA as a positive, dilute to 10 ng/uL with nuclease-free water and RNase inhibitor ( 1 μL of inhibitor +  49 μL of water to create a stock solution of the RNase inhibitor).

5. Place the samples  On ice and add  2 μL of 3' SMARTSeq CDS Primer IIA. Mix well by vortexing, followed by a quick spin.

 10.5 μL total RNA in reaction buffer

🧪 2 μL 3' SMARTSeq CDS Primer IIA

🧪 12.5 μL total volume

6. Incubate the tubes at 🌡️ 72 °C in a pre-heated, hot-lid thermal cycler for

🕒 00:03:00 . During this time, begin the next step (7).

7. Prepare enough Master Mix for all reactions + 10% of the total reaction volume as shown in Table 2. Mix the previously made Buffer+Oligo with RNase Inhibitor:

	A	B	C	D
	Master Mix	1x (μL)	6x (μL) w/o ctrl	8x (μL) w/ ctrl
	Buffer+Oligo	5	33	44
	RNase Inhibitor (40U/ μL)	0.5	3.36	4.4
	Total Volume	5.5	36.3	48.4

Table 2. Buffer+Oligo and RNase Inhibitor Master Mix.

8. Immediately place the tubes from Step 5 🌡️ On ice for 🕒 00:02:00 . Ensure that the thermal cycler lid of another block is preheated to 🌡️ 42 °C .

9. Add 🧪 2 μL per reaction (plus 10%) of the SmartSeqScribe Reverse Transcriptase to the Master Mix from Step 7:

6x reaction (w/o ctrl): 🧪 13.2 μL

8x reaction (w/ ctrl): 🧪 17.6 μL

Note: do NOT vortex the master mix after the enzyme is added; swirl gently with the pipette tip.

IMMEDIATELY PROCEED TO THE NEXT STEP (10).

10. Add  7.5 μL of the Master Mix to each reaction tube. Mix the contents of the tube by gently pipetting followed by a brief spin down.

11. Place the tubes into the second thermocycler block (where the lid was preheated to  42 °C) and run the following program (preprogrammed in step 4.1.3):

 42 °C for  01:30:00

 70 °C for  00:10:00

 4 °C for ∞

Note: the reaction from Step 11 can be stored at  4 °C overnight.

4.3 Step B: cDNA Amplification by LD PCR: Set up

Note: Steps 1-3 below can be done in the PCR workstation, but the remaining steps must be done on the bench space afterward to avoid contamination of the PCR station.

1. Program a thermocycler block to run the following PCR program:

A	B
Temp	Time
95°C	1 min
98°C	10 sec
65°C	30 sec
68°C	3 min
Go to 2	nX
72°C	10 min
4°C	∞

Table 3. cDNA amplification PCR

Note: Regarding cycle number (nX in Step 5), you will want to optimize your PCR depending on the total RNA or number of cells being used in the beginning. For large RNA extractions (such as hundreds of thousands of cells), the maximum RNA volume and concentration should be used with the minimum number of cycles.

	A	B	C
	Input (total RNA)	Input (# of cells)	Typical # of cycles
	10 ng	1,000	7-8
	1 ng	100	10-11
	100 pg	10	14-15
	10 pg	1	17-18

Table 4. Cycling guidelines based on the amount of starting material

4.4 Step B: cDNA Amplification by LD PCR: Begin Procedure

1. Thaw all reagents (except enzymes)  On ice then vortex gently:

2x SeqAmp PCR Buffer
 PCR Primer IIA (12 uM)
 Nuclease-free Water

2. Prepare the PCR Master Mix (+10% of the total volume of each reagent) as shown in Table 5.

	A	B	C	D
	Master Mix	1x (uL)	6x (uL) w/o ctrl	8x (uL) w/ ctrl
	2x SeqAmp PCR Buffer	25	165	220
	PCR Primer IIA (12 uM)	1	6.6	8.8
	SeqAmp DNA Polymerase	1	6.6	8.8
	Nuclease-Free Water	3	19.8	26.4
	Total Volume	30	198	264

Table 5. PCR Master Mix.



Note: vortex the Master Mix before adding the polymerase (enzyme). Add the enzyme right before and swirl gently with the pipette tip to mix, followed by a quick spin down.

3. Add  30 μL of PCR Master Mix to each tube, which contains  20 μL of the first-strand cDNA product from Step A (First Strand cDNA synthesis); for a total volume of  50 μL . Mix well and spin down quickly.

4. Place the tubes into a preheated thermocycler with a heated lid. Run the PCR protocol step in Table 3. *Note: for Step 5, nX cycles, consult Table 4.*

4.5 **Step C: Purification of Amplified cDNA using Agentcourt AMPure XP Kit: Set up**

30m

Note: all of the steps in this section should be done on a lab bench, not in a PCR hood. Amplified DNA should never be introduced to a PCR hood or clean station.

1. Pull out the AMPure XP beads and let sit at  Room temperature for  00:30:00. If you are proceeding directly from Step B, this can be done while the PCR runs. Depending on the number of samples that you are working with, you will need to pull out one or more tubes (generally, the large bottle of beads is aliquoted to  500 μL to prevent constant temperature changes of removing the larger bottle from cold storage).

2. Prepare fresh 80% ethanol; you will need  400 μL per sample, plus a little extra for pipetting error.

3. Thaw the Elution Buffer and maintain at  Room temperature.

4. Fit a 96 well plate (semi-skirted) to hold the PCR tubes in the magnetic stand (depending on what magnetic stand you have).

4.6 **Step C: Purification of Amplified cDNA using Agentcourt Ampure XP Kit: Begin Procedure**

21m

1. Add  1 μL of 10X lysis buffer to each PCR product from Step B (the ENTIRE product).

2. Vortex the AMPure XP beads until they are evenly mixed, then add  50 μL of beads to each PCR tube/sample (*note: the beads are viscous, so work slowly*).

3. Mix by vortexing OR pipetting the entire volume up and down 10 times to mix (ensure that the PCR product is SATURATED with the bead mixture).



4. Incubate at Room temperature for 00:08:00 to let the cDNA bind to the beads.
5. Briefly spin the samples down, then place onto the magnetic separation device for ~ 00:05:00 (or longer) until the supernatant becomes COMPLETELY clear (no beads are left in the supernatant). *To check the supernatant for beads, pipette the liquid out and hold the tip up to the light. If you see any specs or brownish/red coloring, the supernatant still contains beads.*
6. With the samples STILL ON THE STAND: pipette out the supernatant and discard.
7. Then, again with the samples STILL ON THE STAND, add 200 μ L of the fresh 80% ethanol to each sample WITHOUT disturbing the beads. Wait 00:00:30, then carefully pipet off the supernatant (which contains contaminants). The cDNA will stay bound to the beads. **DO NOT DISTURB THE BEADS.**
8. Repeat step 7 again.
9. Briefly spin the samples down and replace to the magnetic stand. Leave on the stand for 00:00:30 then remove the remaining ethanol. (*Note: if the beads are seen in the ethanol, carefully pipette liquid back into the tube and let it sit a little longer.*)
10. Then, with the samples still on the stand, leave the tubes open and allow them to sit at room temperature to allow the ethanol to evaporate for 00:02:00. Ensure the pellet is dry, but only just. The pellet will look matte (without a shine), but it shouldn't be cracked. Keep checking the samples; some will dry faster than others.

*If you **under** dry the pellet, there will still be ethanol in the sample, which will reduce your cDNA recovery; if you **over** dry the pellet, it will take longer than 00:02:00 to rehydrate and will reduce your cDNA recovery.*
11. After the beads dry, with the samples still on the stand, add 17 μ L of Elution Buffer to cover the pellet. Then remove the samples from the magnetic stand, and mix thoroughly to resuspend the beads (light vortex or swirl with the pipette tip).
12. Incubate at Room temperature (with the tubes closed) for 00:02:00.
13. Quickly spin down the tubes. Place the samples back on the magnetic stand for 00:01:00 or longer. Ensure the supernatant is completely clear. *If the beads don't*



pellet, pipette the supernatant and beads together, then pipette down towards the magnet. Continue until no beads remain in the supernatant. This can take a while for some samples (the beads can sometimes have difficulty pelleting).

14. Transfer the completely clear supernatant containing the purified cDNA from each tube into a nuclease-free, low-adhesion tube ( 1.5 mL). Ensure that the tubes are properly labeled. Aliquot out  2 μ L of the cDNA for bioanalysis (place this sample into a  0.2 mL , nuclease-free PCR tube). Samples can then be stored at  -20 °C indefinitely.

Step 5. Bioanalyzer

- 5 Using a bioanalyzer of your choice, determine the size and quantity of the cDNA produced. If the samples are of high enough quality following the previous procedures, proceed with Step 6: Nextera Library Preparation Kit.

Step 6. Nextera Library Preparation Kit

33m 30s

- 6 *Note: Please also reference the original protocol, as well as the pooling guidelines that can be found on Nextera's website.*

6.1 Step A: Tagment Genomic cDNA: Set up

1. Clean your mini-centrifuge and ensure that there is a microplate adapter on your vortexer (replacing the standard vortexing head).

2. Consumables for this step:

- a. ATM (Amplicon Tagment Mix): Thaw  On ice , gently invert the tube 3-5 times, then briefly centrifuge
- b. TD (Tagment DNA buffer): Thaw  On ice , invert 3-5 times, then briefly centrifuge
- c. NT (Neutralize Tagment buffer): Check for precipitates, vortex until resuspended
- d. cDNA from Step 4 (SmartSeq V4 kit): 0.2 ng/uL per sample max concentration; thaw  On ice , invert 3-5 times or vortex lightly; note: some samples will not make this maximum.
- e. Hardshell 96-well plate (skirted)
- f. Microseal B adhesive seals



3. Enter and save the following protocol ("Tagmentation Protocol") into a thermocycler (make sure it can hold a 96 well plate):

1-Preheat Lid

2- 55 °C for 00:05:00

3-Hold at 10 °C

4. Pull out the Index Primers you plan to use (these need to be thawed at

Room temperature for 00:20:00).

6.2 Step A: Tagment Genomic cDNA: Begin Procedure

6m

1. Prepare your cDNA: your target volume is 5 µL with 1 ng total or 0.2 ng/µL. *If your concentration is higher than 0.2 ng/µL, you will need to dilute with RNase-free water or Tris Buffer. If using Tris, it needs to be 10 mM with a pH of 7.5-8.5.*

Then add the following reagents to a hard-shell PCR plate in the order listed:

1- TD: 10 µL

2-cDNA: 5 µL

This ensures that the most valuable item is added last, as well as to reduce potential contamination of the TD reagent.

2. Add 2 µL of ATM to each well. Pipette to mix. Seal with a Microseal B adhesive seal.

3. Centrifuge at 280 x g, 20°C, 00:01:00 (in a plate centrifuge).

4. Place the plate into the pre-programmed thermocycler (Section 6.1.3). Run the Tagmentation Protocol. *Once the program reaches 10 °C , immediately proceed to the next step (step 5), as the transposome is still active.*

5. Quickly spin down the plate. At Room temperature , carefully peel up the seal then add 5 µL of NT to each well. Pipette to mix. Replace the seal (use a new one if necessary).

6. Incubate at Room temperature for 00:05:00 . Prepare for Step B (Amplify Libraries) during this incubation time. Then proceed with Step B.

6.3 Step B: Amplify Libraries: Set up

1. Consumables for this step:

a. NPM (Nextera PCR MasterMix): Thaw  On ice

b. Index Adapters/Primers: Thaw at  Room temperature for  00:20:00 (see above); spin down in a mini centrifuge; you will need i7 primers and i5 primers

c. TruSeq plate fixture

d. Microseal A film

e. Replacement caps for the Index primers

2. Save the PCR program found in Table 6 (below) into a thermocycler that can hold a 96 well plate and has a heated lid

	A	B
	Temp.	Time
	72°C	3 min
	95°C	30 sec
	95°C	10 sec
	55°C	30 sec
	72°C	30 sec
	Go to 3	12x
	72°C	5 min
	10°C	∞

Table 6. Index Adapter PCR protocol

6.4 Step B: Amplify Libraries: Begin Procedure

1. Arrange the Index Primers in the TruSeq Index plate fixture with your S5 (i5) primers in the A-H rows, and N7 (i7) primers in the 1-12 columns.

2. Using a pipette, add  5 µL of each index 1 (i7) adapter down each column.

Replace the cap on each i7 primer with a new **orange cap**. IT IS ALSO CRITICAL THAT YOU CHANGE YOUR PIPETTE TIP EVERY TIME YOU DISPENSE THE  5 µL. THIS WILL PREVENT CONTAMINATION ISSUES.



3. Using a pipette, add  5 μL of each index 2 (i5) adapter to each row. Replace the cap on each i5 primer with a new **white cap**. IT IS ALSO CRITICAL THAT YOU CHANGE YOUR PIPETTE TIP EVERY TIME YOU DISPENSE THE  5 μL . THIS WILL PREVENT CONTAMINATION ISSUES.
4. Add  15 μL of NPM to each well containing the index adapters. Pipette to mix. AGAIN, ALWAYS ENSURE TO CHANGE YOUR PIPETTE TIPS BETWEEN WELLS.
5. Place a new Microseal A film on the plate: remove the purple backing, place on the plate over the wells using gloved fingers and a Kim wipe; make sure the silicone is pressed down tight - this is done by pushing on each well containing a sample; leave the see-through instruction sheet ON. Centrifuge at  280 x g, 00:01:00 .

6. Place the plate into the pre-programmed thermocycler (Table 6). Ensure that the lid of the PCR machine is moderately tight.

Note: you can stop here and store the plate at  2-8 °C for up to 2 days; if you do so, replace the "A" film with the more resilient "B" film.

If you decide to continue in the same day, pull out tubes of AMPure XP beads and let sit at  Room temperature for  00:30:00 minimum.

6.5 Step C: Clean Up Libraries: Set up

Consumables:

- RSB (Resuspension buffer): Thaw at  Room temperature . Store at  4 °C after the initial thaw.
- AMPure beads (aliquoted): Thaw at  Room temperature for  00:30:00 minimum; pull these out during Step B while the thermocycler runs (last step) if continuing in the same day. Depending on the volume you are adding to each sample and the number of samples, you may need to pull out multiple tubes (each should be aliquoted at  500 μL); the AMPure beads should not be reused after leaving at room temperature.
- Prepare fresh 80% ethanol:  400 μL per sample is required, plus some extra for pipetting error.
- V-bottomed or round-bottomed midi plate (v-bottomed is preferred).

6.6 Step C: Clean Up Libraries: Begin Procedure

27m 30s

1. Centrifuge the 96 well plate, taken either directly from the thermocycler or from the refrigerator, at  280 x g, Room temperature, 00:01:00 .
2. Transfer  50 μL of the PCR product from each well to the corresponding well of a v-bottomed (or round-bottomed) midi plate.
3. Add  30 μL or  90 μL of the AMPure XP beads to each well: this is a size selection step. If you wish to grab smaller fragments, use  30 μL ; if you wish to grab larger fragments, use  90 μL . Using  30 μL is best when working with sheared cDNA (see the Nextera Manual). Cover with a Microseal "B" film.
4. Shake the plate at  1800 rpm, Room temperature , 00:02:00 (you can use a 96 well plate adapter, and the vortexer can be set to an approximate RPM; check your vortexer specs).
5. Incubate at room temperature for  00:05:00 . Centrifuge at  280 x g, 00:01:00 , and then briefly at a higher speed.
6. Place on the magnetic stand and wait until the liquid is clear (~  00:02:00).
7. Remove and discard all supernatant from each well. Ensure the liquid is clear by holding the pipette tip containing the liquid up to the light (but ensure no liquid is dripped out of the end or flows back into the pipetter).
8. Wash the samples 2 times:
 - a. Add  200 μL of fresh 80% ethanol to each well
 - b. Incubate on the magnetic stand for  00:00:30
 - c. Remove and discard all supernatant from each well (use the same technique as described above)
9. Using a small pipette tip ( 20 μL), remove any remaining ethanol.
10. Air dry the pellets on the magnetic stand for  00:15:00 (film removed).

Note: you are looking for the pellet to lose its shine (becoming matte), but you do not want the pellet to crack. Over-drying and under-drying can reduce your yield and cause issues downstream. So watch the pellet for the entire time frame.

11. Remove the plate from the magnetic stand.

12. Add  52 μL of the RSB to each well (change tips each time to avoid contamination).

13. Add a new Microseal "B" film to the plate top and secure. Shake at  1800 rpm, Room temperature , 00:02:00

14. Incubate at room temperature for  00:02:00 . Centrifuge at  280 x g, Room temperature, 00:01:00 with a very quick spin at a higher speed to get the excess liquid off of the film.

15. Place the plate back on the magnetic stand and wait until the liquid is clear.

16. Transfer the supernatant to a fresh PCR plate (with wells corresponding to each other to keep track of your samples and indices). Seal with Microseal "B" film. Check the supernatant for any transfer of beads. Aliquot out  2 μL of each sample into individual PCR tubes. Place these tubes and the plate at  -20 $^{\circ}\text{C}$ for storage.

17. Submit the PCR tubes containing the  2 μL aliquots for Bioanalysis and qPCR. This will allow for normalization of the samples manually in Section 7 (below) for an example spreadsheet for normalizing samples. Once the samples are normalized and pooled (see the pooling protocol from Nextera), samples can then be submitted for sequencing.

7. Example Normalization Spreadsheet

7

	A	B	C	D	E	F	G	H	I
	Sample	ng/ μL	avg bp size	dsDNA mw	nM	Final nM	Final volume	Initial volume	RSB to add
	Gpacif_sh	0.487	249	660	2.96336863	10	5	16.87268994	N/A

	A	B	C	D	E	F	G	H	I
	GspRII_sh	0.73 9	278	660	4.02768 694	10	5	12.41407 307	N/A
	Gsilvae_sh	0.32 9	255	660	1.95484 254	10	5	25.57750 76	N/A
	Gruetz_sh	0.811	272	660	4.51760 25	10	5	11.067817 51	N/A
	Gcarib_sh	0.751	343	660	3.31743 087	10	5	15.07190 413	N/A
	Gaust_sh	1.55 3	327	660	7.195811 32	10	5	6.948486 8	N/A
	Gpacif	5.82 9	467	660	18.91181 62	10	5	2.643849 717	2.35615 028
	GspRII	5.94 1	519	660	17.34395 98	10	5	2.882848 005	2.1171519 9
	Gsilvae	10.73	524	660	31.0259 079	10	5	1.6115563 84	3.38844 362
	Gruetz	10.9 3	532	660	31.12895 88	10	5	1.606221 409	3.39377 859
	Gcarib	6.96 4	497	660	21.23041 28	10	5	2.355112 005	2.64488 8
	Gaust	3.46 6	433	660	12.12821 05	10	5	4.1226197 35	0.87738 027
	qPCR results								

	A	B	C	D	E	F	G	H	I
	P1N	1.916	494	660	5.87657 956	10	30	51.05010 438	
	P2S	0.90 8	299	660	4.601195 91	10	30	65.20044 053	58.12527 25
		2.33	494	660	7.146362 41				
		1.3	299	660	6.58761 528				

Table 7. Example Normalization Spreadsheet