

Sep 07, 2018 Version 2

DNA extraction protocol for *Cryptosporidium* spp in stool samples (Adapted from GF-1 Soil Sample DNA extraction Kit) V.2

 Version 1 is forked from [DNA extraction protocol for *Cryptosporidium* spp in stool samples \(Adapted from GF-1 Soil Sample DNA extraction Kit\)](#).

 [PLOS One](#)

DOI

[dx.doi.org/10.17504/protocols.io.tb5eiq6](https://doi.org/10.17504/protocols.io.tb5eiq6)

Asar Khan¹, Sumaira Shams¹, Saima Khan¹, Muhammad Iftikhar Khan², Sardar Khan³, Abid Ali¹

¹Department of Zoology Abdul Wali Khan University Mardan, KPK Pakistan;

²Centre for Biotechnology and Microbiology University of Swat, Pakistan;

³Department of Environmental Science, University of Peshawar, KPK Pakistan



Asar Khan

Department of Zoology Abdul Wali Khan University Mardan 2320...

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.tb5eiq6](https://doi.org/10.17504/protocols.io.tb5eiq6)

Protocol Citation: Asar Khan, Sumaira Shams, Saima Khan, Muhammad Iftikhar Khan, Sardar Khan, Abid Ali 2018. DNA extraction protocol for *Cryptosporidium* spp in stool samples (Adapted from GF-1 Soil Sample DNA extraction Kit). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.tb5eiq6>

Manuscript citation:

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 07, 2018

Last Modified: September 07, 2018



Protocol Integer ID: 15453

Abstract

The GF-1 Soil Sample DNA Extraction Kit is designed for the rapid and efficient purification of bacteria DNA from soil samples without the need for precipitation or organic extraction. The kit uses a high pure specially-treated silica-based material fixed into a column to efficiently bind DNA in the presence of high salt. The kit applies the principle of a minicolumn spin technology and the use of optimized buffers to ensure that only DNA is isolated while cellular protein, humic acid, metabolites, salts and other low molecular weight impurities are removed during the subsequent washing steps. We used this protocol for *Cryptosporidium* spp DNA extraction from stool samples with minor modifications and its working.

Guidelines

All steps should be carried out at room temperature unless stated otherwise.

Materials

MATERIALS

☒ Centrifuges 5810 R **Eppendorf**

☒ Micropipetes tips

GF-1 Soil sample DNA extraction Kit **Vivantis Technologies Sdn. Bhd. Revongen Corporation Center**
☒ **No.12A, Jalan TP 5, Taman Perindustrian UEP, 47600 Subang Jaya, Selangor Darul Ehsan, Malaysia. Catalog #GF-SD-025**

☒ Water bath at 70 °C

☒ Micropipettes (2-20 µl, 50-250µl, 100-1000µl)

☒ Eppendorf tubes (1.5 & 2.0 ml)

☒ Vertex mixer

Safety warnings

! Buffer SB contains guanidine salts which can be harmful when in contact with skin or swallowed. Always wear gloves and practice standard safety precautions. Do NOT disinfect guanidine or extraction waste in solutions containing bleach or any other form of acid. To clean any items contaminated with the reagent, simply soak in detergent and water to remove all traces of guanidine before cleaning with bleach or acidic solution.


Before start


Pre-set the water bath at 70°C.

Lysis

- A.** 1g of stool sample was taken in a 2ml microcentrifuge tube and 500mg (0.5g) of Glass beads was added. Than 1ml of SLX-buffer was added into each tube and vertexd at maximum speed for 5 minutes. Now, 100µl of DS-Buffer was added & vertexd agian for 3 minutes and then incubated at 70°C for 10 minutes.

B. The tubes were then Centrifugated at 5,000xg for 3 minutes at room temperature and 800µl of the supernatant was then moved into a new microcentrifuge tube. A 270µl of P2-Buffer was added into the sample tube and mixed thoroughly by vertexing and incubated on ice for 5 minutes. Then the tubes were Centrifugated again at 14,000xg for 3 minutes to pellet the stool sample and all the supernatant was again transferred into a new 2ml microcentrifuge tube.


 1 g of stool sample

 500 mg of Glass beads (provided in the Kit)

 1 mL SLX-buffer


 100 µL DS-Buffer


 00:05:00 vertex time at max speed

 00:03:00 vertex time, again at max speed

Note

During incubation mix the sample by vertexing twice.


 800 µL supernatant moved to new microcentrifuge tube


 270 µL P2-Buffer was added and vertexed

 00:05:00 vertex time (for ice incubation)

DNA Precipitation

- A 0.7 volume (700 µl) of iso-propanol was added and mixed well in the tubes by inverting the tubes 20-30 times and centrifugated at 14,000xg for 10 minutes. The supernatant was discarded carefully without dislodging the DNA pellet. The tubes were inverted for 1 min on the paper towel to drain the remaining liquid.


 00:10:00 centrifugated at 14,000xg

 700 µL iso-propanol was added into tubes and mixed well

DNA Solubilization


- 3 200µl of EB (Elution Buffer) was added to the tubes and assorted by pulse-vertexing and incubated at 70°C for 10-20 min to dissolve DNA pellet. 100µl of HTR Reagent was added and mixed thoroughly by vertexing for 10 secs and then incubated at room temperature for 2 mins. The tubes were again centrifuged at 14,000g for 2 mins and the supernatant was transferred into a new clean 1.5ml microcentrifuge tube.

 200 µL Elution buffer added and mixed through pulse-vertexed

 70 °C incubation temp.

 00:20:00 for incubation

 100 µL HTR-reagent added and mixed

 00:00:10 vertexing time for HTR reagent


 20 °C RT for incubation

 00:02:00 for incubation

Protein digestion

- 4 2µl of Proteinase-K was poured to each tube (1.5 ml) and vortexed thoroughly and incubated at 37 °C for 10 mins. Then the XP1-Buffer equal to the sample volume (1.5 ml tube volume) was added and vortexed.

 37 °C Incubation

 00:10:00

 2 µL Proteinase-K was added and mixed

DNA wash and purification

- 5 **A.** The columns were inserted into each 2ml collection tube and 600µl of sample was transferred into each column. The tubes were centrifugated at 10,000 x g for 1 minute and the flow-through was discarded. Then 300µl of XP1 Buffer was added and centrifuged at 10,000 x g for 1 min and the flow thorough in the collection tubes was discarded.
- B.** The columns were inserted again into a new 2ml collection tube, and 700µl SPW Wash Buffer was added and centrifugated at 10,000 x g for 1 min, the flow thorough was discarded. The 700µl of SPW Wash Buffer step was repeated and the flow thorough was discarded.

C. The columns were centrifugated at 14,000g for 2 mins to remove all the traces of ethanol and then placed into a clean 2ml Eppendorf tube. 60µl of Elution Buffer was added directly onto the center of the membrane and incubated at 70°C for 5 mins.

🧪 600 µL sample was transferred to column

🧪 300 µL of XP-1 buffer was added and centrifugated at 10,000xg

⌚ 00:01:00 for centrifugation

🧪 700 µL SPW wash buffer was added into columns and centrifuged at 10,000xg

⌚ 00:01:00 Centrifuge time

Note

The 700µl of SPW Wash Buffer step was repeated and the flow thorough was discarded.

🧪 60 µL elution buffer was added onto the center on the cloumn membrane

🌡️ 70 °C Incubation of the columns

⌚ 00:05:00 incubation time

Elution of DNA and storage

- 6 The columns were centrifugated at 14,000g for 1 min to elute DNA, and the elution phase was performed twice with a second 60µl of Elution Buffer. The columns were removed and discarded while the extracted DNA (Eppendorf tube) was stored at -20 °C until amplification.

🧪 60 µL Elution buffer (repeated step)

🌡️ -20 °C for DNA storage