


Dec 11, 2023

EdU Staining Protocol in juvenile *F. hepatica*

 [PLOS Pathogens](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.eq2lyjnrrlx9/v1>



Paul McCusker¹, Duncan Wells¹, Rebecca Armstrong¹, Emily Robb¹, Paul McVeigh¹, Nathan Clarke¹, Erica Gardiner¹, Erin McCammick¹, Aaron Maule¹

¹Queen's University Belfast



Paul McCusker

QUB

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.eq2lyjnrrlx9/v1>

Protocol Citation: Paul McCusker, Duncan Wells, Rebecca Armstrong, Emily Robb, Paul McVeigh, Nathan Clarke, Erica Gardiner, Erin McCammick, Aaron Maule 2023. EdU Staining Protocol in juvenile *F. hepatica*. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.eq2lyjnrrlx9/v1>



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 31, 2023

Last Modified: December 11, 2023

Protocol Integer ID: 85723

Keywords: Fasciola, EdU, Proliferative cells, juvenile fasciola hepatica with the thymidine analogue, juvenile fasciola hepatica, staining proliferative cell, staining protocol, schistosoma mansoni, thymidine analogue, proliferative cell

Funders Acknowledgements:

LIVER FLUKE MOTOR FUNCTION AND PARASITE CONTROL: EXPLOITING A 'TARGET VALIDATION TOOLBOX' AS A DRUG SCREEN-INTERFACE FOR FLUKICIDE DISCOVERY

Grant ID: BB/K009583/1

Probing in vivo parasite biology in vitro

Grant ID: NC/N001486/1

Exploiting stem cell biology for liver fluke control

Grant ID: BB/T002727/1

Abstract

This protocol describes the process of labelling and whole mount staining proliferative cells in *in vitro* juvenile *Fasciola hepatica* with the thymidine analogue 5-ethynyl-2'-deoxyuridine (EdU). The protocol is based on the *Schistosoma mansoni* EdU staining protocol developed by the Collins and Newmark groups. A full copy of that protocol is available here - http://collinslab.org/PDF/Whole_mount_EdU_Smansoni.pdf.

Guidelines

Fixation

While flat fixing is necessary for worms grown for more than three weeks *in vitro* younger worms may be suitable for free fixing in 4%FA.

Materials

CS50 - 50%  Chicken Serum, New Zealand Origin **Thermo Fisher Catalog #16110082** and 50%

 RPMI 1640 Medium, no phenol red **Thermo Fisher Catalog #11835105** with x1

 Antibiotic Antimycotic Solution (100×), Stabilized **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A5955**

PBSTx - 1X PBS + 0.3%  Triton X-100 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T8787-50ML**

4%FA -  500 µL  Formaldehyde solution **Merck MilliporeSigma (Sigma-Aldrich) Catalog #F8775-25ML**

with  4 mL PBSTx

4%PFA - Add  8 g Paraformaldehyde to  100 mL millipore water, stir/heat **in a flow hood** for 1 h between  55 °C and  60 °C on hotplate (**NB. must NOT go above 60°C; if it does, start again**). Following 1 h heating turn off heat and add NaOH dropwise until solution goes clear before adding  100 mL 2X PBS. and mix. Store at  4 °C for 2 weeks or  -20 °C in aliquots for 6 months.

EdU Detection Solution

	A	B	C	D	E
Reagent	200 (µl)	1 ml (µl)	5ml (ml)	10ml (ml)	
PBS	158	789	3.945	7.89	
100mM Copper Protectant in Water (Store 4oC)	2	10	0.05	0.1	
10mM AlexaFluor (488 or 545)/6-FAM azide	0.2	1	0.005	0.01	
500mM Ascorbic Acid (Make fresh – 88 mg in 1 ml PBS)	40	200	1	2	

Copper protectant -  Copper(II) sulfate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #451657**

Alexafluor 488 -  Azide-fluor 488 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #760765**



6-FAM azide -  6-Fam Azide **Metabion**

DAPI solution - 1:1000 of  1 mg/mL stock in PBSTx

Protocol materials

 Triton X-100 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T8787-50ML**

 RPMI 1640 Medium, no phenol red **Thermo Fisher Catalog #11835105**

 Azide-fluor 488 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #760765**

 Copper(II) sulfate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #451657**

 Chicken Serum, New Zealand Origin **Thermo Fisher Catalog #16110082**

 6-Fam Azide **Metabion**

 Antibiotic Antimycotic Solution (100×), Stabilized **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A5955**

 Formaldehyde solution **Merck MilliporeSigma (Sigma-Aldrich) Catalog #F8775-25ML**

 VECTASHIELD Mounting Medium **Vector Laboratories Catalog #H-1000**

Safety warnings

EdU Mutagen

As a thymidine analogue EdU is known for germ cell mutagenicity and reproductive toxicity. Particular care should be taken when disposing of the chemical and around pregnant mothers.



In vitro addition of EdU to juvenile *F. hepatica*

1d

- 1 Juvenile liver fluke are cultured in 50% Chicken serum (CS50) as described in McCusker et al. (2016).
- 2  10 μL of  10 millimolar (mM) 5-Ethynyl-2'-deoxyuridine (EdU) is mixed with  190 μL pre-warmed CS50 and added to worms for between 6 and 24 h.

STEP CASE

Stain immediately following EdU incubation  15 steps

Fixation

- 3 Add  20 μL **4% FA/4% PFA** to a petri dish.
- 4 Pipette worms in into **4% FA/4% PFA** and immediately flatten with coverslip and weight (e.g. half full 15 mL falcon tube) for 10 min.
- 5 Add **4% FA/4% PFA** to float coverslip. Lift coverslip and move worms into 1.5 mL eppendorf with  1 mL **4% FA/4% PFA**. Place tube on a rotator for further 10 min (4%FA)/2 h (4%PFA) at  Room temperature **OR** O/N at  4 °C .

OPTIONAL Storage Step

- 6 Wash in  1 mL **PBSTx** for 10 min.
- 7 Wash in  1 mL 1:1 **PBSTx:Methanol** for 10 min.
- 8 Store in  1 mL methanol at  -20 °C until required.
- 9 Rehydrate in  1 mL 50:50 **PBSTx:Methanol** for 10 min.



Staining

- 10 x2 washes in  1 mL **PBSTx** for 10 min each.
- 11 Incubate in  1 mL fresh **Proteinase K Solution** for 15 min at  Room temperature .
- 12 Postfix in  1 mL **4%FA** for 10 min.
- 13 Wash in  1 mL **PBSTx** for 10 min.
- 14 Incubate in **EdU Detection Solution** for 30 min at  Room temperature . **NB: cover with tinfoil from this point on.**
- 15 x2 washes in  1 mL **PBSTx** for 5 min each
- 16 Incubate in  1 mL **DAPI Solution** for 20 min at  Room temperature or O/N at  4 °C .
- 17 Rinse in **PBSTx** and mount in  VECTASHIELD Mounting Medium **Vector Laboratories Catalog #H-1000** .

Protocol references

McCusker, P., McVeigh, P., Rathinasamy, V., Toet, H., McCammick, E., O'Connor, A., Marks, N. J., Mousley, A., Brennan, G. P., Halton, D. W., Spithill, T. W., & Maule, A. G. (2016). Stimulating Neoblast-Like Cell Proliferation in Juvenile *Fasciola hepatica* Supports Growth and Progression towards the Adult Phenotype In Vitro. *PLoS Neglected Tropical Diseases*, 10(9), [e0004994]. <https://doi.org/10.1371/journal.pntd.0004994>