



Oct 16, 2024

ELISA-based quantitative detection of peptides in plant tissue

 Journal of Experimental Botany

DOI

<https://dx.doi.org/10.17504/protocols.io.n92ld8k67v5b/v1>

Maurice Koenig¹, Zarah Sorger¹, Shania Keh², Gunther Doehlemann¹, Johana Misas Villamil¹

¹Institute for Plant Sciences, Cluster of Excellence on Plant Sciences (CEPLAS), University of Cologne, Cologne, Germany.;

²Institute for Plant Sciences, University of Cologne, Cologne, Germany.



Shania Pin Yin Keh

The Sainsbury Laboratory

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n92ld8k67v5b/v1>

Protocol Citation: Maurice Koenig, Zarah Sorger, Shania Keh, Gunther Doehlemann, Johana Misas Villamil 2024. ELISA-based quantitative detection of peptides in plant tissue. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.n92ld8k67v5b/v1>

Manuscript citation:

Quantitative detection of plant signaling peptides utilizing ELISA. Maurice König, Zarah Sorger, Gunther Doehlemann, Johana C. Misas Villamil. bioRxiv 2024.06.27.600388; doi: <https://doi.org/10.1101/2024.06.27.600388>



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it is working

Created: October 07, 2024

Last Modified: October 16, 2024

Protocol Integer ID: 109259

Keywords: ELISA, Phytocytokines, plant signaling peptides, Zip1, Peptide detection, peptides in plant tissue, quantitative detection of peptide, plant tissue phytocytokine, plant immunity, signaling peptide, custom peptide antibody, containing peptide, peptide, acid phytocytokine, antibody, final bioactive peptide form, plant tissue, based assay, released peptide, detection in complex biological sample, maize leaf, sensitive recognition of zip1

Funders Acknowledgements:

Deutsche Forschungsgemeinschaft

Grant ID: DO 1421/5-2

Deutsche Forschungsgemeinschaft SFB1403

Grant ID: 414786233

Cluster of Excellence on Plant Sciences EXC 2048/1

Grant ID: 390686111

Abstract

Phytocytokines play a critical role in the fine tuning of plant immunity and development. The final bioactive peptide form can vary from 1-5 kDa and is usually expressed in low amounts, making them difficult to detect with conventional methods. Here, we present an ELISA-based assay that allows rapid and cost-effective detection of naturally released peptides in plant tissues. This protocol is optimized for Zip1, a 17-amino-acid phytocytokine derived from *Zea mays* that elicits salicylic acid signaling in maize leaves. Using a custom peptide antibody, we designed an experimental pipeline to achieve a specific, selective and sensitive recognition of Zip1 allowing its detection in complex biological samples. Using this protocol, we are able to quantify in leaf material native Zip1-containing peptides with a detection limit in the nanogram range. Our methodology has a three-day turn-over and allows the use of raw extract from plant tissue. This method can be adapted for the detection and quantification of a variety of plant signaling peptides.



Materials

EQUIPMENTS

- Tecan infinite M200 PRO (with Tecan i-control software)
- Multichannel Pipette (20 - 200µL)

MATERIALS

- 96-well microtiter plate (CORNING COSTAR 9018; Thermofisher scientific, Ref: 44-2504-21)
- Adhesive film for Microplates (VWR, Article No. 391-0620)
- SuperSignalTM West Pico PLUS Chemiluminescence-Substrate (Thermofisher Scientific, Catalog Number: 34580)
- ZIP1 peptide (Sequence: QPWEGESELKLATQGASVRR, GenScript Biotech)
- Anti-ZIP1 (Sequence: EGESELKLATQGASVRR; Davids Biotechnologie, polyclonal antibody from rabbit)
- Anti-rabbit IgG, HRP-linked Antibody (Cell Signaling Technology, Category Number: 7074S)
- Paper towel

BUFFERS

1) Bicarbonate Buffer

3.03g Na₂CO₃

6.0g NaHCO₃

1000ml dH₂O

2) ZIP1 Stocks for ELISA Standard Calibration Curve

In bicarbonate buffer, dilute the synthesized ZIP1 peptide to the following concentrations:

SL1: 20µg/ml

SL2: 10µg/ml

SL3: 5µg/ml

SL4: 0.2µg/ml

SL5: 0.1µg/ml

Store the ZIP1 stocks at -20°C.

3) PBS Buffer, pH 7.4

1.44 g Na₂HPO₄

0.2g KCl

0.246g KH₂PO₄

8.0g NaCl

1000ml dH₂O

4) Blocking Buffer

5% Skim Milk Powder dissolved in PBS Buffer, pH 7.4

5) Covering Agent



70% Ethanol

5% Bromophenol Blue

25% MilliQ H₂O

Designing the 96-well Plate Layout

- 1 The following is an example of an ELISA 96-well microtiter plate layout. The first two columns of the 96-well microtiter plate is used for the ZIP1 standard calibration curve. Remaining spaces is allocated to the experimental group.



It is recommended for a "barrier" of empty wells to be put between each experimental group duplicates.

Duplicates are recommended for every sample group. Whereby if space allows, triplicate will yield a higher (thus better) R-squared value for the standard calibration curve.

	1	2	3	4	5
A	ZIP1 (500ng/μL)	ZIP1 (500ng/μL)		Experimental Sample 1	Experimental Sample 1
B	ZIP1 (250ng/μL)	ZIP1 (250ng/μL)			
C	ZIP1 (100ng/μL)	ZIP1 (100ng/μL)		Experimental Sample 2	Experimental Sample 2
D	ZIP1 (50ng/μL)	ZIP1 (50ng/μL)			
E	ZIP1 (10ng/μL)	ZIP1 (10ng/μL)			
F	ZIP1 (5ng/μL)	ZIP1 (5ng/μL)		Experimental Sample 3	Experimental Sample 3
G	ZIP1 (1ng/μL)	ZIP1 (1ng/μL)			
H	ZIP1 (0ng/μL)	ZIP1 (0ng/μL)		Experimental Sample 4	Experimental Sample 4
	6	7	8	9	10
A			Experimental Sample 5	Experimental Sample 5	

B					
C			Experimental Sample 6	Experimental Sample 6	
D					
E					
F			Experimental Sample 7	Experimental Sample 7	
G					
H			Experimental Sample 8	Experimental Sample 8	
	11	12			
A	Experimental Sample 9	Experimental Sample 9			
B					
C	Experimental Sample 10	Experimental Sample 10			



Coating

30m

2 You can find the recipe to make ZIP1 SL stocks in the materials section.

15m

To create the ZIP1 standard calibration curve, refer to the table below to the specific volumes of bicarbonate buffer and SL stock required for each concentration.

The proportions of bicarbonate buffer to SL stock in the table below is designed for one well.

	ZIP1 (ng/ μ L)	Bicarbonate Buffer (μ L)	SL Stock Volume Added (μ L)	Total Volume (μ L)
	500	47.5	2.5 (SL2)	50
	250	47.5	2.5 (SL3)	50
	100	49.0	1.0 (SL3)	50
	50	37.5	12.5 (SL4)	50
	10	47.5	2.5 (SL4)	50
	5	47.5	2.5 (SL5)	50
	1	49.5	0.5 (SL5)	50
	0	50.0	0.0	50

Table 1. Proportion of bicarbonate buffer to different SL stock concentration to make the ZIP1 standard calibration curve.

3 Load 50 μ L of each of your samples to two separate wells for duplicates.

10m

Note

The samples' protein concentration should be normalised either to cells/mL (protoplast), to fresh weight (leaf samples), or to equal protein abundance with Bradford Assay.

4 Seal the microtiter plate with an adhesive film. Press on the film to ensure a firm seal for each well across the plate.

15s



4.1 Hold the plate at one end, and tap gently on the other end for even mixing of the samples.

5s

5 Rest the microtiter plate on a flat surface  Overnight at  4 °C

16h



Washing

15m

6 The following steps should be done beside a sink.

II

7 Remove all of the liquid in the microtiter plate in one strong swinging motion.

10s

7.1 Stack some paper towel on a flat surface, flip the microtiter plate and tap it repeatedly on the paper towel to fully remove any remaining liquid.

20s

8 Using a multichannel pipette, add  200 µL of PBS buffer to all wells.

3m

8.1 Hold the plate at one end, and tap gently on the other end to ensure the PBS buffer covers the bottom of all wells

5s

8.2 Remove all of the liquid in the microtiter plate in one strong swinging motion.

10s

8.3 Flip the microtiter plate and tap it repeatedly on clean paper towel to fully remove any remaining liquid.

20s

8.4  [go to step #8](#)

10m

Repeat this for a total of three times.

Blocking

5m

9 Using a multichannel pipette, add  200 µL of blocking buffer to all wells.

3m



10 Seal the microtiter plate with an adhesive film. Press on the film to ensure a firm seal for each well across the plate.

15s

10.1 Hold the plate at one end, and tap gently on the other end for even mixing of the samples.

5s

11 Rest the microtiter plate on a flat surface for 02:00:00 at Room temperature

2h



Primary Antibody Incubation

30m

12 [go to step #6](#)

15m

Complete the washing step.

13 Add 100 μL of primary antibody (a-ZIP1, concentration of 1:7 $\mu\text{g}/\text{mL}$, dissolved in blocking buffer) to each wells with samples.

5m

14 Add 200 μL of blocking buffer to the rest of the wells.

5m

15 Seal the microtiter plate with an adhesive film. Press on the film to ensure a firm seal for each well across the plate.

15s

15.1 Hold the plate at one end, and tap gently on the other end for even mixing of the samples.

5s

16 Rest the microtiter plate on a flat surface Overnight at 4 $^{\circ}\text{C}$

16h



Secondary Antibody Incubation

30m

17 [go to step #6](#)

15m

Complete the washing step.

18 Add 100 μL of the secondary antibody (a-rabbit, concentration of 1:1000 , dissolved in blocking buffer) to each wells with samples.

5m



19 Add  200 μ L of blocking buffer to the rest of the wells. 5m

20 Seal the microtiter plate with an adhesive film. Press on the film to ensure a firm seal for each well across the plate. 15s

20.1 Hold the plate at one end, and tap gently on the other end for even mixing of the samples. 5s

21 Rest the microtiter plate on a flat surface for  02:00:00 at  Room temperature 2h



Detection 40m

22  go to step #6 15m
Complete the washing step.

23 Add  50 μ L of SuperSignalTM West Pico PLUS Chemiluminescence-Substrate to each well with samples. 5m

24 Put the microtiter plate into the Tecan infinite M200 PRO machine, ready for detection. 10s

24.1 Open the "Tecan i-control" software, associated with the Tecan infinite M200 PRO machine. 2m



We will use the software available "MycoAlert Assay" template from the "Luminescence" category for detection.

Under "Plate definition", select "Thermofisher Scientific Transparent 96-well"

Under "Part of Plate", select all the wells with samples within for measurement plus an empty well in addition.

24.2 Click "Start" to begin the detection. 5m
Save the excel file after detection has been completed.



**Note**

Make sure you save the excel file after measurement is completed, before you begin the second round of detection or else the data will be lost.

- 25 Take the microtiter plate out and add  100 μL of covering agent to rest of the well

5m

Note

Take extra caution when adding the covering agent to the rest of the plate as the high percentage of ethanol makes the covering agent a volatile liquid that is prone to dripping during pipetting.



- 26 Put the microtiter plate into the Tecan infinite M200 PRO machine, ready for detection.

10s

- 26.1 We will use the same settings as previously used.



For reference:

We will use the software available "MycoAlert Assay" template from the "Luminescence" category for detection.

Under "Plate definition", select "Thermofisher Scientific Transparent 96-well"

Under "Part of Plate", select all the wells with samples within for measurement plus an empty well in addition.

- 26.2 Click "Start" to begin the detection.
Save the excel file after detection has been completed.

5m

**Data Analysis**

15m

- 27 In the detection step, we took two separate measurements in case of the covering agent accidentally drips into any of the sample wells. If the measurements with covering agent is done properly, we will take that data set for further analysis as the covering agent helps eliminate background noise.
- 28 On excel, calculate the average RLU of every duplicates from the standard calibration curve and sample wells.



5m



28.1 Take the average and deduct the RLU value from the empty well.

2m

29 Plot a scatter plot of the ZIP1 concentration (ng/ μ L) against the RLU value for your ZIP1 standard calibration curve.

5m

29.1 From "Chart Elements", select "Trendline" then select "More options" to open "Trendline Options"

1m

29.2 Under "Trendline Options", select "Linear" for the trendline, select "Display equation on plot", and select "Display R squared value on plot"

1m

Note

A R squared value of 0.95 and above is recommended for good practice.

30 Using the displayed equation on plot, you can calculate the concentration of ZIP1 present in your sample wells.

1m

Protocol references

- Alhajj M, Zubair M, & Farhana A. 2023. Enzyme linked immunosorbent assay. StatPearls.
- Alvarez AM. 2004. Integrated approaches for detection of plant pathogenic bacteria and diagnosis of bacterial diseases. *Annual Review of Phytopathology*, 42(1), 339-366.
- Alves TO, D'Almeida CT, Scherf KA, Ferreira MS. 2019. Modern approaches in the identification and quantification of immunogenic peptides in cereals by LC-MS/MS. *Frontiers in Plant Science*, 10, 1470.
- Arispe N, Diaz JC, Flora M. 2008. Efficiency of histidine-associating compounds for blocking the Alzheimer's A β channel activity and cytotoxicity. *Biophysical Journal*, 95(10), 4879-4889.
- Aydin S. 2015. A short history, principles, and types of ELISA, and our laboratory experience with peptide/protein analyses using ELISA. *Peptides*, 72, 4-15.
- Boonham N, Kreuze J, Winter S, van der Vlugt R, Bergervoet J, Tomlinson J, Mumford R. 2014. Methods in virus diagnostics: from ELISA to next generation sequencing. *Virus research*, 186, 20-31.
- Depotter JRL, Misas Villamil JC, Doeblemann G. 2022. Maize immune signalling peptide ZIP1 evolved de novo from a retrotransposon. *BioRxiv*, 2022-05. doi: <https://doi.org/10.1101/2022.05.18.492421>
- Diaz-Amigo C, Popping B. 2013. Accuracy of ELISA detection methods for gluten and reference materials: a realistic assessment. *Journal of Agricultural and Food Chemistry*, 61(24), 5681-5688.
- Gorris MT, Sanz A, Peñalver J, López MM, Colomer M, Marco-Noales E. 2020. Detection and diagnosis of *Xylella fastidiosa* by specific monoclonal antibodies. *Agronomy*, 11(1), 48.
- Gust AA, Pruitt R, Nürnberger T. 2017. Sensing danger: key to activating plant immunity. *Trends in Plant Science*, 22(9), 779-791.
- Hanada K, Tamai M, Yamagishi M, Ohmura S, Sawada J, & Tanaka I. 1978. Isolation and characterization of E-64, a new thiol protease inhibitor. *Agricultural and Biological Chemistry*, 42(3), 523-528.
- Hou S, Liu D, He P. 2021. Phytocytokines function as immunological modulators of plant immunity. *Stress Biology*, 1, 1-14.
- Hou S, Wang X, Chen D, Yang X, Wang M, Turrà D., Di Pietro A, Zhang W. 2014. The secreted peptide PIP1 amplifies immunity through receptor-like kinase 7. *PLoS Pathogens*, 10(9), e1004331.
- Hou S, Liu Z, Shen H, Wu D. 2019. Damage-associated molecular pattern-triggered immunity in plants. *Frontiers in Plant Science*, 10, 646.

- Huffaker A, Pearce G, Ryan CA. 2006. An endogenous peptide signal in Arabidopsis activates components of the innate immune response. *Proceedings of the National Academy of Sciences*, 103(26), 10098-10103.
- Huffaker A, Dafoe NJ, Schmelz EA. 2011. ZmPep1, an ortholog of Arabidopsis elicitor peptide 1, regulates maize innate immunity and enhances disease resistance. *Plant Physiology*, 155(3), 1325-1338.
- Huo T, Wang CT, Yu TF, Wang DM, Li M, Zhao D, Xiu-Ting L, Jin-Dong F, Zhao-Shi X, Song XY. 2021. Overexpression of ZmWRKY65 transcription factor from maize confers stress resistances in transgenic Arabidopsis. *Scientific Reports*, 11(1), 4024.
- Katikireddy KR, O'Sullivan F. 2011. Immunohistochemical and immunofluorescence procedures for protein analysis. *Gene Expression Profiling: Methods and Protocols*, 155-167.
- Koenig M, Moser D, Leusner J, Depotter JRL, Doehlemann G, Misas Villamil JC. 2023. Maize Phytocytokines Modulate Pro-Survival Host Responses and Pathogen Resistance. *Molecular Plant-Microbe Interactions*, 36(9), 592-604.
- Lin F, Jiang L, Liu Y, Lv Y, Dai H, Zhao H. 2014. Genome-wide identification of housekeeping genes in maize. *Plant Molecular Biology*, 86, 543-554.
- Liu B, Teng D, Yang Y, Wang X, Wang J. 2012. Development of a competitive ELISA for the detection of soybean α subunit of β -conglycinin. *Process Biochemistry*, 47(2), 280-287.
- Martinelli F, Scalenghe R, Davino S, Panno S, Scuderi G, Ruisi P, Villa P, Stroppiana D, Boschetti M, Goulart LR, et al. 2015. Advanced methods of plant disease detection. A review. *Agronomy for sustainable development*, 35, 1-25.
- Monaghan EK. 2003. Enzyme linked immunosorbent assay (ELISA) for detection of sulfur-rich protein (SRP) in soybeans (*Glycine max* L.) and certain other edible plant seeds. PhD Thesis, Florida State University. Accessed August 2024 <https://repository.lib.fsu.edu/islandora/object/fsu:180557/datastream/PDF/view>
- Posthuma-Trumpie GA, Korf J, van Amerongen A. 2009. Lateral flow (immuno) assay: its strengths, weaknesses, opportunities and threats. A literature survey. *Analytical and Bioanalytical Chemistry*, 393, 569-582.
- Qian G, Wang L, Wu Y, Zhang Q, Sun Q, Liu Y, Liu F. 2009. A monoclonal antibody-based sensitive enzyme-linked immunosorbent assay (ELISA) for the analysis of the organophosphorous pesticides chlorpyrifos-methyl in real samples. *Food Chemistry*, 117(2), 364-370.
- Rampitsch C, Ames N, Storsley J, Marien L. 2003. Development of a monoclonal antibody-based enzyme-linked immunosorbent assay to quantify soluble β -glucans in oats and barley. *Journal of Agricultural and Food Chemistry*, 51(20), 5882-5887.

- Rhodes J, Yang H, Moussu S, Boutrot F, Santiago J, Zipfel C. 2021. Perception of a divergent family of phytocytokines by the Arabidopsis receptor kinase MIK2. *Nature Communications*, 12(1), 705.
- Rzemieniecki J, Stegmann M. 2022. Regulation of pattern-triggered immunity and growth by phytocytokines. *Current Opinion in Plant Biology*, 68, 102230.
- Sakamoto S, Putalun W, Vimolmangkang S, Phoolcharoen W, Shoyama Y, Tanaka H, Morimoto S. 2018. Enzyme-linked immunosorbent assay for the quantitative/qualitative analysis of plant secondary metabolites. *Journal of Natural Medicines*, 72, 32-42.
- Venbrux M, Crauwels S, Rediers H. 2023. Current and emerging trends in techniques for plant pathogen detection. *Frontiers in Plant Science*, 14, 1120968.
- Wang XiangHong WX, Liu Tao LT, Xu Na XN, Zhang Yan ZY, Wang Shuo WS. 2007. Data from: Enzyme-linked immunosorbent assay and colloidal gold immunoassay for ochratoxin A: investigation of analytical conditions and sample matrix on assay performance. CABI Databases DOI: 10.1007/s00216-007-1506-6
- Wu J, Meng C, Shang H, Rong S, Zhang C, Hong J, Zhou X. 2011. Monoclonal antibody-based triple antibody sandwich-enzyme-linked immunosorbent assay and immunocapture reverse transcription-polymerase chain reaction for *Odontoglossum* ringspot virus detection. *Journal of Virological Methods*, 171(1), 40-45.
- Yoshimatsu K, Sawada JI, Jaziri M, Shimomura K. 1995. Detection of cardenolides by ELISA in plant sciences. *Studies in Natural Products Chemistry*, 15, 361-380.
- Ziemann S, van der Linde K, Lahrmann U, Acar B, Kaschani F, Colby T, Kaiser M, Ding Y, Schmelz EA, Huffaker A, Holton N, et al. 2018. An apoplastic peptide activates salicylic acid signalling in maize. *Nature Plants*, 4(3), 172-180.