


Jun 11, 2018

Establishment of axenic sea-ice diatom cultures, modified from Jaeckisch et al. (2011)

DOI

dx.doi.org/10.17504/protocols.io.qvzdw76



Anders Torstensson¹, Jody Deming¹

¹University of Washington



Anders Torstensson

University of Washington

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.qvzdw76

Protocol Citation: Anders Torstensson, Jody Deming 2018. Establishment of axenic sea-ice diatom cultures, modified from Jaeckisch et al. (2011). [protocols.io https://dx.doi.org/10.17504/protocols.io.qvzdw76](https://dx.doi.org/10.17504/protocols.io.qvzdw76)

Manuscript citation:

Torstensson, A., Young, J. N., Carlson, L. T., Ingalls, A. E. and Deming, J. W. (2019), Use of exogenous glycine betaine and its precursor choline as osmoprotectants in Antarctic sea-ice diatoms1. J. Phycol. doi:10.1111/jpy.12839

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

This protocol has been successful for the three tested species, *Nitzschia lecointei*, *Navicula perminuta*, *Fragilariaopsis cylindrus* CCMP1102. Cultures have been stable and uncontaminated for > 18 months (DAPI screening and no growth on BD Difco 2216 agar plates).

Created: June 11, 2018

Last Modified: June 11, 2018

Protocol Integer ID: 12953

Abstract

Table 1. Antibiotic concentrations used for establishing axenic sea-ice diatom cultures.

Antibiotic	Final concentration ($\mu\text{g/ml}$)
Ampicillin	50
Gentamycin	15
Streptomycin	125
Chloramphenicol	10
Ciprofloxacin	10
Penicillin	100

References

Guillard RRL. 1975. Culture of phytoplankton for feeding marine invertebrates. In: Smith WL, Chanley MH eds. *Culture of marine invertebrate animals*. New York: Plenum, 29-60.

Jaeckisch N, Yang I, Wohlrab S, Glöckner G, Kroymann J, Vogel H, Cembella A, John U. 2011. Comparative Genomic and Transcriptomic Characterization of the Toxigenic Marine Dinoflagellate *Alexandrium ostenfeldii*. *PLoS ONE* **6**(12): e28012.

Attachments



[Axenic diatom protoc...](#)

105KB

- 1 Pre-treat an optically dense diatom culture with 3 min of vortex-mixing (< 8°C) and 30 sec of ultrasonication (80% in 5-s pulses) on ice to remove attached bacteria from diatoms and EPS aggregates.
- 2 Centrifuge cells at 1000 × g for 10 min (< 8°C) and wash the pellet twice with sterile F/2 medium (Guillard, 1975) to remove the majority of loose bacterial cells.
- 3 Add six antibiotics to the medium according to Table 1 (see abstract), and incubate for 5 days at -1°C and under light (20–45 µmol photons m⁻² s⁻¹).
- 4 Centrifuge cells at 1000 × g for 10 min (< 8°C) wash pellet in fresh F/2 medium.

- 5 Pick < 10 cells to start the axenic cultures (pipetting or plate on F/2 agar to pick single colonies).
- 6 After recovery and growth (~2 months), check for bacterial cells using DAPI staining and after agar plating on half-strength BD Difco™ Marine Agar 2216 (< 8°C), kept at seawater salinity using artificial seawater (411 mM NaCl, 9.39 mM KCl, 26.1 mM MgCl₂·6H₂O, 28.4 mM MgSO₄·7H₂O, 5.01 mM TAPSO buffer).
- 7 If bacteria are still present, repeat protocol. Axenic cultures are usually established within the first or second attempt.