

May 31, 2025

Version 2

Experimental procedures for TF Perturb-Seq (UT Southwestern) V.2

DOI

<https://dx.doi.org/10.17504/protocols.io.14egn635pl5d/v2>

Lei Wang¹, Sushama Sivakumar¹, Yihan Wang¹, Sean C Goetsch¹, Vrushali Pandit¹, Huan Zhao¹,
Anjana Sundarajan¹, Daniel Armendariz², Chikara Takeuchi¹, Gary C Hon¹, Nikhil V Munshi¹, Mpathi Nzima¹

¹University of Texas Southwestern Medical Center; ²Yale University

IGVF @ UTSW



Mpathi Nzima

University of Texas Southwestern Medical Center at Dallas

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egn635pl5d/v2>

Protocol Citation: Lei Wang, Sushama Sivakumar, Yihan Wang, Sean C Goetsch, Vrushali Pandit, Huan Zhao, Anjana Sundarajan, Daniel Armendariz, Chikara Takeuchi, Gary C Hon, Nikhil V Munshi, Mpathi Nzima 2025. Experimental procedures for TF Perturb-Seq (UT Southwestern). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.14egn635pl5d/v2> Version created by **Mpathi Nzima**



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 01, 2025

Last Modified: May 31, 2025

Protocol Integer ID: 119423

Keywords: seq data production, experimental procedures for tf perturb, best practices for perturb, seq, tf perturb, ut southwestern, perturb, experimental procedure, optimized procedure

Abstract

Optimized procedures and best practices for Perturb-Seq data production, developed at the UT Southwestern Characterization Center

Part 1: sgRNA plasmid library construction

1

1.1 Synthesize sgRNA oligo pool using **Twist Bioscience**

sgRNA sequences were derived from this study: Replogle JM, Bonnar JL, Pogson AN, Liem CR, Maier NK, Ding Y, Russell BJ, Wang X, Leng K, Guna A, Norman TM, Pak RA, Ramos DM, Ward ME, Gilbert LA, Kampmann M, Weissman JS, Jost M. Maximizing CRISPRi efficacy and accessibility with dual-sgRNA libraries and optimal effectors. Elife. 2022 Dec 28;11:e81856. doi: 10.7554/eLife.81856.

Lenti-Guide(10X) backbone:

GTGGAAAGGACGAAACACCGNNNNNNNNNNNNNNNNNNNNgtttaagagctaagctggaaacagc

1.2 Resuspend sgRNA oligo pool with low EDTA TE buffer. [0.5 ng/ul. ~ 0.024uM]

1.3 Perform PCR to amplify double-stranded sgRNAs using **library-fwd** and **library-rev** primers to yield 139-bp fragments.

Primers (synthesized from IDT):

library-fwd: taacttgaaagtatttcgatttcttgctttatatcttGTGGAAAGGACGAAACACCG

library-rev: gttgataacggactagccttatttaaacttgctatgctgttccagcttagctcttaaac

Set up 4 tubes of 50 ul PCR reactions shown below (total 3.2 ng sgRNA oligos):

2.5 ul 10 uM **library-fwd** + **library-rev** primers

1.6 ul sgRNA oligo pool (0.5 ng/ul. ~ 0.024uM)

20.9 ul H₂O

25 ul 2X KAPA HotStart Readymix (Roche, Cat. # 7958935001)

PCR program:

95C for 30s

5 cycles(98C for 20s, 59C for 15s, 72C for 10s)

5 cycles(98C for 20s, 65C for 15s, 72C for 10s)

72 for 60s

4 forever

Combine reactions and run 1.2% Agarose gel to purify 139 bp products with Qiagen MinElute Gel Extraction Kit. Purify with 2 MinElute Spin Columns and elute with 12 ul EB buffer for each column (Total 24 ul). Measure DNA concentration with Qubit kit.

- 1.4 Digest 15 ug LentiGuide(10X)-BFP-Puro (LW203) in 200 ul reaction volume with 15 ul NEB Esp3I for 3h and run 0.7% agarose gel to purify the vector.

Plasmid map of LentiGuide(10X)-BFP-Puro (LW203): <https://benchling.com/s/seq-cp12NBXnoIGQY6NSXY5b?m=slm-GTRPOQU0bCAenT5kHEgB>

Purify with three QIAquick Spin Columns (for DNA size between 70 bps to 10 kbs) and elute with 35 ul EB buffer for each column. Combine three elutions together(105 ul total). Measure DNA concentration with Qubit kit.

- 1.5 Assemble the double-stranded sgRNA library into the lentiviral backbone with the NEBuilder HiFi DNA Assembly(NEB, Cat. # E2621X).

Set up 2 tubes of 80 ul reactions shown below:

500 ng Esp3I digested LentiGuide(10X)-BFP-Puro
20 ng sgRNA PCR library
Add water to 50 ul
Add 50 ul NEB HIFI DNA assembly Master Mix

Incubate at 50 degrees for 60 min.

Purify assembled DNA with 0.8X SPRIselect beads(Beckman Coulter, Cat. # B23318). Pipette mix 160 ul (0.8X) SPRIselect beads with 200 ul PCR reactions. Incubate at room temperature for 5 min. Keep the tube on the magnet until the supernatant is clear. Wash SPRIselect beads three times with 500 ul 80% ethanol. Elute in 20 ul EB buffer. Measure DNA concentration with Qubit kit.

- 1.6 Electroporation with Lucigen Endura Duo competent cells.

Electroporate 25 ul E.coli + 2 ul ligated DNA (>50 ng) with 0.1 cm cuvette with BioRAD MicroPulser Electroporator. Add 2 ml SOC recovery media, transfer the E.coli to a 14 ml tube. (Electroporate 100 ng plasmid library instead of ligated DNA for pre-made libraries).

Electroporate another 25 ul E.coli + 2 ul ligated DNA (>50 ng) with 0.1 cm cuvette with BioRAD MicroPulser Electroporator. Add 2 ml SOC recovery media, transfer the E.coli to a 14 ml tube.

Culture at 37 degree shaker for 1 h.

Spread 0.05, 0.01, 0.001 ul Ecoli containing media together with 50 ul LB to Agarose plates and culture at 37 degree overnight. The next day, count the plated E.coli colonies



to ensure that there are >1000 clones per sgRNA in 2 ml of culture.

Transfer 4 ml of E.coli culture (stored at 4 degree O/N) to a flask with 300 ml fresh LB + 150 ug/ml ampicillin and culture at 37 degree shaker for 16 hours. Pellet the E.coli. Extract plasmids with Zymo Midiprep kit. Elute plasmids with 250 ul + 200 ul elution buffer. Measure DNA concentration.

Part 2: Bulk sequencing of sgRNAs

2 Goal: This section describes procedures to construct next-generation sequencing libraries of sgRNAs, either from the plasmid library described in Part 1 or from genomic DNA where sgRNAs have been integrated into cells. Typically, we sequence sgRNAs from 3 samples:

- a. the plasmid library from Part 1;
- b. genomic DNA from iPSCs/hESCs before differentiation;
- c. genomic DNA from differentiating cells, ~5 days before harvest

2.1 For genomic DNA samples, prepare genomic DNA from >2M cells (or 1000X sgRNA number) using the Qiagen DNeasy Blood & Tissue Kit (Qiagen, Cat. # 69504).

2.2 Prepare PCR amplicon for NGS bulk sequencing (v2 with 10X TT index):

Primers needed:

Fwd-v2: ACACTCTTTCCCTACACGACGCTCTTCCGATCTgaaagtatttcgatttcttggt

Rev : GTGACTGGAGTTCAGACGTGTGCTCTTCCGATCTaagtgataacggactagcc

10X TT set A index or IDT-synthesized TT set A index (A1 ~ A10)

Map of sgRNA PCR amplicon v2 Benchling file: <https://benchling.com/s/seq-ev18tEY07pX59CNZ075A?m=slm-2MGRsOieiH7dX2Ms3Ecn>

Prepare 1st round of PCR

Prepare master mix of 4 ~ 12 tubes of PCR reactions(100 ul per reaction) (Calculate PCR reaction volume according to the size of sgRNA library. Set up enough amount of genomic DNA for PCR to make sure the complexity of the sgRNA library.)

5 ul 10 uM Fwd-v2 primer

5 ul 10 uM Rev primer

0.125 ug sgRNA plasmid library (or 0. 5 ug genomic DNA),add water to 40 ul.

50 ul 2X KAPA HIFI HS RM (Roche, Cat No.: 7958935001)

PCR protocol:

95°C for 2 min,

21 cycles for genomic DNA or 8 cycles for plasmids (27 cycles for genomic DNA or 15 cycles for plasmids if PCR products from 1st round were sequenced with Nanopore):

98°C for 20 s, 55°C for 20s ,72°C for 20s

72°C for 1 min

4°C forever

Mix all 100 ul PCR reactions together. Take 100 ul and add 60 ul (0.6X) SPRIselect to remove plasmids. After 5 minute incubation, keep the tube on the magnet until the supernatant is clear. Transfer 150 ul supernatant to a new tube and add 60 ul SPRIselect. Incubate at room temperature for 5 min. Keep the tube on the magnet until the supernatant is clear. Wash the SPRIselect beads three times with 300 ul 80% ethanol. Elute PCR product in 40 ul EB buffer.

Prepare 2nd round PCR

Set up 2nd round master mix of 1 PCR reaction.

20 ul 10X plate TT set A primer (or **5ul** 10 uM IDT-synthesized TT set A 1 ~ 10 primers + **15 ul** H2O)

30 ul PCR1 elution

50 ul 2X KAPA HotStart ReadyMix

PCR protocol:

98°C for 45s

8 cycles: 98°C for 20s, 62°C for 10s, 72°C for 10s

72°C 1 min

4°C forever

Post PCR SPRIselect purification

Mix 70 ul(0.7X) SPRIselect with 100 PCR products to remove large fragments.

After 5 minute incubation, keep the tube on the magnet until the supernatant is clear.

Transfer 150 ul supernatant to a new tube and add 30 ul SPRIselect. Incubate at room temperature for 5 min. Keep the tube on the magnet until the supernatant is clear. Wash the SPRIselect beads three times with 300 ul 80% ethanol. Elute with 35 ul EB for 275 bps PCR products. Spec DNA using a Qubit fluorometer and an Agilent TapeStation.

- Expected Qubit fluorometer concentration: 1 to 30 ng/ul.
- Expected the library size on a TapeStation: 275 bps



2.3 Submit the libraries for NGS illumina sequencing or do Nanopore sequencing (Oxford Nanopore Technologies) shown below

Nanopore sequencing according to the modified nanopore protocol (Ligation sequencing amplicons - Native Barcoding Kit 24 V14 (SQK-NBD114.24)). Briefly, repair and prepare the input amplicon DNA, ligate native barcodes and sequencing adaptors, load the libraries into either MinION or PromethION flowcell. The three modifications are: 1. Changing "Add 15 ul of resuspended AMPure XP Beads" to "Add **19** ul of resuspended AMPure XP Beads " in End-prep part 11, 2. Changing "0.4X AMPure XP beads (AXP)" to "**1.0X** AMPure XP beads(AXP)" in Native barcode ligation part 11, 3. Changing "Add 20 ul of resuspended AMPure XP beads (AXP)" to "Add **50** ul of resuspended AMPure XP beads (AXP)" in Adapter ligation and clean-up part 9.

Part 3: Transfection of Adherent Cells for Virus Packaging

3 Goal: Package sgRNA lentivirus.

3.1 The day before transfection, seed 293T (WT) cells at 5.5×10^6 cells per 10cm dish.

3.2 Check cell health. Do not proceed if cell density or health is not correct for transfection.

3.3 Mix the DNA in a sterile tube. Scale according to the number of plates to be transfected.

	A	B
	DNA	Amount for 10-cm dish
	transgene	8 ug
	psPAX2	6ug
	pMD2G	2ug

Transfection

3.4 Add 1ml Opti-MEM per plate to the tube, mix/vortex.



- 3.5 Add 64ul Transporter 5 transfection reagent per plate to the tube, vortex 5 seconds, gently spin down.
- 3.6 Incubate at room temperature for 20 minutes.
- 3.7 Gently pipette up and down 3 times to mix, then gently add 1ml per plate to cells, spreading around the plate.
- 3.8 Gently mix the dish.
- 3.9 Change the media the next day.

Part 4: Virus collection and concentration

- 4 **Goal:** Generate, harvest, and concentrate sgRNA lentivirus.
- 4.1 **Collect virus 48 hours after the media change:** For ultracentrifugation (UC), pre-chill Beckman Optima XL-100K centrifuge to 4C.
- 4.2 Draw media with virus up into 10 ml or larger syringe.
- 4.3 Attach a 0.45um filter, filter the virus into sterile tube(s) of appropriate size and type.
- 4.4 If concentrating:
 1. Filter in ultracentrifuge tubes at 25-30ml per tube. Supplement with DMEM if needed.
 2. Balance the ultracentrifuge buckets and contents.
 3. Concentrate virus with the SW28 Rotor at 25000 rpm @ 4C for 90 minutes.
 4. Pour out supernatant
 5. Resuspend in appropriate volume of complete mTeSR+CloneR2.
- 4.5 Aliquot virus if needed.



4.6 Bleach everything that has contacted virus using 50% bleach.

4.7 If titering, use the virus fresh for titer infection (Part 6).

4.8 Store remaining virus at 4C for large-scale infection (Part 7).

Part 5: Single-cell suspension of ESC/iPSCs

5 **Goal:** Detach cultured ESCs and iPSCs into a single-cell suspension for lentiviral infection.

5.1 Remove media from cells, wash with 2ml PBS, remove PBS.

5.2 Add 1ml TrypLE Select and incubate at 37C for up to 15min. Check every 3min and proceed to the next step when cells are rounded and able to detach with a firm tap to the plate.

5.3 Pipette the cells to dissociate them, then add 4ml PBS to neutralize.

5.4 Spin at 300 g for 5min. Carefully remove supernatant.

5.5 Resuspend pellet in media+Rock Inhibitor. Pipette thoroughly to break up remaining cell clusters and to fully separate cells from Matrigel.

5.6 Spin at 300 g for 5min. Carefully remove supernatant. When aspirating, watch for the "pop" from the top of the pellet that indicates the Matrigel layer of the pellet being removed.

5.7 Resuspend pellet in an appropriate volume of media depending on usage.

Part 6: Measure lentiviral infection in hESCs/iPSCs

6 **Goal:** Measure the titer of a lentivirus infected into H9/WTC11 cells.



- 6.1 Culture target cells to 80% confluency (~2M cells per well). Check cell health and density.
Do not proceed if cells look abnormal or overgrown.
- 6.2 4-5 hours before infection, replace culture media with media containing CloneR2.
- 6.3 Prepare the cells per the Single Cell Suspension protocol (Part 5). Resuspend in ~1ml media+Rock1 per 2 wells collected.
- 6.4 Count and aliquot appropriate number of cells - 1M cells per condition. Spin at 300*g for 5 minutes at room temperature. Remove supernatant completely.
Check cell viability. Repeat wash if viability is <80%. Do not proceed if still <80% live.
- 6.5 Resuspend cell pellet in media+CloneR2 at 1M cells per 1ml. Transfer to clean wells of a 6well Ultra-Low Attachment Tissue Culture plate.
- 6.6 Add virus - 10ul, 50ul, 200ul is usually a good starting point. Gently tap/shake the plate to mix.
If using concentrated virus, spike in at 10% with a different fluorophore to help estimate MOI.
- 6.7 Incubate cells in humidified 37C/5% CO2 for 3 hours, either on a plate shaker or tapping/shaking the plate every ~30 minutes to keep the cells in suspension.
- 6.8 After 3 hours, plate the cells to new 6well Matrigel coated plates with an additional 1ml media+CloneR2. Incubate the cells in a humidified 37C/5% CO2 incubator.
- 6.9 Change the media the next day.
- 6.10 Measure titer 2 days later:
 - a. Collect cells per the Single Cell Suspension protocol. Resuspend cell pellet at ~10M cells per ml (~100-200ul per well).
 - b. Measure BFP+ cells using a flow cytometer.
 - c. Calculate the appropriate amount of virus for full-scale infection, taking into account a ~10% loss in infectivity due to storing the virus in 4C for 2 days.

Part 7: Large-scale lentiviral infection of hESCs/iPSCs

- 7 **Goal:** Infect a sgRNA lentivirus into H9/WTC11 cells.

Notes: The below protocol is essentially the same as the titer infection, but with more cells and a single virus quantity. This infection is 2 days after the titer infection, using the same batch of virus.

- 7.1 Culture target cells to 80% confluency (~2M cells per well).
Check cell health and density. Do not proceed if cells look abnormal or overgrown.
- 7.2 4-5 hours before infection, replace culture media with media containing CloneR2.
- 7.3 Prepare the cells per the Single Cell Suspension protocol (Part 5). Resuspend in ~1ml per 2 wells collected.
- 7.4 Count and aliquot appropriate number of cells. Spin at 300*g for 5 minutes at room temperature. Remove supernatant completely.
Check cell viability. Repeat wash if viability is <80%. Do not proceed if still <80% live.
- 7.5 Resuspend cell pellet in 1ml media+CloneR2 per 1 M cells.
- 7.6 Add the appropriate amount of virus. Invert the tube several times to mix.
- 7.7 Transfer to a clean 6well Ultra-Low Attachment Tissue Culture plate at 1M cells per well.
- 7.8 Incubate cells in humidified 37C/5% CO2 for 3 hours, either on a plate shaker or tapping/shaking the plate every ~30 minutes to keep the cells in suspension.
- 7.9 After 3 hours, combine all cells and add media+CloneR2 to 500k cells per ml.
- 7.10 Plate 2ml (1M cells) per well to new 6-well Matrigel coated plates. Incubate the cells in a humidified 37C/5% CO2 incubator.
- 7.11 Change the media the next day.
- 7.12 Incubate the cells in a humidified 37C/5% CO2 incubator until sorted. Media change every 1-2 days. Split as needed, maintaining complexity of the cell library.
 - a. Terminate experiment if:

- i. Viability of cells is <50% or is too low to maintain library complexity
- ii. Cells appear unhealthy
- iii. Expression of marker fluorophore is much lower than expected

b. Take a representative picture before the first split to estimate integration efficiency

Part 8: Sort for BFP+ hESCs/iPSCs

8 **Goal:** Enrich for lentivirally infected cells by sorting for BFP+ cells

8.1 **Cell Sorting:** Prepare the cells per the Single Cell Suspension protocol (Part 5). Resuspend cells in media+Rock Inhibitor to ~10M cells/ml.

8.2 Meanwhile, prepare 15ml conical tubes to sort into by adding 2ml media+Rock Inhibitor and thoroughly coating the inside of the tube.

8.3 Strain the cells through a 40um strainer into the tubes to sort from.

8.4 Using a FACS instrument, sort for the desired population of BFP+ cells. Collect into the prepared 15ml tubes.

8.5 Spin at 300 g for 5min. Carefully remove supernatant.

8.6 Resuspend pellets in the appropriate amount of media+CloneR2.

8.7 Plate cells in 6-well plate - 250k to 1M cells per well.

8.8 Change the media 2 days after the sort.

8.9 Incubate the cells in a humidified 37C/5% CO2 incubator until differentiation. Media change every 1-2 days. Split as needed, maintaining complexity of cell library.

8.10 Two weeks after infection, differentiate the cells to cardiomyocytes using the appropriate protocol.

Part 9: Differentiation into cardiomyocytes

- 9 Goal: Use GiWi protocol to efficiently differentiate ES/iPSCs into cardiomyocytes
- 10 **CDM3 supplement:** Slowly add 1.056 g of L-ascorbic acid 2-phosphate to a 50 ml WFI water, inverting intermittently. Mix until clear. Add 2.5 g of recombinant human albumin, invert intermittently, mix until dissolved. Filter sterilize. Make 5 ml aliquots in 15 ml Falcon tubes and store at -20°C.
- 11 **CDM3 complete media:** RPMI 1640 0.5 L + CDM3 supplement 5 ml. CDM3 is stable at 4 °C for >4 weeks. Optional: supplementation with 4ug/ml insulin and 1 mg/ml polyvinyl alcohol (PVA) from d0-d8 greatly improves cell viability, cardiomyocyte yield and consistency. Before adding CDM3 supplement, add 0.5 g PVA to 500 ml RPMI 1640 and allow to solubilize at room temperature for a few hours, sterile filter and add 10 mM HEPES, CDM3 supplement and 4 ug/ml insulin.
- 12 **Start differentiation: Day 0:** CDM3-C: prepare sufficient CDM3-C by adding CHIR99021 at a final concentration of 3-4 µM to CDM3 complete medium and warm at 37°C. Aspirate medium and add 5 ml of CDM3-C.
- 13 NOTE: both plating density and CHIR99021 concentration should be optimized for maximal differentiation efficiency and consistency.
- 14 NOTE: Without the addition of insulin, significant cell death will be noted; hPSCs that undergo epithelial to mesenchymal transition will result in higher cardiomyocyte yield than those that maintain an epithelial morphology.
- 15 **Day 2:** CDM3-C59: prepare sufficient CDM3-C59 media by adding Wnt-C59 at a final concentration of 2 µM to CDM3 complete medium and warm at 37°C. At 48h, aspirate medium and add 5 ml of CDM3-C59.
- 16 NOTE: for optimal differentiation, be sure to replace the media at the same time each day.
- 17 **Day 4 and 6:** At 96h and day 6, aspirate medium from wells and add 5 ml of CDM3 complete medium.
- 18 **Day 8 and 10:** Selection media: RPMI 1640 0.5 L + B27 supplement 10 ml. Selection media is stable at 4 °C for >2 weeks. Aspirate medium from wells and add 2 ml of selection media.
- 19 NOTE: Metabolic purification can be started at day 8 or at first onset of spontaneous beating. At least 4 days of metabolic selection may be required. Some non-cardiomyocyte cell types can survive metabolic selection. Supplementation with 4ug/ml



insulin and 1 mg/ml polyvinyl alcohol greatly improves differentiation purity and typically does not require selection.

20 **Day 12:** Harvest CMs for scRNA-seq

Part 10: Protocol for cell hashing followed by 10X single-cell RNA-Seq

21 **Goal:** Harvest cells and perform cell hashing for superloading of 10X RNA-Seq libraries

21.1 Harvest differentiated cardiomyocytes into single-cell suspension.

21.2 Resuspend 10M - 13M harvested cells in 1ml staining buffer (1 % BSA in 50% PBS / 50% culture media) in 1.5 ml Eppendorf tube

21.3 Dilute 1ul (0.5 ug) TotalSeqTM-C anti-human Hashtag Antibody (BioLegend) to 5 ul staining buffer in 0.2 ml PCR strip on ice. Repeat for a total of 12 unique Hashtag antibodies.

21.4 Aliquot 50 ul of cell suspension to each Hashtag-Ab PCR tube, mix and incubate at 4°C with rotation for 20 min.

21.5 Add 150 ul staining buffer to each tube and centrifuge at 275 x g for 5 min.

21.6 Remove supernatant and add 200 ul staining buffer to each tube. Mix and transfer cell suspensions to a 15 ml tube with 5 ml staining buffer on ice. Centrifuge the 15 ml tube at 275 g for 5 min.

21.7 Wash the cells in the 15 ml tube with a 7 ml ice-cold staining buffer once.

21.8 Wash the cells in the 15 ml tube with a 7 ml ice-cold 1% BSA in PBS once. Count cell concentration to determine how much volume should be added in the next step.

21.9 Resuspend cells in ice-cold 1% BSA in PBS at a concentration of 5 ~ 6 M cells/ml.

21.10 Filter the cell suspension with 30 um Pre-Separation filters (Miltenyi Biotec).



- 21.11 Count the cells. Ensure cell viability is >70% and cells are in a single-cell suspension.
- 21.12 There are no super-loading 10X 5' HT kits for CRISPR Screening & Cell Multiplexing available in the market. 5' HT libraries were prepared with 10X 5' HT reagent kits with Feature Barcode technology for CRISPR Screening & Cell Surface Protein. But super-loading the cells in protocol CG000513 part **1.2 Load Chromium Next GEM Chip N** (shown below):
1. For cardiomyocytes, load 173k cells to target 60k recovered cells
 2. For other stem-cell derived cells, load 115k cells to target 60k recovered cells

Part 11: Construction of scRNA-Seq libraries

- 22 **Goal:** Construct sgRNA and HTO libraries during scRNA-Seq.
- 22.1 Prepare transcriptome, sgRNA, and HTO libraries according to 10X CG000513 protocol with some modifications
- 22.2 For transcription libraries follow the protocol CG000513 **Step 5: 5' Gene Expression Library Construction**.
- 22.3 For sgRNA libraries follow the protocol CG000513 **Step 6: CRISPR Screening Library Construction**, except 4 times of the input sgRNA cDNA recommended were used:
1. CG000513 step 6.1.j ~ 6.1.m: Elute with 20.5 ul instead of 40.5 ul.
 2. CG000513 part 6.2.b: Add 10 ul (instead 5 ul) Guide RNA cDNA Cleanup for each PCR reaction. Add 40 ul (instead 45 ul) Feature S1 primers 4 for each PCR reaction.
 3. CG000513 part 6.2.e and 6.4.d Increase the PCR cycles from 10 to 11.
- 22.4 For HTO libraries follow the protocol CG000513 **Step 7: Cell Surface Protein Library Construction** without modification. TotalSeq™ anti-human Hashtag reagents (BioLegend) are designed to label most human cells, using a combination of two clones that recognize cell surface protein CD298 and β 2 microglobulin, respectively.
- 22.5 Measure the library concentration with Qubit reagents(Thermo Fisher Scientific) and analyze the peak size of the libraries with TapeStation System(Agilent Technologies).