

Mar 03, 2025

Extracellular Vesicle Isolation

DOI

<https://dx.doi.org/10.17504/protocols.io.n92ldr4eng5b/v1>

Emily Tabaie¹

¹University of California, Riverside



Emily Tabaie

University of California, Irvine

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n92ldr4eng5b/v1>

Protocol Citation: Emily Tabaie 2025. Extracellular Vesicle Isolation. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.n92ldr4eng5b/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 03, 2025



Last Modified: March 03, 2025

Protocol Integer ID: 123711

Keywords: extracellular vesicle isolation these step, extracellular vesicle isolation, extracellular vesicles from primary neuronal cell culture, extracellular vesicle, primary neuronal cell culture, exoeasy kit from qiagen, isolation, cell, qiagen, exoeasy kit

Abstract

These steps are used to isolate extracellular vesicles from primary neuronal cell cultures. It uses a combination of centrifugation and a ExoEasy Kit from Qiagen for isolation.

Extracellular Vesicle Isolation

- 1 Do all steps in the hood and be sterile
- 2 Grow Neurons and ME49 B7 parasites separately
- 3 Change Neuron media every 3 days and save at -80°C
- 4 At Day 9 of the neurons purify the parasites and infect the neurons with an MOI of 0.5
- 5 After infection wait 3 hours and then change half of the neuronal media
- 6 Continue change and saving the neuron media until the neurons die
- 7 Once all the supernatant is collected and filtered centrifuge at 500 x g for 15 minutes at 4°C
- 8 Take out the supernatant and ultracentrifuge it at 15,000 x g for 20 minutes at 4°C
- 9 Use the exoEasy Maxi Kit (Qiagen 76064)
- 10 Add 1 volume buffer XBP to 1 volume of sample and mix well by gently inverting the tube 5 time
- 11 Let the mixture warm up to room temperature
- 12 Add the sample/XBP mix onto the exoEasy spin column and centrifuge at 500 x g for 1 min at room temperature
- 13 Discard the flow-through and place the column back into the same collection tube



- 14 Add 10 mL buffer XWP and centrifuge at 4,000 rpm for 5 min at room temperature to remove residual buffer from the column
- 15 Discard the flow-through together with the collection tube
- 16 Transfer the spin column to a fresh collection tube
- 17 Add 400 µl of Buffer XE to the membrane and incubate for 1 min
- 18 Centrifuge at 500 x g for 5 min at room temperature to collect the eluate
- 19 Re-apply the eluate to the exoEasy spin column membrane and incubate for 1 min
- 20 Centrifuge at 4,000 rpm for 5 min at room temperature to collect the eluate and transfer to an appropriate tube (not supplied)
- 21 Aliquot out the EVs and store at -80°C
- 22 Run NTA on the EVs the same day before freezing EVs