

Apr 03, 2017

🌐 Fluorescence activated cell sorting (FACS) of *Perkinsus marinus* transformants

📖 [Nature Methods](#)

DOI

dx.doi.org/10.17504/protocols.io.hh2b38e

Protist Research to Opti...



Imen Lassadi

OPEN ACCESS



DOI: dx.doi.org/10.17504/protocols.io.hh2b38e

External link: <https://doi.org/10.1038/s41592-020-0796-x>

Protocol Citation: Imen Lassadi: Fluorescence activated cell sorting (FACS) of *Perkinsus marinus* transformants. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.hh2b38e>

Manuscript citation:

Faktorová D, Nisbet RER, Robledo JAF, Casacuberta E, Sudek L, Allen AE, Ares M, Aresté C, Balestreri C, Barbrook AC, Beardslee P, Bender S, Booth DS, Bouget F, Bowler C, Breglia SA, Brownlee C, Burger G, Cerutti H, Cesaroni R, Chiurillo MA, Clemente T, Coles DB, Collier JL, Cooney EC, Coyne K, Docampo R, Dupont CL, Edgcomb V, Einarsson E, Elustondo PA, Federici F, Freire-Beneitez V, Freyria NJ, Fukuda K, García PA, Girguis PR, Gomaa F, Gornik SG, Guo J, Hampl V, Hanawa Y, Haro-Contreras ER, Hehenberger E, Highfield A, Hirakawa Y, Hopes A, Howe CJ, Hu I, Ibañez J, Irwin NAT, Ishii Y, Janowicz NE, Jones AC, Kachale A, Fujimura-Kamada K, Kaur B, Kaye JZ, Kazana E, Keeling PJ, King N, Klobutcher LA, Lander N, Lassadi I, Li Z, Lin S, Lozano J, Luan F, Maruyama S, Matute T, Miceli C, Minagawa J, Moosburner M, Najle SR, Nanjappa D, Nimmo IC, Noble L, Vanclová AMGN, Nowacki M, Nuñez I, Pain A, Piersanti A, Pucciarelli S, Pyrih J, Rest JS, Rius M, Robertson D, Ruaud A, Ruiz-Trillo I, Sigg MA, Silver PA, Slamovits CH, Smith GJ, Sprecher BN, Stern R, Swart EC, Tsaousis AD, Tsy-pin L, Turkewitz A, Turnšek J, Valach M, Vergé V, Dassow Pv, Haar Tvd, Waller RF, Wang L, Wen X, Wheeler G, Woods A, Zhang H, Mock T, Worden AZ, Lukeš J, Genetic tool development in marine protists: emerging model organisms for experimental cell biology. *Nature Methods* 17(5). doi: [10.1038/s41592-020-0796-x](https://doi.org/10.1038/s41592-020-0796-x)

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: March 31, 2017

Last Modified: March 16, 2018

Protocol Integer ID: 5402



Cell Culture and Electroporation

- 1 See protocol "Oyster parasite *Perkinsus marinus* transformation using Amaxa electroporator and non-proprietary electroporation buffer"
Monitor transfection efficiency, by testing for the presence of fluorescent cells 5 to 6 days post transfection.

Cell recovery

- 2 Once the transfection success is confirmed, transfer the cells to a larger volume in a T75 flask (or equivalent) to allow cell number to increase (up to one week).

Fluorescence activated cell sorting protocol

- 3
The experiment should be undertaken in sterile conditions

Use *Perkinsus marinus* wild type cells as a control to set gating for non-fluorescent, single, live cells.

Set up a template that includes a bivariate plot to display forward scatter (FSC) and side scatter (SSC), and one histogram for each fluorophore that will be used (e.g. eGFP, mCherry)

Gate cells for fluorescence above that seen for the untransformed control cells, and sort these cells either as population or as single cells to a 96 well plate.

Sorted single cells in 96 well plates can then be cultured for 2 months at 25°C in the dark.

Inspect cells by microscopy to confirm fluorescence status.