



Jul 11, 2024

Fluorescent in situ hybridisation for juvenile *Fasciola hepatica*

 [PLOS Pathogens](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.3byl4q7xjvo5/v1>

Duncan Wells¹, Paul McCusker², Rebecca Armstrong², Emily Robb²

¹QUB; ²Queen's University Belfast



Duncan Wells

QUB

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl4q7xjvo5/v1>

Protocol Citation: Duncan Wells, Paul McCusker, Rebecca Armstrong, Emily Robb 2024. Fluorescent in situ hybridisation for juvenile *Fasciola hepatica*. [protocols.io](https://dx.doi.org/10.17504/protocols.io.3byl4q7xjvo5/v1) <https://dx.doi.org/10.17504/protocols.io.3byl4q7xjvo5/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: July 03, 2023

Last Modified: July 11, 2024

Protocol Integer ID: 84389

Keywords: juvenile fasciola hepatica this fish protocol, use in fasciola hepatica, juvenile fasciola hepatica, fasciola hepatica, fluorescent in situ hybridisation, fluorescent, fish protocol, lab, situ hybridisation, newmark lab

Abstract

This FISH protocol has been adapted for use in *Fasciola hepatica*, specifically 4-week-old juveniles, grown *in vitro* in 50% Chicken Serum, 50% RPMI culture medium. The majority of this protocol comes from those of the Newmark Lab (King and Newmark, 2013) and Collins Lab (Collins and Collins, 2017). We are grateful to members of these labs for several helpful discussions on this technique.

Guidelines

Following a long period of relative success with this technique, we had noticed difficulties emerged. After trying various tweaks, we decided to purchase the vast majority of reagents again. This seemed to have helped significantly.

For all steps prior to and including the riboprobe incubation, all tips and tubes should be RNase free. Reagents should be made up with DEPC-treated double distilled H₂O - for 1 litre of DEPC-treated water, add 1ul Diethyl pyrocarbonate (D5758, Sigma) per 1ml of H₂O. Shake vigorously for 15 seconds. Leave at room temperature overnight, or 2 hours at 37°C, then autoclave.

We have found that longer riboprobes (800-1000 bases) have generally worked better than shorter (150-300 bases) ones.

Alongside every positive test, we run the sense probe as a negative control. This is especially helpful when we see probes with otherwise ambiguous/questionable staining.

Highlighted solutions have recipes at the bottom of the protocol.

Safety warnings

! Various reagents through this protocol have hazards associated with them. Consider these before use.

Before start

Ensure in situ oven is at the desired temperature. We usually switch the oven on at the start of the day 1 and allow the oven to reach a stable temperature.



Riboprobe Preparation

- 1 Amplify target of interest with an end-point PCR using FastStart Taq polymerase, dNTPack (4738381001, Sigma). Primers should have a 5'-TAATACGACTCACTATAGGGT-3' promoter sequence. For the antisense probe (positive) use the regular forward primer and the T7-tagged reverse primer. For the sense probe (negative control) use the T7-tagged forward and the regular reverse. 2h 30m
 - 2 Purify resulting PCR products with the ChargeSwitch PCR clean-up kit (CS12000, Thermo) according to the manufacturers guidelines. Only 15µl of elution buffer is used to maximise amplified target concentration. 20m
 - 3 Perform a transcription reaction containing 2µl transcription buffer, 1µl T7 polymerase (both from T7 polymerase kit, EP01111), 6µl purified PCR product, 1µl DIG-RNA labelling mix (11277073910, Sigma). Incubate Overnight at 27 °C . 20m
 - 4 DNase treat the probes by adding 0.25µl DNase I (18047019, Thermo), 0.25µl transcription buffer (from previous step) and 2µl RNase-free H₂O. Incubate at 37 °C for 00:15:00 . 15m
 - 5 Add 1.25µl **4M LiCl** and 37.6µl molecular-grade ethanol (51976, Sigma). Precipitate probes Overnight at -80 °C . 5m
 - 6 16000 x g, 4°C, 00:20:00 20m
 - 7 Remove supernatant, wash pellet with 70% ethanol, 16000 x g, 4°C, 00:05:00 , discard supernatant and allow to air-dry for 00:15:00 . 20m
- Expected result**

Pellet is barely visible to the eye. We would recommend removing supernatants while viewing the pellet under a microscope.
- 8 Resuspend pellet in 20µl RNase-free H₂O. Check purity and concentration on a spectrophotometer/fluorometer (we use a DeNovix DS-11 FX). 20m



- 9 Adjust concentration to 50ng/μl with **Hybridisation Buffer**, store probes at -20 °C until use.

10m

Preparation of yeast RNA (mainly from Jiang [2012])

- 10 Add 1g of torula RNA (RNA from torula yeast, type VI [R6625, Sigma]) to a 50ml conical tube.

5m

- 11 Add 10ml phenol (P4557, Sigma) and 15ml Tris-EDTA buffer (12090015, Thermo). Vortex tubes until the yeast RNA has dissolved. 16000 x g, 4°C, 00:10:00

10m

- 12 Collect the aqueous phase into a new tube (50ml), add 1:1 phenol/chloroform (P2069, Sigma), mix well, 16000 x g, 4°C, 00:10:00 .

10m

- 13 Collect aqueous phase into a new tube, add 1:1 chloroform, mix well, 16000 x g, 4°C, 00:10:00

10m

- 14 Collect aqueous phase into a new 50ml tube, add 1ml 3M sodium acetate (71196, Sigma) and 30ml molecular-grade ethanol. Leave at -20 °C Overnight to precipitate.

- 15 16000 x g, 4°C, 00:20:00 , discard supernatant, wash pellet in 70% ethanol, 16000 x g, 4°C, 00:10:00

30m

- 16 Remove supernatant, allow pellet to airdry for 00:15:00 , add 5ml formamide pre-heated to 50 °C

15m

- 17 Vortex until the pellet is resuspended

- 18 Measure concentration and purity on a spectrophotometer/fluorometer. Adjust concentration to 50mg/ml and store in 1ml aliquots at -20 °C .

Worm Preparation

1d 0h 40m



19 Carry out excystment (see McVeigh *et al.*, 2014), transfer newly excysted juveniles to 96-well round-bottomed plates, in groups of <25 NEJs per well. Add 200µl of pre-warmed, pre-mixed 50:50 chicken serum (New Zealand Origin, 16110082, Thermo) and RPMI (11835030, Thermo). Replace this every other day. Worms are maintained in a  37 °C incubator with 5% CO₂.

1d

20 At 4 weeks old, wash juveniles out of 50:50 and add RPMI. Transfer worms into a small petri dish (Starstedt – 83.3925) containing **0.25% Tricaine** diluted in RPMI. Leave worms in this for  00:03:00 .

3m

Expected result

Worms should initially curl ventrally but this should be greatly reduced following this incubation.

21 Transfer worms in a minimal amount of tricaine/RPMI into a 20µl droplet of **4% Formaldehyde** and set a glass coverslip on top. Leave for  00:10:00 . Remove the coverslip and collect worms into 1.5ml hydrophobic tubes (1210-10, SSIBio). Add 1ml 4% formaldehyde. Rotate tubes for  00:10:00 .

20m

22 Remove 4% formaldehyde and wash once with **PBST**, then replace this with a pre-mixed 1:1 ratio of PBST and methanol (32213, Sigma) rotating for  00:10:00 . Worms can then be stored in 100% methanol at  -20 °C until use. We are yet to test how long worms can be stored for and still be suitable. We would be reluctant to use worms stored for longer than 6 months.

10m

Day 1 FISH (each step begins with removing the solution from the last) (volumes are 1ml unless otherwise stated)

23 Remove 100% methanol from worms and replace with 1:1 PBST:Methanol for  00:10:00 rotating.

10m

24 Add PBST, incubate for  00:10:00 rotating

10m



- 25 Add **1xSSC** and incubate for  00:10:00 with no rotation 10m
- 26 Incubate worms in 200µl **Bleaching solution** under a bright light for  01:00:00 1h
- 27 Add 1xSSC and incubate for  00:10:00 with no rotation 10m
- 28 Add PBST, incubate for  00:10:00 rotating. 10m
- 29 Add **Proteinase K solution**, incubate for  00:15:00 with no rotation. 15m
- 30 Add 4% formaldehyde, incubate for  00:10:00 rotating. 10m
- 31 Add PBST
- 32 Transfer worms in 500µl PBST into in situ baskets (12.444, Intavis) in a 24-well plate. Add 500µl **Prehybridisation buffer**.
- 33 Add plate to a shaker for  00:10:00 set at 100rpm (any future shaking steps are at 100rpm) 10m
- 34 Add Prehybridisation buffer (pre-warmed to  52 °C) to the next well and transfer the baskets to these. Place the plate in an in situ oven pre-warmed to  52 °C . Have the plate rocking. Our in situ oven (Boekel Scientific Shake 'N Bake™ Rocking Hybridization Oven/Incubator, 136400) only allows for gentle rocking.
- 34.1 Approximately 10 minutes before the end of the above 2 hours, prepare the probe as follows - heat 20µl of 50ng/µl probe to  80 °C for  00:05:00 then place on ice. Add 980µl Hybridisation buffer to the probe. 5m
- 35 Add the probe mix to the next well, transfer baskets to this well and incubate  Overnight rocking at  52 °C .



Day 2 FISH (volumes are 1ml unless otherwise stated) (each basket wash is carried out in a new well unless otherwise stated)

- 36 Remove 500µl of probe mix from the well with the worms and add 500µl pre-heated **2xSSC**, incubate for 00:20:00 rocking at 52 °C . 20m
- 37 Wash in 2xSSC, incubate for 00:20:00 rocking at 52 °C . 20m
- 38 Repeat step 36) 2 more times, moving baskets into new wells each time.
- 39 Wash with **0.2xSSC**, incubate for 00:20:00 rocking at 52 °C 20m
- 40 Repeat step 38) 3 more times, moving baskets into new wells each time.
- 41 Wash with **TNT**, incubate for 00:10:00 shaking at Room temperature 10m
- 42 Repeat step 41) once
- 43 Incubate in **Blocking solution** for 01:00:00 shaking at Room temperature 1h
- 44 Incubate in **Antibody solution** Overnight shaking at 4 °C 1h

Day 3 FISH (volumes are 1ml unless otherwise stated)

- 45 Wash with TNT for 00:05:00 , then 00:10:00 , and then 6x 00:20:00 shaking at Room temperature . 35m
- 46 Protect plates from light from this point on. 4h 10m

**FOR FLUORESCENT ISH**

Incubate in freshly-made **Tyramide solution** for  00:10:00 , shaking at

 Room temperature . Go to step 46).

FOR CHROMOGENIC ISH

Transfer worms to small petri dishes and incubate in **Exposure buffer** for up to

 04:00:00 with no rotation, checking every 20 minutes or so. If no signal has

developed after 4 hours, incubate the dishes at  4 °C . This slows any subsequent development but this allows you to leave it overnight. When a sufficient signal has developed, halt the reaction by washing worms with TNT. Wash the worms in 100% ethanol to reduce any background/non-specific staining. Wash worms with TNT. Go to step 47).

- 47 Incubate worms in **DAPI solution** for  00:20:00 , followed by two 10 minute washes in TNT.

20m

- 48 Mount worms in 10µl Vectashield (H-1000, Vector Labs) on standard microscope slides and seal coverslips with clear nail polish.

Recipes (presented in order of first appearance) (DEPC-treated H₂O presented as DT-H₂O, double distilled water presented as ddH₂O)

2h

- 49 **4M LiCl** (100ml)
Add 16.96g LiCl (L9650, Sigma) to 50ml DT-H₂O. Add this slowly as this is an exothermic reaction. Once dissolved, bring volume to 100ml with ddH₂O.

2h

Hybridisation Buffer (50ml)

25ml deionised formamide (F9037, Sigma)

12.5ml 20xSSC (S6639, Sigma)

100µl 1mg/ml yeast RNA

5ml 10% (v/v) Tween-20 (5ml Tween-20 [P1379, Sigma] + 45ml DT-H₂O)

5ml 50% Dextran sulfate (S4030, Sigma)

2.4ml DT-H₂O

0.25% Tricaine

50mg Ethyl 3-aminobenzoate methanesulfonate (E10521, Sigma) in 2ml RPMI

4% Formaldehyde (4.5ml)

500µl 36.5% formaldehyde (F8776, Sigma)
4ml PBST (see below)

PBST (200ml)

1 tab PBS (P4417, Sigma)
200ml DT-H₂O
600µl Triton-X (T8787, Sigma)

1xSSC (20ml)

1ml 20xSSC
19ml DT-H₂O

Bleaching solution (2ml)

1.77ml DT-H₂O
100µl deionised formamide
50µl 20xSSC
80µl H₂O₂ (H1009, Sigma)

Proteinase K solution (~5ml)

4.95ml PBST
50µl 10% SDS (dissolve 5g SDS [L4509, Sigma] in 50ml DT-H₂O)
2.5µl Proteinase K (03115887001, Sigma) (check the bottle for exact stock concentration as this varies, adjust calculation accordingly. End concentration is to be 10ug/ml)

Prehybridisation buffer (50ml)

25ml deionised formamide
12.5ml 20xSSC
100µl 50mg/ml yeast
5ml 10% (v/v) Tween-20
7.4ml DT-H₂O

2xSSC (20ml)

17.98ml ddH₂O
2ml 20xSSC
20µl Triton-X

0.2xSSC (20ml)

19.78ml ddH₂O
200ul 20xSSC
20µl Triton-X

**TNT (~1L)**

1L ddH₂O

12.11g Tris Base (T1503, Sigma)

8.77g NaCl (27800.291, VWR)

3ml Triton-X

pH to 7.5 and filter sterilise

Blocking Solution (8ml)

400µl Horse serum (H0146)

40µl Western Blocking Solution (11921673001, Sigma)

7.56ml TNT

Antibody solution (4ml)

2µl Anti-DIG-AP (Chromogenic ISH) (11093274910, Sigma) / 3µl Anti-DIG-HRP

(Fluorescent ISH) (11207733910, Sigma)

3998µl Blocking Solution

Tyramide Solution (~5ml)

5ml TSA Buffer (see below)

1µl 20mg/ml 4-iodophenylboronic acid (see below)

6µl 5% H₂O₂ (492µl TSA Buffer + 8µl H₂O₂)

2µl FAM-Tyramide (see below)

TSA Buffer (200ml)

58.44g NaCl

3.09g Boric acid (B6768, Sigma)

100ml ddH₂O, bring to 200ml when dissolved

pH to 8.5 and filter sterilise

20mg/ml 4-iodophenylboronic acid (2ml)

40mg 4-iodophenylboronic acid (IPBA, 471933, Sigma) in 2ml dimethylformamide (DMF, 227056, Sigma)

FAM-Tyramide

- Tyramine solution = Dissolve 10mg tyramine hydrochloride (T2879, Sigma) in 1ml DMF, add 10µl triethylamine (T0886)

- FAM-DMF = Dissolve 10mg 5(6)-FAM, SE (5-(and-6)-Carboxyfluorescein, Succinimidyl Ester) (C1311, Thermo) in 1ml DMF

Add 342.5µl tyramine solution to 1ml FAM-DMF

Incubate at for  02:00:00 , rotating at room temperature

Add 8.6ml ethanol for a 1mg/ml stock



Split into 20µl aliquots and store at  -20 °C

Exposure Buffer (~10ml)

Prepare 50ml tubes of each of;

- 1M Tris base = 6.06g
- 1M NaCl = 2.92g
- 1M MgCl₂ = 4.76g

Bring each to 50ml with ddH₂O

Prepare PVA - bring 80ml ddH₂O to  75 °C while stirring and slowly add 10g PVA (P8136, Sigma). If you add this too quickly it will form clumps and not dissolve.

1ml 1M Tris base

1ml 1M NaCl

500µl 1M MgCl₂

100µl 10% Tween-20

7.4ml PVA

45µl 4-Nitro blue tetrazolium (NBT, 11383213001, Sigma)

35µl 5-bromo-4-chloro-3-indolyl-phosphate (BCIP, 11383221001, Sigma)

DAPI Solution (10ml)

10µl 1mg/ml DAPI (dissolved in ddH₂O) (D9542, Sigma)

9.99ml TNT